



Own brands

range

Excellent value for money



Alexandra Takacs
Marketing & Administration Manager

Welcome!

Dear HENRY SCHEIN customer, dear HENRY SCHEIN customer,

I am pleased to present you our new Henry Schein private label catalog 2023/2024 to introduce. Discover our extensive own-brand range on over 65 product pages, which we continually add new products for you.

You can recognize our own brand Henry Schein by its quality seal. This stands for a large selection, best quality and consistently low prices.

This time we paid particular attention to our expanded disinfection range, with our new EuroSept Xtra product line. In times when hygiene and cleanliness are top priority, we have carefully expanded our range of disinfectants, cleaning products and protective items to make your work environment as safe as possible so that you can concentrate on what is important: providing the best possible care for your patients.

Browse through our new catalog and be inspired by the diverse products. Thank you for your trust in Henry Schein.

Sincerely, yours

Alexandra Takacs
Marketing & Administration Manager



At henryscheinmed.at you will constantly find new products and offers!

- > Access to over 10,000 articles
 - > Everything at a glance with your customer account (Invoice, order history and delivery tracking)
 - > No waiting time on the phone
- > View invoice as PDF
 - > Register the return and contact the person responsible
 - Find field staff
 - > Create watch lists

PRODUCT OVERVIEW

Injection/infusion	04	Instruments	42
Examination gloves	08	Instrument disinfection/processing	43
Surgical requirements	12	dressing materials	46
Cleaning/disinfection	15	Practical / laboratory items	53
Sanitary papers	33	Rapid tests	56
Scalpels	36	Diagnosis	58
Scissors	37	ECG/ultrasound	64
Needle holder	39	X-ray/dental care	65
Tweezers	40	Nursing/elderly care	67
Arterial clamps	41		

SERVICE

About Henry Schein	69	Technical service/repairs/ Maintenance	78
team	70-71	Device checks	79
Webshop	72	Probe analyses	80
Newsletter	73	Leasing/Financing	81
Practice Green	74-75	Practice planning	82
Henry Schein own brand	76-77	All from a single source	83



100 pieces from

3.29

PRICE-PERFORMANCE

Single-use syringes

With Luer connector, 2 parts

Particularly inexpensive standard syringes, easy to vent, pronounced piston stop, with Luer connection or Luer lock connection. Manufactured according to ISO standards, pyrogen-free, individually sterile packed. With an easy-to-read and smudge-proof scale. Pack (100 pieces)

Size	Order no.	Price
2ml/3ml	900 3017	3.60
from 10 packs		3.49
each from 25		3.29
packs each 5ml/6ml	900 3018	5.60
from 10 packs of		5.49
20 packs of 10ml/		5.19
12ml each	900 3019	7.95
from 12 packs each		7.85
from 24 packs each		7.69
20ml/ 24ml	900 3020	10.95
from 8 packs each		10.75
from 24 packs each		9.95



Disposable insulin syringes

The ideal syringe for dosing small units. The smooth-gliding piston allows even extremely slow drawing and injection. Tested according to DIN 13098.

Easy-to-read, smudge-proof scale. Minimum remaining volume. Individually sterile packed. Pack (100 pieces)

Size	Order no.	Price
With cannula U10030 G x1/2" (0.3 x12mm) 1m	900 8529	15.95
from 24 packs each		14.99
With cannula U40 30 G x1/2" (0.3 x12mm) 1ml	900 8530	15.95
from 24 packs each		14.99



Single-use syringes

3 parts with rubber pistons

Comfortable handling with even easier running and better readability.

- > Excellent quality
- > Smudge-proof scale and transparent cylinder
- > Extremely easy to run and absolutely tight
- > Sterile in individual blister packaging
- > Latex-free, pyrogen-free, non-toxic and PVC-free
- > Environmentally friendly packaging

Pack (100 pieces)

1 disposable syringe with Luer connector

Size	Order no.	Price
2 ml / 3 ml (centr.)	900 3302	5.45
from 5 packs of 5 ml		4.95
each (eccentric)	900 3311	8.25
from 5 packs each		7.95
10 ml (eccentric)	900 3304	10.50
from 5 packs of 20		9.95
ml each	900 3305	15.50
(eccentric) from 5 packs each		14.95

2 disposable syringes with Luer lock attachment

Size	Order no.	Price
3ml/0.1ml	900 3759	13.90
5ml/0.2ml	900 3760	16.90
10ml/0.5ml	900 3761	19.50



100 pieces from

9.55

PRICE-PERFORMANCE

3-part fine dosing syringes

1 ml, 3 parts, for tuberculin/heparin/allergy examinations. Luer approach, centric. Highly transparent cylinder. Black graduation. Smooth-running piston plug with double sealing ring. Minimum remaining volume. PVC free. Latex free. Pack (100 pieces)

Pack of	900 3310	10.95
18 packs each		9.55



100 pieces from

12.95

PRICE-PERFORMANCE


Tuberculin disposable syringes

With Spardorn, 3 parts

1ml; 100 pitch with central Luer connection made of plastic. No loss of medication due to Spardorn. Individually sterile packed.

In the donor box. Pack (100 pieces)

Cardboard	900 8276	13.75
from 20 packs each		12.95



100 pieces from

9.95

PRICE-PERFORMANCE

Tuberculin disposable syringes

With Spardorn, 2 parts

1ml; 100 pitch with central Luer connection made of plastic. No loss of medication due to displacer piston (saving mandrel).

Individually sterile packed. In the donor box. Pack (100 pieces)

Cardboard	900 3016	10.95
from 18 packs each		9.95



Piece off

1.69

PRICE-PERFORMANCE

Disposal boxes

Particularly robust safety containers for disposing of used cannulas, syringes, pipette tips and scalpel blades.

With twisting and pulling devices for cannulas in the lid (for all cannula types). Once sealed in the end cap, the container can no longer be opened. Burns without pollutant emissions.

Size	Order no.	Price
1 0.25 l	988 1898	1.69
2 0.5L	988 1899	1.99
3 1 liter	988 1900	2.29
4 1.5L	988 1901	2.49
5 2 liter	988 1907	2.49
6 2.2L	988 1902	2.49
7 3 liter	988 1908	2.99
8 3.5L	988 1903	3.29
9 4 liter	988 1904	3.49
10 6 l	988 1905	4.49
11 10 liter	988 1906	5.99



Disposable cannulas

Standard sizes

- > Excellent quality
 - > Stainless steel according to ISO 9626/1991
 - > Precision tip according to ISO 9626 with precise, smooth cannulas cut
 - > Base made of transparent polypropylene
 - > Polypropylene protective cap
 - > Thin-walled siliconized cannula tube
 - > For Luer and Luer lock syringes
 - > Sterile, pyrogen-free and non-toxic
 - > 5 years expiry time for sterile goods
 - > Individually sterile packed
- Pack (100 pieces)

100 pieces from

2.69

PRICE-PERFORMANCE

Size	Ø x length mm	Color	Order no.	Price
iv, for highly viscous solutions, im				
Gr. 1	0.90x40mm	yellow	900 3333	2.99
iv, large volume, in, aqueous solutions				
Gr. 2	0.80x40mm	green	900 3334	2.99
iv, in				
Gr. 12	0.70x30mm	black	900 3335	2.99
iv, in small quantities				
Gr. 14	0.60x30mm	blue	900 3336	2.99
iv, in small quantities				
Gr. 16	0.60x25mm	light blue	900 3337	2.99
iv, sc, hand vein, pediatrics				
Gr. 17	0.55x25mm	purple	900 3338	2.99
iv, sc, hand vein, pediatrics				
Gr. 18	0.45x25mm	brown	900 3339	2.99
Insulin, in				
Gr. 20	0.40x20mm	Gray	900 3340	2.99
from 10 packs per variety				2.69

Disposable cannulas

Special sizes

Pack (100 pieces)


Ø x length mm	Color	Order no.	Price
0.40x40mm	Gray	900 6444	4.55
0.50x16mm	orange	900 3629	4.55
0.80x16mm	green	900 3628	4.55
0.80x25mm	green	900 3627	4.55
0.80x50mm	green	900 6443	4.55
0.90x25mm	yellow	900 3626	4.55
1.10x25mm	beige	900 3633	4.55
long tip 1.20 x 40 mm	pink	900 3625	4.55
short tip 1.20 x 40 mm	pink	900 3624	4.55
1.60x25mm	transparent	900 3631	4.55
1.60x40mm	transparent	900 3632	4.55
from 10 packs per variety			3.95



100 pieces from

3.95

PRICE-PERFORMANCE



100 pieces from

16.90

PRICE-PERFORMANCE

Protect yourself safely from
needlestick injuries!

Safety disposable cannulas

The disposable safety cannulas are closed with one hand; an audible "click" indicates whether the safety mechanism is closed. The cannulas can be used with standard Luer and Luer Lock syringes.

Pack (100 pieces)

Gauge / Ø x length mm	Color	Order no.	Price
G 21 / 0.8 x 25 mm	green	900 7650	17.95
G 23 / 0.6 x 25 mm	blue	900 7654	17.95
G 25 / 0.5 x 16 mm	orange	900 7657	17.95
G 25 / 0.5 x 25 mm	orange	900 7655	17.95
G 25 / 0.5 x 38 mm	orange	900 7656	17.95
G 30 / 0.30 x 12.5 mm	yellow	900 7658	17.95
from 10 packs per variety			16.90

Pack off

75.90

PRICE-PERFORMANCE

Latex free



Infusion devices for
Gravity infusion

Without intermediate piece

According to DIN 58362 for gravity infusion. With ISO mandrel. Free from pyrogen and DEHP. With bacteria-proof air filter. Hose length 180cm, Luer lock attachment. Individually sterile packed.


Pack (100 pieces)

Pack	988 3618	78.90
from 5 packs		77.90
each from 10 packs each		75.90

Pack off

78.90

PRICE-PERFORMANCE



Infusion devices for
Gravity infusion

With flexible intermediate piece

According to DIN 58362 for gravity infusion. With ISO mandrel and flexible intermediate piece. Free from pyrogen and DEHP. With bacteria-proof air filter. Hose length 180cm, Luer lock connection. Individually sterile packed. Pack (100 pieces)

Pack	988 3617	83.90
from 5 packs		79.90
each from 10 packs each		78.90

1 cellulose swab from the roll. Fits into every dispenser

box. Ready to use. Soft and absorbent quality. Can be sterilized in steam. Swab size: 4x5cm

Size	Order no.	Price
Pack (2 rolls of 500 pieces)	900 0853	5.75
from 12 packs each		5.45
Swab dispenser		
Piece	420 499	14.90

2 alcohol wipes

Soaked with 70% isopropyl alcohol. Separately packaged.

Dimensions: 65x 30 mm

Pack (200 pieces)	104 6352	4.45
from 20 packs		4.29
each from 40 packs each		3.79

3 alcohol wipes

Alcohol swabs in the practical dispenser box.

The empty dispenser box can be used as a cannula collector.

Pack (150 pieces)	900 6440	4.95
from 10 packs each		4.75



1,000 pieces

5.45

PRICE-PERFORMANCE




200 pieces



2

3



With quick clip opening

Latex free

piece only

5.79

PRICE-PERFORMANCE

tourniquet

- > Easy to put on
- > Can be operated with one hand
- > To decongest in a measured manner
- > Skin-friendly dust band
- > Without metal — also for radiological environments

Piece	900 6441	6.19
from 5 pieces each		5.89
from 10 pieces each		5.79



100 pieces from

8.90

PRICE-PERFORMANCE



Premium latex examination

gloves High-quality examination gloves with a very good fit and excellent tactile sensitivity. Made of particularly tear-resistant latex. Corresponds to DIN EN 455.

Pack (100 pieces)

Size	Order no.	Price
Powder free		
XS (5-6)	101 9924	9.90
S (6-7)	101 8787	9.90
M (7-8)	101 3806	9.90
L (8-9)	101 0445	9.90
XL (9-10)	900 1715	9.90
from 10 packs per variety from		9.50
30 packs per variety each		8.90
Lightly powdered		
XS (5-6)	101 6863	9.50
S (6-7)	101 2710	9.50
M (7-8)	101 7034	9.50
L (8-9)	101 5682	9.50
XL (9-10)	900 1748	9.50
from 10 packs per variety from		8.90
30 packs per variety each		8.50

100 pieces from

7.90

PRICE-PERFORMANCE



Economy Latex Examination

Gloves Inexpensive latex examination gloves.

Corresponds to DIN EN 455.
Pack (100 pieces)

Size	Order no.	Price
Powder free		
XS (5-6)	102 4184	8.90
S (6-7)	102 8304	8.90
M (7-8)	102 5289	8.90
L (8-9)	102 9168	8.90
from 10 packs per variety from		8.50
30 packs per variety each		7.90
Lightly powdered		
XS (5-6)	101 0954	8.90
S (6-7)	101 5137	8.90
M (7-8)	101 9020	8.90
L (8-9)	101 4477	8.90
from 10 packs per variety from		8.50
30 packs per variety each		7.90



Maxima

New Generation latex gloves

Powder free

Particularly suitable for hypersensitive users who have developed allergies to chemical additives. The gloves fit left and right. Polymer covered, fully textured. Reinforced cuffs protect against rolling up and increase the tear resistance of the glove. The glove has a good fit, offers high sensitivity and excellent tactile sensitivity. Corresponds to DIN EN 455 and 374.

Pack (100 pieces)

Size	Order no.	Price
XS (5-6)	900 9548	8.90
S (6-7)	900 9549	8.90
M (7-8)	900 9550	8.90
L (8-9)	900 9551	8.90
XL (9-10)	900 9552	8.90
from 10 packs per variety from		8.50
30 packs per variety each		7.90



Nitrile examination gloves

Powder free

- > The HS nitrile gloves have a very good material strength, which offers optimal protection
- > They reduce hand fatigue, promote comfort and support excellent tactile sensitivity with soft Comfort
- > Ambidextrous, straight fingers and structured surface
- > Latex-free, powder-free
Pack (100 pieces)

Size	Order no.	Price
Color blue, scent: grape		
xs	988 3755	9.90
s	988 3756	9.90
M	988 3757	9.90
L	988 3758	9.90
XL	988 3759	9.90

Color orange, scent: orange		
XS	988 3750	9.90
S	988 3751	9.90
M	988 3752	9.90
L	988 3753	9.90

Color green, scent: green tea		
XS	988 3760	9.90
S	988 3761	9.90
M	988 3762	9.90
L	988 3763	9.90
XL	988 3764	9.90

Color pink, scent: bubble gum		
XS	988 3765	9.90
S	988 3766	9.90
M	988 3767	9.90
L	988 3768	9.90
XL	988 3769	9.90

from 10 packs per variety from	8.90
30 packs per variety each	8.50



With microstructured
Surface for optimal grip

Premium grip

Latex examination gloves

Powder free

Non-sterile, powder-free gloves with a water-soluble protein content of less than 50 µg/g, high wearing comfort thanks to good fit, thin but tear-resistant, tightness value AQL 1.5, fits left and right, in an optimized dispenser box. White.

Pack (100 pieces)

Size	Order no.	Price
XS	988 3706	9.50
S	988 3707	9.50
M	988 3708	9.50
L	988 3709	9.50
XL	988 3710	9.50
from 10 packs per variety from		8.90
30 packs per variety each		8.50

100 pieces from

9.50

PRICE-PERFORMANCE

With aloe vera
Internal coating



Aloe Vera Latex

Examination gloves

Powder free

Powder-free latex examination gloves with excellent tactile sensitivity. Thanks to the special aloe vera inner coating, your hands are cared for even while you work. Corresponds to DIN EN 455.

Pack (100 pieces)

Size	Order no.	Price
XS (5-6)	900 2871	10.50
S (6-7)	900 2872	10.50
M (7-8)	900 2873	10.50
L (8-9)	900 2874	10.50
XL (9-10)	900 2875	10.50
from 10 packs per variety from		9.90
30 packs per variety each		9.50



Pack of 200 pieces

Nitrile examination gloves

Powder free

- > Soft, thin and tear-resistant with high tactile sensitivity
- > Excellent comfort and fit
- > EN 455 and EN 374 (PPE category III)
- > Powder and latex free,
Pack (200 pieces)

Size	Number	Order no.	Price
XS	200 pieces	988 4725	11.50
S	200 pieces	988 4726	11.50
M	200 pieces	988 4727	11.50
L	200 pieces	988 4728	11.50
XL	180 pieces	988 4729	11.50
from 10 packs per variety from			10.50
30 packs per variety each			9.90



100 pieces from

7.90

PRICE-PERFORMANCE



Nitrile

Examination gloves
Powder free

Non-sterile latex-free examination gloves. Fully structured surface. Powder-free nitrile glove with rolled edge, good fit, high wearing comfort. This glove does not contain any natural latex proteins. Color white. Complies with DIN EN 455, 420 and 374. Pack (100 pieces)

Size	Order no.	Price
XS (5-6)	988 0926	9.50
S (6-7)	988 0927	9.50
M (7-8)	988 0928	9.50
L (8-9)	988 0929	9.50
XL (9-10)	988 0930	9.50
from 10 packs per variety		8.90
from 30 packs per variety each		7.90



NEW

Black nitrile examination gloves

Powder free

Non-sterile, powder-free, black nitrile glove with rolled edge, seamless. Tightness value AQL 1.5, good fit, high wearing comfort, fits right and left. This glove does not contain any natural latex proteins. The glove complies with EN standards 420 and 374-1-2. EPI / PPE Cat. 1. Pack (100 pieces)

Size	Order no.	Price
S	957 833	11.90
M	957 834	11.90
L	957 835	11.90
XL	957 836	11.90
from 10 packs per variety		10.90
from 30 packs per variety each		9.90



CRITERION vinyl
Examination gloves
Powder free

Non-sterile latex-free gloves made of vinyl, powder-free, seal value AQL 1.5, high wearing comfort thanks to good fit, thin but tear-resistant, seamless, in an optimized dispenser box. Corresponds to DIN EN 455, 420 and 374.

Size	Number	Order no.	Price
S (6-7)	100 pieces	900 8166	9.50
M (7-8)	100 pieces	900 8165	9.50
L (8-9)	100 pieces	900 8164	9.50
XL	90 pieces	900 8163	9.50
(9-10) from 10 packs per			8.90
variety from 30 packs per variety each			7.90



Pack off

7.90

PRICE-PERFORMANCE

CRITERION vinyl
Examination gloves
Lightly powdered

Non-sterile latex-free gloves made of vinyl, lightly powdered, seal value AQL 1.5, high wearing comfort thanks to good fit, thin but tear-resistant, seamless, in an optimized dispenser box. Corresponds to DIN EN 455, 420 and 374.

Size	Number	Order no.	Price
S (6-7)	100 pieces	900 8171	9.50
M (7-8)	100 pieces	900 8170	9.50
L (8-9)	100 pieces	900 8169	9.50
XL	90 pieces	900 8168	9.50
(9-10) from 10 packs per			8.90
variety from 30 packs per variety each			7.90



Extra long

Length 400mm



Nitrile examination gloves, extra long

Powder free

Non-sterile, powder-free nitrile glove with extra long shaft (400 mm glove). Pack (50 pieces)

- > Seamless
- > Tightness value: AQL 1.5
- > Complies with EN 455 part 1,2,3 + EN 374-3 (chemical protection) standards
- > High wearing comfort thanks to good fit, fits right and left
- > Very tear-resistant due to higher material thickness

Size	Order no.	Price
S	988 0159	31.50
M	988 0160	31.50
L	988 0161	31.50
XL	988 0162	31.50
from 5 packs per variety		29.90



NEW

Nitrile examination gloves without additives

Powder free

Non-sterile, powder-free nitrile glove with rolled edge. Accelerator-free nitrile gloves - low dermatitis potential, suitable for people who may be hypersensitive to chemical additives. Ambidextrous. Fingertips structured, seamless, tightness value AQL 1.5, good fit, high wearing comfort, fits left and right. This glove does not contain any natural latex proteins. White.

Pack (100 pieces)

Size	Order no.	Price
XS	988 3069	10.50
S	988 3070	10.50
M	988 3071	10.50
L	988 3072	10.50
XL	988 3073	10.50
from 10 packs per variety		9.50
from 30 packs per variety each		8.90



Easy to put on even with wet hands

CRITERION
Surgical gloves

Sterile, powder-free

- > Powder-free, no skin irritation, very good quality
 - > Very comfortable to wear, good tactile feel
 - > Easy to put on, even with wet or damp hands through polymeric inner coating
 - > Micro-roughened surface — very easy to grip
- Corresponds to DIN EN 455, 420 and 374.
Pack (100 pieces)

Size	Order no.	Price
5.5	104 6125	39.90
6	104 5073	39.90
6.5	104 9218	39.90
7	104 3502	39.90
7.5	104 3503	39.90
8	104 3563	39.90
8.5	104 7315	39.90
9	104 7227	39.90
from 4 packs per variety		36.90



Easy to put on thanks to light powdering

CRITERION
Surgical gloves

Sterile, lightly powdered

- > Light powdering for easy putting on of the gloves
 - > Good wearing comfort
 - > Good tactile feeling
- Corresponds to DIN EN 455, 420 and 374.
Pack (100 pieces)

Size	Order no.	Price
5.5	104 7721	34.90
6	104 3500	34.90
6.5	104 9242	34.90
7	104 3501	34.90
7.5	104 5920	34.90
8	104 4696	34.90
8.5	104 5787	34.90
9	104 9192	34.90
from 4 packs per variety		31.90

**Maxima mouthguard**

surgical for binding

Mouthguard to tie. Latex-free and fiberglass-free, 3-layer. For excellent wearing comfort and the best breathability. Microorganism-free product. All of our masks are produced in an ISO 7 cleanroom and are compatible in accordance with the EN ISO 11737-1 standard requirements. Type IIR masks with splash protection, comply with the EN 14683:2014 standard. Bacterial filtration efficiency: up to 99%. Pack of 50 pieces.

Size	Order no.	Price
blue	900 9068	7.90
green	900 2164	7.90
pink	900 2167	7.90
	900 2165	7.90
white from 5 packs per variety		6.90

**Maxima mouthguard**

with elastic ear loops Mouth mask with elastic ear loops

(Earloop). Latex-free and fiberglass-free, 3-layer. For excellent wearing comfort and the best breathability. Microorganism-free product. All of our masks are produced in an ISO 7 cleanroom and are compatible in accordance with the EN ISO 11737-1 standard requirements. Type IIR masks with splash protection, comply with the EN 14683:2014 standard. Bacterial filtration efficiency up to 99%. Pack of 50 pieces

Size	Order no.	Price
blue	900 9061	6.90
pink	900 9063	6.90
	900 9064	6.90
	900 2168	6.90
	900 9062	6.90
yellow green white from 5 packs per variety		5.90

**Skin stapler**

The sterile, disposable skin stapler impresses with its ergonomic design and user-friendly handling. The flattened device head gives the user an ideal view of the wound. The skin stapler contains 35 staples. The remaining number of brackets is easy to read.

Format 5.9 x 3.9 mm pieces from	900 7010	7.50
5 pieces each		6.50

**Disposable surgical gowns**

Sterile disposable surgical gowns. Made of three-layer, liquid-repellent fleece. Overlapping at the back for tying. With elastic cuffs. Colour blue. Pack (25 pieces)

Size	Order no.	Price
M (length: 125 cm, circumference: 150 cm)	988 3773	159.90
L (length: 135 cm, circumference: 155 cm)	988 3774	159.90
XL (length: 145 cm, circumference: 164 cm)	988 3775	159.90
XXL (length: 155 cm, circumference: 164 cm)	988 3776	159.90

**Overshoes**


Non-slip overshoes in blue, one size. The product is non-sterile and latex-free. Pack (100 pieces)

Pack (50 pairs)	988 9864	7.90
from 10 packs each		7.49

**Disposable razors**


> No risk of injury thanks to blade protection, ideal angle of the razor
> Ready to use

Pack (100 pieces)	988 2214	16.90
from 5 packs each		15.50

**Staple remover**

The disposable staple remover can be used in various surgical procedures. It is a high-quality medical product and the ergonomic design makes it easy and quick to use. Packed sterile.

Piece	900 7011	2.59
from 10 pieces each		1.99

**Visitor smock**

Closing at the back, made of non-woven fabric, with tie belt and rubber cuffs, germ and moisture repellent.

Length: approx. 100 cm.

Pack (10 pieces)

green	988 0937	29.90
from 5 packs each		25.90

**Acrylic glove dispenser**

Size	Order no.	Price
Single box	979 5054	18.90
from 3 pieces each		17.90
Double box	979 5055	20.90
from 3 pieces each		19.90
Triple box	979 5056	25.90
from 3 pieces each		24.90



Fig. similar

Surgical cap

Soft quality in beret shape with elastic band, made of green fleece fabric. Protects the hairstyle. Fiberglass free and breathable. Pack (100 pieces)

Color	Order no.	Price
Diameter approx. 53 cm		
white	988 9860	5.90
blue	988 9862	5.90
from 5 packs per variety		5.49
Diameter approx. 61 cm		
white	988 9861	6.90
blue	988 9863	6.90
from 5 packs per variety		6.49



Pack off

72.90

PRICE-PERFORMANCE

Surgical perforated surgical drape without adhesive strips

Blue sterile non-adhesive perforated drape for surgical procedures.

Size	Quantity	Order no.	Price
50x70cm	ø 7 cm pack (90 pieces)	988 2244	79.90
from 5 packs each			77.90
75x90cm	ø 10 cm pack (30 pieces)	988 2245	74.90
from 5 packs each			72.90



Pack off

72.90

PRICE-PERFORMANCE

Surgical drape with adhesive strips

Blue sterile drape with a 5 cm adhesive strip for surgical procedures.

Size	Crowd	Order no.	Price
50x50cm	Pack (110 pieces)	988 2239	72.90
50x70cm	Pack (60 pieces)	988 2237	74.90
75x90cm	Pack (50 pieces)	988 2238	79.90



Surgical surgical perforated drape with adhesive strips

Blue sterile perforated drape. The oval window has an adhesive strip on the outside. Pack (60 pieces)

Size	Order no.	Price
Two parts		
45x75cm	ø 7 cm 988 2246	149.90
from 5 packs each		139.90



200 pieces from

165.00

PRICE-PERFORMANCE

Surgical tape

Blue sterile adhesive strip 10 x 50 cm.

Pack (200 pieces)	988 2247	179.00
from 5 packs each		165.00



Surgical drape without adhesive strips

Blue sterile drape without adhesive strips for surgical procedures.

Size	Crowd	Order no.	Price
40x50cm	Pack (100 pieces)	988 2240	58.90
50x50cm	Pack (90 pieces)	988 2241	59.90
50 x 70 cm	Pack (90 pieces) 75 x	988 2242	63.90
90 cm	Pack (50 pieces)	988 2243	69.90



NEW

EuroSept Xtra Washing Lotion Washing lotion, pH

skin-neutral, for stressed skin and hands.

The pleasantly fine-pored, creamy foam cleanses gently and thoroughly. The lanolin it contains as well as another moisturizing agent and the skin-friendly cleansing surfactants prevent the skin from drying out during washing. The washing lotion is free of dyes and has a pleasantly subtle scent.



NEW

NEW

EuroSept Xtra Washing Lotion Sensitive

Creamy, delicate, pH skin-neutral washing lotion for regular cleaning of skin and hands. Moisturizing ingredients and glycerin prevent the skin from drying out during washing.

Cleans thoroughly and is economical to use. In addition, the mild washing lotion is free of perfumes and dyes and is therefore suitable for use on healthy skin in order to prevent the development of allergies due to the substances mentioned.

Size	Order no.	Price
500ml bottle	988 4925	3.99
from 5 pieces each		3.49
1 liter bottle	988 4926	5.99
from 5 pieces each		5.39
5 liter canister	988 4928	20.90
from 3 pieces each		18.90

Size	Order no.	Price
500ml bottle	988 4920	3.19
from 5 pieces each		2.99
1 liter bottle	988 4921	5.49
from 5 pieces each		4.99
5 liter canister	988 4923	18.90
from 3 pieces each		16.90

HS washing lotion	Wash lotion	Washing lotion sensitive
cleansing,	ÿ	ÿ
nourishing	ÿ	ÿ
moisturizing	ÿ	ÿ
paraben free	ÿ	ÿ
soap-free	ÿ	ÿ
alkali-free	ÿ	ÿ
Color	blue	-
Perfume	ÿ	ÿ
Shelf life after opening	12 months	12 months
Container sizes	500ml/1L/5L	500ml/1L/5L



Piece off

47.90

PRICE-PERFORMANCE

Wall dispenser

For 500 / 1,000 ml dispenser bottles. Dosage quantity easily adjustable from 1.0 - 3.0 ml in 0.5 ml increments. Pump delivers gels and liquids drip-free. Pre-punched drill holes for easy installation.

Size	Order no.	Price
Wall dispenser for 500 ml bottle		
Piece	979 3828	49.90
from 3 pieces each		47.90
Replacement pump for 500ml bottle		
Piece	979 5121	5.99
Wall dispenser for 1 liter bottle		
Piece	979 3843	54.90
from 3 pieces each		52.90
Replacement pump for 1 liter bottle		
Piece	979 5133	5.99



Very productive & very caring

EuroSept Plus Hand cream

The very high quality cream is quickly absorbed without leaving a greasy film. Soothes irritated skin and provides long-lasting moisture. Even small amounts are sufficient. Suitable for all skin types. The hand cream has a very low perfume additive of just 0.5% to prevent allergies.

Tube (100ml)

Piece	900 8925	4.99
from 5 pieces each		4.49



EuroSept Xtra Lotion

> Rich, long-lasting care lotion with high-quality care complexes for smooth, supple skin (particularly suitable for dry skin)

> Has a moisturizing and moisturizing effect
> Protective and skin-soothing high-quality care

Bottle (500ml)

Piece	988 2694	8.90
from 5 pieces each		7.90



NEW

Pumice mouse

Pumice stone for manual cleaning. The pumice mouse removes all stains on the skin such as ink, paint, oil and other strongly sticking substances.

Piece	900 2127	7.50
from 5 pieces each		6.90



Hand washing brush Autoclavable.

Piece	900 8916	2.59
from 10 pieces each		1.99



Piece off

1.99

PRICE-PERFORMANCE

Autoclavable

Application aids for disinfectants

Size	Order no.	Price
① dosing pump for easy removal and decanting of 5 liter canisters.		
Delivery: 20ml per stroke.		
Piece	979 6272	7.50
Delivery: 30ml per stroke.		
Piece	580 01	8.90
② dosing pump for 500 ml bottles	979 6274	3.89
for 1 liter bottles	979 6275	4.89
③ spray heads for 1 liter bottles		
Piece	979 2975	3.99
④ foam spray heads		
Piece	979 3194	4.49
⑤ Canister tap for filling from 5, 6 or 10 liter canisters - with stopcock.		
Piece	979 6276	8.79
⑥ canister keys for easy opening of the canister caps.		
Piece	979 6278	3.89



Gel-like

SafeSept Max®

Hand disinfection gel

Ready-to-use alcoholic disinfectant for hygienic hand disinfection.

> Fragrance-free, moisturizing
> Aldehyde and phenol free

Spectrum of action: Bactericidal, yeasticidal, tuberculocidal, limited virucidal (incl. HBV, HIV, HCV)

Size	Order no.	Price
1 liter bottle	988 2756	8.50
from 5 bottles each		7.99
5 liter canister	988 2757	36.90
from 3 canisters each		35.90



Virucidal

SafeSept Max®

Hand disinfection liquid

Ready-to-use, virucidal alcoholic hand disinfectant for hygienic hand disinfection. Bactericidal, yeasticidal, tuberculocidal, virucidal (incl. HBV, HIV, HCV). Effective against Noro virus within 15 seconds (EN 14476, MNV). Hygienic hand disinfection (EN 1500): 3 ml - 30 seconds. Particularly skin-friendly. Due to its user-friendly combination of active ingredients, SafeSept Max® Hand Disinfectant Liquid is able to achieve a comprehensive spectrum of effects while simultaneously providing skin care and skin protection.

Composition: 100 g contains: 80 g of ethanol (CAS: 64-17-5); 0.1 g alkyl(C12-16)dimethylbenzyl ammonium chloride.

> Virucidal according to EN 14476
> Effective against norovirus (MNV) within 15 seconds
> Contains care ingredients
> Included in the VAH's current disinfectant list

Size	Order no.	Price
150ml bottle	988 4848	3.29
from 5 bottles each		2.99
500ml bottle	988 2752	5.99
from 5 bottles each		5.49
1 liter bottle	988 4851	8.50
from 5 bottles each		7.99
5 liter canister	988 2754	36.90
from 3 canisters each		35.90



NEW

EuroSept Xtra hand disinfection gel

Alcoholic hand disinfectant gel for use on the skin, particularly skin-friendly and particularly practical to use thanks to the gel form. It contains high-quality moisturizing ingredients that prevent the skin from drying out. The product is compatible with all washing lotions, it is aldehyde-free, phenol-free and does not contain any quaternary ammonium compounds.

Spectrum of action Bactericidal, yeasticidal, limited virucidal plus. For hygienic hand disinfection. Included in the VAH's current list of disinfectants. Use biocidal products carefully. BAuA Reg. No.: N-102481

Size	Order no.	Price
250ml bottle	988 4837	3.59
from 5 bottles each		2.99
500ml bottle	988 4838	5.99
from 5 bottles each		5.49
1 liter bottle	988 4839	8.50
from 5 bottles each		7.99
5 liter canister	988 4840	36.90
from 3 canisters each		35.90
Dosing pump for 500ml and 1 liter	979 2558	5.99

NEW

EuroSept Xtra hand disinfection liquid

Virucidal alcoholic hand disinfectant for use on the skin. The EuroSept Xtra hand disinfectant is particularly skin-friendly. Due to its user-friendly combination of active ingredients, it is able to achieve a comprehensive spectrum of effects with simultaneous skin care and skin protection effects. Effective against norovirus within 15 seconds.

(EN 14476, MNV), bactericidal, yeasticidal, tuberculocidal, virucidal (with HBV, HIV, HCV). Included in the VAH's current disinfectant list. Registered as a biocide product with the Federal Institute for Occupational Safety and Health in accordance with the Biocide Act BAuA Reg. No. N-66697.

Size	Order no.	Price
500ml bottle	988 4844	5.99
from 5 bottles each		5.49
1 liter bottle	988 4845	8.50
from 5 bottles each		7.99
5 liter canister	988 4846	36.90
from 3 canisters each		35.90
Dosing pump for 500ml and 1 liter	979 2558	5.99



Quick disinfection with a ready-to-use cloth.

EuroSept® Xtra
Surface disinfectant wipes

application areas
Ready-to-use, low-alcohol wipes made from lint-free fleece material. Available with neutral or lemon scent.

Advantages

- > Broad antimicrobial activity spectrum ready for use
- > Formula with very good Material compatibility
- > Fast disinfection without quaternary Ammonium compounds (QAC)
- > Innovative booster additives

Physical Properties

Conditions	Concentration and exposure time
Bacteria according to EN 13727, EN 1276, EN 13697 (dirty conditions)	Conc. 1 min. all sizes
Candida albicans according to EN 13624, EN 1650, EN 13697 (dirty conditions)	Conc. 1 min. all sizes
TB (tuberculosis) pathogens (M. terrae) and mycobacteria (M. avium) according to EN 14348 (dirty conditions)	Conc. 1 min. all sizes
"Limited virucidal" (enveloped viruses such as HBV, HIV, HCV according to EN 14476 and EN 16777) (dirty conditions)	Conc. 1 min. all sizes
"Limited virucidal plus" (enveloped viruses, e.g. HBV, HIV, HCV) and noroviruses, adenoviruses and rotaviruses according to EN 14476 and EN 16777) (dirty conditions)	Conc. 1 min. all sizes
Adenoviruses (EN 14476, EN 16777) (pure conditions)	Conc. 1 min. all sizes
Polyoma SV 40 (RKI/DVV) (pure conditions)	Conc. 1 min. all sizes
Exposure times for the impregnated cloths: Bacteria and Candida albicans according to EN 16615 and VAH (dirty conditions)	Conc. 5 min. XS, S Conc. 2 min. L Conc. 1 min. M, XL

Ingredients: 100 g liquid contains: 20.0 g 1-propanol; 5.0 g 2-propanol

EuroSept® Xtra surface disinfection wipes XL wipes (30 x 30 cm) 70 pieces Item no. 9884937 One Step (15 x 30 cm) 140 pieces Item no. 9884777
L Flowpack (20 x 30 cm) 80 pieces Item no. 9884949 S Container (20 x 22 cm) 200 pieces Item no. 9884002
M cloths (15 x 30 cm) 140 pieces. Item no. 9884936 S Refill pack (20 x 22 cm) 200 pieces Item no. 9884003
XS container (13 x 20 cm) 120 pieces Item no. 9884197 XS refill pack (13 x 20 cm) 120 pieces Item no. 9884196

EuroSept® Xtra surface disinfectant wipes lemon scent
S container (20 x 22 cm) 200 pieces item no. 9884951 S Refill pack (20 x 22 cm) 200 pieces Item no. 9884950

Supplementary product: Practical dosage container Item no. 9796544

BRAND QUALITY AT A TOP PRICE



EuroSept Xtra

Quick disinfection wipes M/XL

Ready-to-use wipes soaked in an alcoholic disinfectant solution for disinfecting non-invasive medical products and surfaces in areas close to patients. It is not necessary to wipe the surface again when using EuroSept Xtra Wipes.

Size	Crowd	Order no.	Price
15x30cm	140 pieces	988 4936	18.90
from 3 packs each			17.90
30x30cm	70 pieces	988 4937	20.90
from 3 packs each			18.90



EuroSept Xtra

Quick disinfectant wipes S Lemon

Ready-to-use alcoholic wet wipes for prophylactic surface disinfection of alcohol-resistant non-invasive medical devices. Carefully wipe the surface to be disinfected and allow the disinfectant to take effect.

Only for use on alcohol-compatible materials (e.g. not suitable for acrylic glass). Contains fragrances.

Crowd	Order no.	Price
200 pieces refill pack	988 4950	13.90
from 5 packs each		11.90
200 pieces containers	988 4951	14.90
from 12 packs each		12.90



EuroSept Xtra quick disinfection wipes L

Ready-to-use, alcohol-based wet wipes for prophylactic surface disinfection of non-invasive medical devices in the medical sector. Soaked with EuroSept Xtra Surface Disinfection. Stable for 28 days after opening. The opening and expiry date must be documented manually on the sticker or in the space provided before the bag is opened for the first time. Very good material compatibility, even with materials that are otherwise incompatible with alcohol.

> No separate dispenser system required anymore, therefore practical to use

20x30cm	80 pieces	988 4949	8.90
from 3 packs each			7.90



EuroSept Xtra

Quick disinfectant wipes S

Ready-to-use alcoholic wet wipes for prophylactic surface disinfection of alcohol-resistant, non-invasive medical devices. Carefully wipe the surfaces to be disinfected and allow the disinfectant to take effect. Attention: Not suitable for treating alcohol-sensitive surfaces such as:

Acrylic glass. Composition: 100 g of soaking liquid contains 20 g of 1-propanol

Size	Crowd	Order no.	Price
22x20cm	200 wipes tin	988 4002	15.50
from 3 packs each			13.90
22x20cm	200 wipes, refill	988 4003	12.50
from 3 packs each			10.90



NEW

EuroSept Plus
Disinfectant wipes XS

Ready-to-use, alcohol-based wet wipes for prophylactic surface disinfection of alcohol-resistant, non-invasive medical devices. Carefully wipe the surfaces to be disinfected and allow the disinfectant to take effect. Attention: Not suitable for treating alcohol-sensitive surfaces such as acrylic glass.

Size	Crowd	Order no.	Price
14 x 20 cm	120 wipes tin	988 4197	6.90
from 3 packs			6.50
each 14 x 20 cm	120 wipes refill pack	988 4196	7.90
from 3 packs each			7.50

NEW



EuroSept Plus
One step disinfectant wipes

Ready-to-use cloths (fleece) soaked in an alcoholic disinfectant solution for disinfecting non-invasive medical products and surfaces in areas close to patients. In a practical stand-up pouch

Pack (120 wipes)		
Format 17.5 x 36	988 2785	26.90
cm from 3 packs each		24.90



EuroSept Xtra One Step
disinfectant wipes

EuroSept Xtra System Surface Dispenser Wet Wipes Ready-to-use wet wipes for prophylactic surface disinfection of non-invasive medical devices in the medical sector. Very good material compatibility, even with materials that are otherwise alcohol incompatible.

Pack (140 wipes)		
Format 15 x 30 cm	988 4777	24.90
from 3 packs each		21.90

Protecting the planet for a healthier future

EuroSept Xtra Organic
Quick disinfectant wipes

Advantages

- > 100% biodegradable cloths made from vegetable fibers
- > cloths with a higher weight mean a greater reach with less cloths
- > Fast and effective formula*
- > Aldehyde and phenol free
- > *Tested and found to be effective against a broad spectrum of viruses and Bacteria found – including MRSA, TB, HIV-1 and H1N1

Physical Properties

Conditions

	exposure time
Bactericidal, levurocidal (4-field test) according to EN 16615 (low and high organic load)	60 seconds
Bactericidal, yeasticidal according to EN 13727, EN 13697, EN 13624	30 seconds
Bactericidal, levurocidal (including MRSA) according to EN 1276 EN 1650 EN13697 (low organic load)	60 seconds
Tuberculocide according to EN 14348 (low organic contamination)	60 seconds
Limited virucidal according to EN 14476	30 seconds
Limited virucidal plus, virucidal according to EN 14476 (low organic load)	5 minutes

Ingredients: 100 g liquid contains: 20.0 g 1-propanol; 5.0 g 2-propanol

Further information: Shelf life: 2 years in original packaging; 28 days after first opening, if it stays tightly closed. MD: Class IIa



EuroSept Xtra Bio rapid disinfection wipes

Made from 100% biodegradable vegetable fibers. Aldehyde-free, phenol-free. Pack of 80 pieces

Piece	988 9964	10.90
from 3 packs each		9.90

Instructions for use

Ready-to-use surface disinfectant wipes. Open the safety film under the cap and remove the first wipe from the pack. Use the removed cloth as quickly as possible and do not use it a second time. Completely wet the surfaces with the pre-soaked wipes and spread the disinfectant solution. Keep the surface moist throughout the entire exposure time. In food areas, surfaces must be treated with water of drinking quality. After disinfection, wipe the surfaces with a hygienically clean cloth if necessary. Intended for professional use only.

Keep the package lid and safety film tightly closed.

Highly effective wipes for alcohol-sensitive surfaces.

EuroSept® Xtra wipes for alcohol-free surface disinfection

Advantages

- > Suitable for use on alcohol-sensitive surfaces
- > Broad spectrum of action against Bacteria, yeast and viruses
- > Ready-to-use, pH-neutral formula
Dries practically without leaving any residue

Physical Properties

Conditions	Concentration and exposure time
Bacteria and Candida albicans, according to EN 13727, EN 13624, EN 1276, EN 1650, EN 13697 (dirty conditions)	Conc. 2 min. S
Bacteria and Candida albicans, according to EN 13727, EN 13624, EN 1276, EN 1650 (dirty conditions)	Conc. 1 min. M, XL
Bacteria and Candida albicans according to EN 13697 (dirty conditions)	Conc. 2 min. M, XL
Bacteria and Candida albicans according to EN 13727, EN 13624 (dirty conditions)	Conc. 1 min. L
"Enveloped viruses (e.g. HBV, HIV, HCV) according to RKI/DVV and EN 14476 (dirty conditions)	Conc. 1 min. S, M, XL
Rotaviruses according to EN 14476 (clean conditions)	Conc. 30 min. All sizes
Polyoma SV 40 viruses (RKI/DVV) (dirty conditions)	Conc. 1 min. S, M, L, XL
Exposure times of the impregnated cloths: Bacteria and Candida albicans according to EN 16615 and VAH (dirty conditions) Exposure times of the impregnated cloths: Bacteria and Candida albicans according to EN 16615 and VAH (dirty conditions)	Conc. 1 min. S, M, L, XL

Ingredients: 100 g contains: 0.24 g benzyl-C12-16-alkyldimethylammonium chloride, 0.24 g N, N-didecyl-N, N-dimethylammonium chloride, 0.24 g C12-14 alkyl [(ethylphenyl)methyl] dimethylammonium chloride

EuroSept® Xtra wipes for alcohol-free surface disinfection

S Refill pack (20 x 22 cm) 200 pieces Item no. (9884362) M cloths (15 x 30 cm) 140 pieces Item no. (9884942)

S containers (20 x 22 cm) 200 pieces Item no. (9884363) XL towels (30 x 30 cm) 70 pieces

L flow pack (20 x 30 cm) 80 pieces Item no. (9884948)

Article no. (9884943)

BRAND QUALITY AT A TOP PRICE



NEW

EuroSept Xtra quick disinfection wipes M/XL, without alcohol

Alcohol-free, ready-to-use disinfectant wipes for quick disinfection and cleaning of medical inventory and surfaces of all kinds, including in the medical sector. Based on quaternary ammonium compounds. EuroSept Xtra System Wet Wipes non-alcoholic are alcohol, aldehyde and phenol-free. High, lint-free material quality, guaranteed release of active ingredients and good cleaning power. Included in the VAH's current list of disinfectants. Registered as a biocide product with the Federal Institute for Occupational Safety and Health in accordance with the Biocide Act BAuA Reg. No. N-66695 / N-68687.

Complies with EU Directive 93/42/EEC for medical devices. Neutral smell. Note: Suitable for acrylic glass! For use with empty tissue dispenser 979 6544.

Crowd	Size	Order no.	Price
140 wipes refill bags	15x30cm	988 4942	15.90
from 3 packs each			13.90
70 wipes refill bags	30x30cm	988 4943	17.50
from 3 packs each			15.50

NEW

EuroSept Xtra

Quick disinfection wipes L, without alcohol

Alcohol-free wet wipes for prophylactic surface disinfection of non-invasive medical devices. Particularly suitable for use on alcohol-sensitive surfaces such as acrylic glass. Spectrum of action and exposure times: surface disinfection (wipe disinfection), bactericidal, yeasticidal according to VAH/EN 16615: 5 min.; tuberculocide (M.terrae): 1 min; limited virucidal (including HBV, HIV, HCV), rota viruses: 30 sec.; Noro viruses (MNV):1 min. Composition: 100 g of soaking liquid contains: 45 g of ethanol (CAS: 64-17-5). Contains fragrances.

- > No separate dispenser system required, making it practical to use
- > Alcohol-free wipes for quick disinfection of large areas
- > Quickly effective, particularly gentle on materials, also suitable for alcohol-sensitive surfaces
- > Simple, quick "to the point" application
- > Suitable for all areas close to patients
- > Included in the VAH's current disinfectant list
- > Stable for 28 days after opening

Pack (80 pieces)

20x30cm	80 pieces	988 4948	8.50
from 3 packs each			7.50

NEW

EuroSept Xtra

Quick disinfection wipes S, without alcohol

Ready-to-use alcohol-free wet wipes for prophylactic surface disinfection of non-invasive medical devices in the medical sector. Particularly suitable for use on alcohol-sensitive surfaces such as acrylic glass.

Size	Crowd	Order no.	Price
Tin 20 x 22 cm	200 wipes	988 4363	15.90
from 3 packs each			13.90
Refill pack 20 x 22 cm 200 wipes		988 4362	14.90
from 3 packs each			12.90

Practical lint-free cloths for soaking with Henry Schein disinfectant solutions.

Improved disinfection for alcohol-sensitive surfaces.



Physical Properties

Weight	Size	Sheets per roll	Rolls per box	Liquid per roll
50g/m2	15x30cm	140	4	1.4L
50g/m2	30x30cm	70	4	1.4L

Composition: PES 90% / VIS 10%



EuroSept Xtra cloth dispenser (empty)

Tightly closing dispenser bucket system for hygienic filling of fleece wipes and disinfectants/cleaners or before soaked wipes. The wipes can be removed quickly and conveniently from the practical, reusable bucket.

> Lid diameter: 22.5cm, base diameter: 19.5cm	
> Height including lid: 23cm, height without lid: 21.5cm	
Piece	979 6544
from 3 pieces each	9.50

EuroSept Xtra cleaning cloths (non-woven, dry)

Dry fleece cloths for soaking with VAH/DGHM-listed disinfectants. Application: Place the fleece cloth roll in the dispenser bucket and slowly pour 1.4 liters of ready-to-use disinfectant solution over it. Let the soaked roll stand for about 30 minutes until all cloths are completely evenly moistened. Then write the label with the relevant information. The soaked cloths may be used for a maximum of 28 days! Recommendation according to EN 16615: HS EuroSept Xtra surface disinfection; The exposure times of the respective soaking liquid compatible with the corresponding cloth (according to EN 16615) must be observed.

Size	Crowd	Order no.	Price
15x30cm	3 rolls x 140 wipes	988 4930	33.90
from 3 packs each			29.90
30x30cm	3 rolls x 70 wipes	988 4931	35.90
from 3 packs each			33.90

Physical Properties

Conditions

	concentration and Exposure time*
Bacteria according to EN 13727, EN 1276 and EN 13697 (dirty conditions)	Conc. 1 min.
Candida albicans according to EN 13624, EN 1650 and EN 13697 (dirty conditions)	Conc. 1 min.
Enveloped viruses (e.g. HBV, HIV and HCV) according to RKI/DVV (dirty conditions)	Conc. 1 min.
Rotaviruses according to DVV (low burden)	Conc. 1 min.
Noroviruses according to EN 14476 (low burden)	Conc. 30 min.
Polyoma SV 40 viruses (RKI/DVV) (impurity)	Conc. 30 min.

Ingredients: 100 g contains: 0.24 g benzyl-C12-16-alkyldimethylammonium chloride, 0.24 g N,N-didecyl-N, N-dimethylammonium chloride, 0.24 g C12-14 alkyl[(ethylphenyl)methyl]dimethylammonium chloride

*Do not wipe dry before the end of the application time.



EuroSept Xtra Surface disinfection, without alcohol

Alcohol-free, ready-to-use wipe disinfectant for prophylactic surface disinfection of non-invasive medical products. Particularly suitable for use on alcohol-sensitive surfaces such as acrylic glass. Application: To avoid the formation of streaks, the treated surfaces should be wiped at regular intervals with a damp (water) disposable cloth.

Size	Order no.	Price
1 liter bottle with spray head	988 4905	8.50
from 3 bottles each		6.90
5 liter canister	988 4907	25.90
from 3 canisters each		19.90

Fast-acting disinfection

with a broad spectrum of effects.

Fast-acting, fully virucidal solution

for sensitive surfaces.



EuroSept® Xtra
Surface disinfection

- Advantages
- > Broad spectrum of antimicrobial activity
 - > Fast-acting, ready-to-use solution
 - > Alcohol-based formula evaporates completely for a streak-free surface*
 - > Without quaternary ammonium compounds (QAC)
 - > Available in neutral scent or lemon scent



EuroSept® Xtra alcohol-free, fully virucidal
Solution for surface disinfection

- Advantages
- > Suitable for use on alcohol-sensitive surfaces
 - > Wide spectrum of effects
 - > Fast-acting, ready-to-use formula
 - > Fully virucidal
 - > Suitable for terminal disinfection when thermal sterilization is not possible

Physical Properties

Conditions	Concentration and exposure time
Bacteria according to EN 13727, EN 1276, EN 13697, EN 16615 (high load)	Conc. 1 min.
Candida albicans according to EN 13624, EN 1650, EN 13697, EN 16615 (high load)	Conc. 1 min.
TbB (tuberculosis) pathogens (M. terrae) and mycobacteria (M. avium) according to EN 14348 (high load)	Conc. 1 min.
Enveloped viruses (e.g. HBV, HIV, HCV) according to EN 14476 and EN 16777 (high burden)	Conc. 1 min.
Noroviruses (EN 14476) (clean conditions)	Conc. 1 min.
(EN 16777) (low load)	Conc. 2 min.
Adenoviruses (EN 14476, EN 16777) (low load)	Conc. 1 min.
Polyoma SV 40 (RKI/DVV) (low load)	Conc. 1 min.
“Limited virucidal plus” (noroviruses, adenoviruses and rotaviruses) (EN 14476) (low burden)	Conc. 1 min.
(EN 16777) (low load)	Conc. 2 min.

Ingredients: 100 g contains: 20.0 g 1-propanol; 5.0 g 2-propanol



EuroSept Xtra
Surface disinfection

Ready-to-use rapid disinfection for prophylactic surface disinfection of non-invasive medical devices in the medical sector. Very good material compatibility even with materials that are otherwise incompatible with alcohol (due to the low alcohol content the material can be used on many surfaces).

Size	Order no.	Price
1 liter bottle with spray head	988 4881	8.50
from 3 bottles each		6.90
5 liter canister	988 4883	25.90
from 3 canisters each		19.90



EuroSept Xtra surface
disinfection Lemon

Ready-to-use rapid disinfection for prophylactic surface disinfection of non-invasive medical devices in the medical sector. Very good material compatibility even with materials that are otherwise incompatible with alcohol (due to the low alcohol content the material can be used on many surfaces).

Size	Order no.	Price
1 liter bottle with spray head	988 4885	8.50
from 3 bottles each		6.90
5 liter canister	988 4887	25.90
from 3 canisters each		19.90

Physical Properties

Conditions	Concentration and exposure time
Bacteria according to EN 13727 and EN 13697 (high conditions)	Conc. 1 min.
Candida albicans according to EN 13624 and EN 13697 (high conditions)	Conc. 1 min.
TB pathogen according to EN 14348 (M. terrae) (low burden)	Conc. 5 min.
Virucide according to EN 14476 (low exposure)	Conc. 2 min.
Spores (C. difficile) according to EN 13704 (low load)	Conc. 60 min.

Ingredients: 100 g contains: 0.48 g didecyldimethylammonium chloride, 0.2 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine

*Do not wipe dry before the end of the application time.



EuroSept Xtra
Surface disinfection, without alcohol

Virucidal, ready-to-use alcohol-free wipe disinfectant for mechanical surface disinfection of medical devices. Also suitable for terminal disinfection when thermal sterilization is not possible. Bacteria according to EN 13727, VAH and EN 16615 (high load): conc. 2 min. Candida albicans according to EN 13727, VAH and EN 16615 (high load): conc. 2 min. 100 g of soaking liquid contains: 0.48 g didecyldimethylammonium chloride; 0.2 g N-(3-Aminopropyl)-N-dodecyl-propane-1,3-diamine. VO 648/2004/EC: Contains: < 5% cationic surfactants, fragrances.

Size	Order no.	Price
1 liter bottle with spray head	988 4889	9.90
from 3 bottles each		7.90
5 liter canister	988 4891	28.90
from 3 canisters each		24.90



Effective disinfection for

sensitive surfaces.

EuroSept® Xtra alcohol-free
Concentrate for surface disinfection

- Advantages
- > Suitable for use on alcohol-sensitive surfaces
 - > Broad spectrum of effectiveness without sticky residue
 - > Suitable for dosing systems
 - > Aldehyde and phenol free
 - > Concentrated formula for less packaging and more sustainability

Physical Properties

Conditions	concentration and Exposure time*
Bacteria and Candida albicans (EN 13727, EN 13697, EN 13624, VAH, EN 16615) (high load)	0.5% 15 min.
TB pathogen (M. terrae) according to EN 14348 (low burden)	2.0% 15 mins
(low load)	0.5% 60 mins
(high load)	1.0% 60 mins
“Limited virucidal” (enveloped viruses according to RKI/DVV, e.g. HBV, HCV, HIV) (high conditions)	0.5% 15 mins
Rotavirus (EN 14476) (low strain)	0.5% 15 mins
Norovirus (EN 14476) (high load)	1.5% 15 mins 1.0% 60 mins
Adenovirus (EN 14476) (low strain)	2.0% 30 mins
“Limited virucidal Plus” (EN 14476) (low exposure)	2.0% 30 mins
Polyoma SV 40 viruses (EN 14476) (high load)	1.0% 15 mins

Ingredients: 100 g contains: 10.8 g didecyltrimethylammonium chloride, 6.0 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine

*Do not wipe dry before the end of the application time.

NEW



EuroSept Xtra
Surface disinfection
Concentrate,
without alcohol

Surface disinfectant for prophylactic surface disinfection of medical devices. To avoid the formation of streaks, the treated surfaces should be wiped at regular intervals with a damp (water) disposable cloth. Only to be used as a wipe disinfectant! Unused, ready-to-use solution must be changed daily.

Size	Order no.	Price
5 liter canister	988 4903	61.90
from 3 canisters each		58.90



Maxima surface disinfection is virucidal, without alcohol

Ready-to-use, fast-acting surface disinfectant for cleaning and disinfecting non-invasive medical devices. Well suited for alcohol-sensitive surfaces such as acrylic glass. The foam has a bactericidal (with MRSA), fungicidal and limited virucidal effect (with HIV, HCV, HBV, H1N1 and H5N1) according to VAH/DGHM. Also effective against enveloped viruses, rota and noro viruses. VAH/DGHM listed.

750ml bottle with foam spray head	900 6137	9.90	Piece	988 2739	7.49
-----------------------------------	----------	------	-------	----------	------



SafeSept Max®
Quick disinfectant wipes, without alcohol

- > Quickly effective, particularly gentle on materials
- > Practical to use
- > Alcohol-free wipes for quick disinfection of large areas Surface disinfection (wiping disinfection): Bactericidal, yeasticidal according to VAH: 1 min.; effective against all enveloped viruses. RKI recommendation including HBV/HIV/HCV/Vaccinia/BVDV/Influenza/Flu viruses: 30 seconds. Dimensions: 14 x 20 cm

Dispenser container (150 wipes)

Maxima
Surface disinfection, without alcohol

Pre-soaked disinfectant wipes in a practical dispenser bucket. Aqueous wet wipes for prophylactic surface disinfection of non-invasive medical devices. Exposure times of the soaking liquid: Bacteria and Candida albicans according to EN 13727, EN 1276, EN 1650, EN 13624 and EN 13697 (high load): Conc. 2 min., Limited virucidal (enveloped viruses e.g. HBV, HIV, HCV) according to RKI/DVV with high load: Conc. 1 min., Rota viruses according to EN 14476 (low load): Conc. 1 min., Noro viruses according to EN 14476 (low load): Conc. 30 min., Polyoma-SV 40 viruses (RKI/DVV with exposure): Conc. 30 min. Exposure times for the soaked cloths: Bacteria and Candida albicans according to EN 16615 and VAH under high exposure: Conc. 1 min. Soaked with Note: Suitable for acrylic glass! Dimensions: 28 x 28 cm

Crowd	Order no.	Price
Dispenser bucket with 75 pre-soaked wipes	988 4345	25.90
Refill pack with 75 pre-soaked wipes	988 4346	22.90





NEW



EuroSept Xtra instrument disinfection concentrate

Liquid concentrate for the manual disinfectant cleaning of medical instruments using the immersion bath process. Spectrum of action: bactericidal, yeasticidal, limited virucidal. Also suitable for use in an ultrasonic bath at room temperature. Very good material compatibility, extremely gentle to use (aldehyde and phenol free). Does not fix proteins (blood) and has very good cleaning performance. Not suitable for final disinfection of semi-critical and critical medical devices.

Size	Order no.	Price
2 liter bottle	988 4863	26.50
5 liter canister	988 4864	51.90
from 3 canisters each		49.90

NEW



EuroSept Xtra Enzyme cleaner for instruments

Offers good material compatibility and is therefore the ideal cleaning concentrate for the manual reprocessing of surgical instruments. Does not foam and is suitable for use in an ultrasonic bath. The enzyme cleaner offers excellent cleaning performance by combining non-ionic surfactants with a stabilized complex of three enzymes. The self-acting formula effectively penetrates and removes all types of organic contaminants such as clotted blood and mucus, even in hard-to-reach areas. Suitable for ultrasonic bath. Instrument cleaning: 1% (10 ml/L). Cleaning in an ultrasonic bath: 0.5% (5 ml/L). Temperature: Prepare solution with warm water (30 - 40 °C). Exposure time: at least 5 minutes. Complies with EU Directive 93/42/EEC for medical devices. Neutral smell.

Size	Order no.	Price
1 liter bottle	988 4858	16.90
5 liter canister	988 4236	64.90
from 3 canisters each		59.90

NEW



EuroSept Xtra Instrument disinfection concentrate virucidal

Liquid concentrate for manual, regular preparation and terminal disinfection of non-critical and semi-critical medical and surgical instruments using the immersion bath process. Not suitable for the terminal disinfection of critical medical devices!

5 liter canister	988 4868	69.50
from 3 canisters each		67.90

EuroSept Xtra Tray cleaners

Liquid concentrate for cleaning as part of the correct preparation of medical devices such as impression trays and mixing aids. The innovative formulation offers very good cleaning performance due to the mildly alkaline recipe and is gentle on materials at the same time.

EuroSept Plus Tray Cleaner Concentrate is suitable for cleaning all corrosion-resistant metals as well as sensitive materials such as anodized aluminum. Also ideal for use in ultrasonic baths. Complies with EU Directive 93/42/EEC for medical devices. Neutral smell. Bottle (1 liter)

1 liter bottle	988 4916	25.90
from 3 bottles each		22.90

Instrument tub

With sieve insert and dosage marking.
Dimensions: 34 x 21 x 12 cm (W x H x L)

4 liter tub	979 6282	54.90
-------------	----------	-------



NEW



EuroSept Plus Orange Oil Cleaner

Special cleaning solution for practices and laboratories with natural orange oil to remove zinc oxide cement, alginate and plaster residues as well as impression pastes, blood residues and other contaminants.

The Orange Oil Cleaner cleans impression trays, instruments, glass plates and ceramics. It is therefore ideal for use on stainless steel and other surfaces. Please check for compatibility in an inconspicuous area beforehand. Pleasant orange scent, high concentration of orange oil. Attention: very high cleaning effect!

Compared to the HS-Orange Solvent, the HS-Orange Oil-Cleaner is made from pure, high-quality orange oil and therefore has a higher cleaning effect. Not suitable for painted surfaces. Do not use on plastic surfaces.

250ml bottle	988 2716	8.90
from 5 bottles each		7.90



EuroSept Plus Foam cleaner

Ready-to-use cleaning spray for quick and intensive cleaning of plastic surfaces, artificial leather surfaces, devices and objects.

- > Very compatible with materials, even with sensitive ones
- Surfaces such as artificial leather or plexiglass
- > Excellent cleaning performance
- > Ready to use, without dyes

The ready-to-use solution is applied undiluted to objects and surfaces from a distance of approx. 20-30 cm as foam (use a suitable foam spray head) and distributed evenly with a lint-free, undyed disposable cloth and allowed to take effect briefly. Make sure it is completely wetted.

Before large-area applications on sensitive materials such as artificial leather, test for material compatibility in a small, hidden area. To avoid the formation of streaks, the treated surfaces should be wiped at regular intervals with a disposable cloth moistened with water. VO 648/2004/EC: Contains < 5% cationic surfactants, perfume. Spray bottle (1 liter)

1 liter bottle with spray head	988 3996	11.90
from 3 bottles each		9.90

NEW



EuroSept Plus Orange solvent

Special cleaning solution for practices and laboratories with natural orange oil to remove zinc oxide cement, alginate and plaster residues as well as impression pastes, blood residues and other contaminants. The Orange Solvent cleans impression trays, instruments, glass plates and ceramics. Ideal for use on stainless steel and other surfaces. Please check for compatibility in an inconspicuous area beforehand. Pleasant orange scent.

Attention: Not suitable for painted surfaces. Do not use on plastic surfaces. Very high cleaning effect!

Size	Order no.	Price
250ml bottle	988 2719	8.90
from 3 bottles		7.90
each 500 ml bottle	988 2718	15.90
from 3 bottles each		14.90

400ml can only

9.90

PRICE-PERFORMANCE



Silicone spray For the care of medical instruments.

Increases the gliding properties of catheters and probes and prevents them from sticking. Spray can (400 ml)

Piece	900 2402	10.90
from 5 cans each		9.90

NEW



Maxima® oil spray

Universal cleaning and care spray for handpieces and contra-angle handpieces, air- and ball-bearing turbines. CFC-free. Without spray attachment. Spray can (500 ml)

Pieces (without adapter)	900 1570	8.50
from 6 pieces each		6.90



EuroSept Plus liquid
cleaner for thermal disinfectant Foam-free, alkaline cleaning concentrate for mechanical cleaning of surgical instruments, medical utensils, as well as laboratory glass in hospitals, doctor's offices and laboratories. Anodized aluminum surfaces and titanium alloys should be checked for durability as color changes may occur. This must be taken into account, for example, for implants with color coding. Dosing must be carried out using a suitable dosing device and depends on the degree of contamination of the product. It is approx. 2 ml - 5 ml. Labeling according to Regulation (EC) No.

648/2004: Contains < 5 - 15% amphoteric surfactants, enzymes.

5 liter canister	988 3990	41.90
from 3 canisters each		39.90



EuroSept Xtra powder
cleaner for thermal disinfectant
Cleaner (granules) for the mechanical preparation of dental instruments, with the exception of handpieces and contra-angle handpieces. After the cleaning step, neutralization with Thermo Neutralizer is recommended. Rinse aid is required to dry the instruments stain-free. Instruments made of anodized aluminum must be checked for alkali resistance. Labeling according to VO (EC) 648/2004: 30% phosphates, 5% non-ionic surfactants. Drying in the last rinse cycle. Neutral smell.
Set: 5 liter cleanser , 5 l neutralizer, 1 l rinse aid.

Buckets (4.5 kg)	988 4364	40.90
from 3 buckets each		38.90



EuroSept Xtra liquid
cleaner for thermal disinfectant
Mildly alkaline liquid cleaner for the reprocessing of surgical, microsurgical and MIS (minimally invasive surgery) instruments, rigid endoscopes, anesthesia accessories, containers. Cleaning with a gentle pH value. Can be used both for mechanical cleaning of instruments (RDG) and for cleaning endoscopes in the RDG-E. Removes residues of dried and denatured blood. Supports the removal of biofilms. No neutralization step required for mechanical processing, resulting in short program sequences.

5 liter canister	988 9743	41.90
from 3 canisters each		39.90



EuroSept Plus neutralizer for thermal disinfectant

Neutralizer based on citric acid for the mechanical preparation of dental instruments after alkaline cleaning to optimize the cleaning result. Neutral smell.

5 liter canister	988 3989	40.90
from 3 canisters each		38.90

EuroSept Plus rinse aid for thermosinfectant

Low-foaming rinse aid (concentrate) for the mechanical reprocessing of medical instruments and utensils in automatic cleaning and disinfection machines (RDG). The rinse aid is automatically added in the final rinse using a dosing device. The dosage depends, among other things, on the items to be washed and the water quality and is 0.3 ml – 1.0 ml/l.

1 liter bottle	988 3987	10.90
from 3 bottles each		9.90



Towel rolls
Soft, bright white pulp wipes. Internal roll, 20 cm wide.

Size	Order no.	Price
Cardboard (6 rolls of 280m each) 1-layer	900 3315	59.90
Cardboard (6 rolls of 160m each) 2-layer	900 3316	59.90
Cardboard (12 rolls of 120m each) 1-layer	900 3314	56.90
Cardboard (12 rolls of 60m each) 2-layer	988 4366	56.90



Disposable towels

1-ply disposable towels, recycled quality.

Recycled white. ZZ fold.	Dimensions: 25 x 23 cm	
Box (5,000 pieces)	900 3324	53.90

2-layer disposable towels, particularly soft.

White. Z-fold.	Dimensions: 21 x 21.5 cm	
Box (1,000 pieces)	900 9964	57.90

White. ZZ fold.	Dimensions: 21 x 24 cm	
Box (3,200 pieces)	900 3325	56.90

White. C fold.	Dimensions: 25 x 31cm	
Carton (2,432 pieces)	900 3323	52.90

Recycled white. C-fold box (1,920	Recycled quality. Dimensions: 25 x 41cm	900 6451	56.90
-----------------------------------	---	----------	-------

pieces) 3-layer disposable towels, particularly soft and absorbent.

Pure white. ZZ fold.	Dimensions: 22 x 32 cm	
Box (2,500 pieces)	988 0905	74.90

Pure white. ZZ fold.	Dimensions: 22 x 42 cm	
Box (2,000 pieces)	900 3327	76.90



Suitable for papers with C and ZZ folds



Practical dispensers for towel rolls

Paper towel dispenser
White towel dispensers.

Size	Order no.	Price
1 large paper towel dispenser (for C and Z folds, suitable for order no. 900 3324 / 900 3325 / 900 3326 / 900 3327) Dimensions: 38 x 28 x 14 cm		
Piece	900 8446	31.90
from 3 pieces each		30.90
2 small paper towel dispensers (for C and ZZ folds, suitable for the Be position number 900 3320 / 900 3321/ 900 3323 / 900 6451) Dimensions: 38 x 28 x 11cm		
Piece	900 8447	29.90
from 3 pieces each		28.90
3 large paper towel dispensers (suitable for towel rolls 900 3315 / 900 3316) Dimensions: 37 x 25 x 23 cm		
Piece	900 8448	32.90
from 3 pieces each		31.90
4 small paper towel dispensers (suitable for towel rolls 900 3314 / 900 3318) Dimensions: 35 x 18 x 17 cm		
Piece	900 8449	25.90
from 3 pieces each		24.90
Spare key for small paper towel dispenser (not shown)		
Piece	988 0657	3.50



30 pieces from

6.90

PRICE-PERFORMANCE

Medical records

Absorbent cellulose flakes with soft textile fleece seal and underside made of liquid-tight PE film.

Pack (30 pieces)

Size	Quantity	Order no.	Price
40x60cm	22g	900 8132	7.90
from 9 packs each			6.90
60x60cm	35g	900 8133	8.90
from 4 packs each			7.90
60x90cm	55g	900 8134	10.90
from 4 packs each			9.90



Pack off

37.50

PRICE-PERFORMANCE

cellulose

Unbleached

Size	Order no.	Price
20x30cm	3 x 5kg 900 8866	38.90
40x60cm	3 x 5kg 900 8865	38.90
from 3 pieces per variety		37.50



Carton off

37.50

PRICE-PERFORMANCE

New quality: extra soft

Doctor roles

2-ply. White. Particularly tear-resistant. Very absorbent. Perforations every 40cm. Length: 50 m.

Carton (9 rolls)

Size	Order no.	Price
39cm wide	900 8496	37.50
45cm wide	900 8497	39.90
50cm wide	900 8498	45.90
55cm wide	900 8499	49.90
59cm wide	900 8500	54.90

Protecting the planet for a healthier future



Medical rolls made from

recycled paper

Very durable. Perforations every 40cm

Carton (9 rolls)

Size	Order no.	Price
Roll width 50 cm	Total length 50m 900 9962	47.90
Roll width 59 cm	Total length 50m 900 9963	56.50

100 wipes

19

1,

PRICE-PERFORMANCE

Without color and
Fragrances

Tissue towels

Particularly soft, white, 2-ply, low-lint tissues. The tissue cloths can be used in a variety of ways. Without dyes and fragrances.

Dispenser box (100 wipes)	900 2377	1.29
from 30 pieces each		1.19



Waste bags

Size	Order no.	Price
1 waste bag made of low-pressure film. Transparent, rolled, low-pressure film (wrinkles). Pack (20 x 20 pieces)		
30 liters	900 9403	37.90
2 Solid quality		
Color: black, rolled. Pack (20 x 20 pieces)		
30 liters	900 9400	33.90
50 liters	900 9401	59.90
120 liters	900 9402	59.90

NEW

Pack only

36.90

PRICE-PERFORMANCE



Wet crepe napkins

> Absorbent and wet-resistant

> Size 34 x 38 cm

> Color: White

Pack (1,000 pieces)	900 1391	36.90
---------------------	----------	-------

50 pieces from

9.90

PRICE-PERFORMANCE

Particularly absorbent
and soft

Multi-purpose cloths


Ideal for cleaning instruments, laboratory tables and devices. Food safe. Dimensions: 32 x 36 cm.

> 8 times absorbency

> Sterilizable up to 134 °C

Box (50 wipes)	900 2938	10.90
from 10 pieces each		9.90

NEW



cellulose napkins

Soft and absorbent cellulose napkins.

Size	Crowd	Order no.	Price
40 x 40 cm 2-ply	125 pieces	900 3753	5.50
from 12 packs			4.90
each 40 x 40 cm	12 x 125 pieces	900 3785	49.90
2-ply 40 x 40 cm	50 pieces	900 3754	4.90
4-ply from 16			4.59
packs each 40 x 40 cm 4-ply 50 pieces		900 3786	51.90

Premium toilet paper

Bright white, embossed tissue.

Size	Order no.	Price
2 ply	Pack (48 rolls of 200 sheets) 900 3328	19.90
3 ply	Pack (56 rolls of 250 sheets) 900 6481	36.90



Pack off


19.90

PRICE-PERFORMANCE



Scalpel handle No. 3
Suitable for blades No. 10-15. Length: 12.5 cm, rustproof.

Made of metal	Piece	100 7520	5.99
---------------	-------	----------	------



Carbon scalpel blades
Sterile scalpel blades made of carbon-reinforced steel.

figure	Crowd	Order no.	Price
10	10 pieces	112 6183	5.99
11	10 pieces	112 6175	5.99
12	10 pieces	112 6176	5.99
12D/12B	10 pieces	112 6180	5.99
15	10 pieces	112 6181	5.99
15C	10 pieces	112 6186	5.99
18	10 pieces	572 0864	5.99
20	10 pieces	112 6182	5.99
from 5 packs each			4.99
10	100 pieces	572 0852	11.90
11	100 pieces	572 0853	11.90
12	100 pieces	572 0854	11.90
12D/12B	100 pieces	572 0855	11.90
15	100 pieces	572 0856	11.90
15C	100 pieces	572 0857	11.90
18	100 pieces	572 0858	11.90
20	100 pieces	572 0859	11.90
21	100 pieces	572 0860	11.90
22	100 pieces	572 0861	11.90
23	100 pieces	572 0862	11.90
from 5 packs each			9.90



Scalpel handle No. 4
Suitable for blades No. 20-24. Length: 12.5 cm, rustproof.

Made of metal	Piece	100 3973	5.99
---------------	-------	----------	------

figure



Disposable scalpels

Sterile disposable scalpels made of stainless steel with plastic handle. Pack (10 pieces)

figure	Order no.	Price
10	112 6196	4.50
11	112 6189	4.50
12	112 6194	4.50
12B	112 6188	4.50
15	112 6187	4.50
15C	572 0866	4.50
18	572 0867	4.50
20	112 6193	4.50
21	572 0868	4.50
22	572 0869	4.50
23	572 0870	4.50
24	572 0871	4.50
from 5 packs each		3.75

piece only

5.99

PRICE-PERFORMANCE



Scalpel handle Air
"Air instruments" made of stainless steel, extremely light and ergonomically shaped. Perfect grip. Easy to clean. Figure 3 is suitable for blade shapes 10, 11, 12, 12B, 15 and 15C. Figure 4 fits blades 20, 21, 22 and 23.
Length 12.5cm.

Size	Order no.	Price
Figure 3	900 8098	5.90
Figure 4	900 8099	5.90



default

Size	Order no.	Price
Pointed/pointed, curved		
115mm	900 3101	11.90
130mm	900 3100	12.90
Pointed/pointed, straight		
115mm	900 3109	11.90
130mm	900 3108	12.90



default

Size	Order no.	Price
Pointed/blunt, curved		
115mm	900 3105	11.90
130mm	900 3104	12.90
Pointed/blunt, straight		
115mm	900 3113	9.90
130mm	900 3112	12.90
165mm	900 3110	14.90
Pointed/pointed, straight		
145mm	100 4549	17.90
Slim model (not shown)		
Pointed/blunt		
straight, 145 mm	900 3111	13.90
curved, 145 mm	900 3103	13.90



default

Size	Order no.	Price
Blunt/blunt, curved		
130mm	900 3107	12.90
145mm	900 3106	14.90
Blunt/blunt, straight		
130mm	900 3096	12.90
145mm	900 3095	14.90



Metzenbaum dissecting scissors

Size	Order no.	Price
Bent		
145mm	100 0511	11.90
145mm - Maxima Line	100 8335	19.90
180mm	100 5632	13.90
Straight		
145mm	100 9093	11.90



Fine dissecting scissors

Size	Order no.	Price
Bent		
90mm	900 3091	10.90
105mm	900 3090	11.90
115mm	900 3093	12.90
Straight		
90mm	900 3092	10.90
115mm	900 3094	11.90
150mm	900 3097	14.90



Piece off
11.90
PRICE-PERFORMANCE

Lister bandage scissors

Size	Order no.	Price
90mm	100 9243	11.90
115mm	900 3123	12.90
140mm	100 1087	14.90
140mm - Maxima Line	100 6589	19.90
180mm	900 3087	18.90
180mm - Maxima Line	101 5267	24.90
185mm	100 3334	19.90
200mm	900 3086	20.90



Lister bandage scissors with clip


140mm	100 3465	17.90
-------	----------	-------



Multi-purpose scissors

Made of stainless steel. With plastic handles and serrated cutting edges.
Can be sterilized up to 140 °C.

Size	Order no.	Price
145mm	900 3124	7.50
190mm	100 7898	8.50



Kelly gum scissors

Art	Size	Order no.	Price
Straight	160mm	900 2072	10.50
Bent	160mm	900 2025	10.50




Littauer thread scissors

Size	Order no.	Price
115mm	100 6306	14.90
140mm	100 5115	16.90



Surgical scissors
according to Kelly Bent

115mm	100 7266	15.90
-------	----------	-------



Scissors for mayo

Art	Size	Order no.	Price
Straight	140mm	100 1106	12.90
Straight	145mm	100 4018	14.90
Straight	170mm	100 1237	16.90



Spencer ligature scissors

90mm	900 3089	14.90
------	----------	-------



Needle holder according to Derf

120mm	100 7035	9.50
-------	----------	------

Hegar-Olsen needle holder

140mm	100 1125	15.90
-------	----------	-------

Mayo-Hegar needle holder

Size	Order no.	Price
130mm	100 0959	13.90
140mm - Maxima Line	100 8670	22.90
150mm	100 2570	14.90
180mm	100 8192	15.90
200mm	100 6577	17.90
With carbide insert		
150mm - Maxima Line	100 6116	41.90

Needle holder according to Mathieu


Size	Order no.	Price
140mm	900 2079	24.90
165mm	900 2080	26.90
170mm	900 3597	27.90
200mm	900 3120	29.90
With carbide insert		
140mm	900 6214	49.90
170mm	900 6215	53.90
200mm	900 6216	56.90

Castroviejo needle holder

Size	Order no.	Price
With carbide insert		
140mm	900 8637	52.90
170mm	900 8638	56.90

Crile Wood needle holder

Size	Order no.	Price
150mm	100 2492	13.90
150mm - Maxima Line	101 7128	24.90
With carbide insert		
150mm - Maxima Line	100 7943	47.90





Anatomical tweezers

Size	Order no.	Price
Straight		
115mm - Maxima Line	102 3437	14.90
130mm	900 2011	6.99
140mm	100 5152	7.19
145mm - Maxima Line	101 1272	16.90
150mm	100 6424	7.99
200mm	100 9558	8.99
Narrow design		
120mm	900 3150	10.90
130mm	900 3148	11.90
145mm	900 3147	12.90
160mm	900 3146	13.90
180mm	900 3149	14.90



NEW

College tweezers length 15 cm, stainless steel.

Art	Size	Order no.	Price
Short tip, smooth	Tip length 12 mm	900 2059	7.50
Short tip, grooved	Tip length 12 mm	900 1646	7.90
Long tip, grooved	Tip length 16 mm	100 3313	8.90

Tissue tweezers
Adson, 1:2 teeth

Size	Order no.	Price
120mm	100 5891	10.50



Surgical forceps

Size	Order no.	Price
1:2 teeth		
130mm	900 2009	7.99
140mm	100 4877	8.99
150mm	100 9871	9.29
180mm	900 2010	9.99
2:3 teeth		
140mm	100 0389	8.99
3:4 teeth		
140mm	100 9796	8.99
4:5 teeth		
140mm	100 1179	8.99




NEW

Meriam tweezers with stop pin


Length 16cm	fluted	100 5709	8.50
-------------	--------	----------	------

Adson tweezers
Anatomical

Size	Order no.	Price
120mm	100 5659	9.90
120mm - Maxima Line	101 6279	14.90

Rochester-Pean artery forceps
Matted, without teeth

Size	Order no.	Price
Bent		
130mm	900 3073	13.90
145mm	900 3072	14.90
160mm	100 2005	15.90
180mm	900 3071	16.90
200mm	100 0226	17.90
Straight		
130mm	900 3077	13.90
145mm	900 3076	14.90
160mm	100 0629	15.90
180mm	900 3075	16.90
200mm	100 2838	17.90




piece only

19.90

PRICE-PERFORMANCE

Baby mosquito
Anatomical

100mm	bent	102 4346	19.90
-------	------	----------	-------




Piece off

8.49

PRICE-PERFORMANCE

Cloth clamp according to Backhaus

Size	Order no.	Price
90mm	100 0107	8.49
140mm	100 0341	9.49


Arterial clamp
Halsted Mosquito
Anatomical

Size	Order no.	Price
Bent		
90mm	100 7530	11.90
125mm	100 4862	13.90
125mm Maxima Line	100 2376	23.90
Straight		
90mm	100 5499	11.90
125mm	100 4662	13.90
120mm Maxima Line	100 0846	23.90

Halsted-Mosquito artery forceps (not shown)

Surgical, 1x2 teeth

Bent		
125mm	900 3078	13.90
Straight		
125mm	900 3079	13.90

Kocher artery forceps
Matted, with teeth

Size	Order no.	Price
Bent		
130mm	900 3067	13.90
145mm	900 3066	14.90
Straight		
130mm	900 3070	13.90
145mm	900 3069	14.90
160mm	900 3068	15.90



Skin punch

Makes taking skin and tissue samples easier. Individually sterile packed. Pack (10 pieces)

Size	Order no.	Price
ø 2 mm	548 176	27.90
ø 3 mm	548 178	27.90
ø 4 mm	548 179	27.90
ø 5 mm	548 181	27.90
ø 6 mm	548 182	27.90
ø 7 mm	624 641	27.90
ø 8 mm	548 183	27.90



Volkman sharp spoon

Size	Figure	Order no.	Price
2x3mm	Fig. 0000	900 3590	23.90
3x2mm	Fig. 000	900 3591	23.90
4x6mm	Fig. 00	900 3592	23.90
5x8mm	Fig. 0	900 3593	23.90
7x10mm	Fig.	900 3595	23.90
8x12mm	1 Fig. 2	900 3594	23.90



Martini Sharp Spoon

Size	figure	Order no.	Price
Round, length: 140mm			
2.5x3mm	Fig. 00-0	900 3117	16.90
3 x 3.5mm	Fig. 0-1	900 3118	16.90
Willing Sharp Spoon (not pictured)			
Oval, length: 140mm			
2.5x3mm	Fig. 00-0	900 3116	15.90
3 x 3.5mm	Fig. 0-1	900 3115	15.90



Sharp spoon Hemingway

6 mm round handle, oval spoon shape. Length 18cm.

Size	Order no.	Price
small	900 2065	7.50
medium	900 2066	8.50
large	900 2064	9.50



Raspatorium Williger

length 16cm; Width 6.5mm

900 1745 12.90



Disposable vaginal speculum

Specula made of medical plastic. Individually sterile packed.
Size information clearly visible on the packages.
Pack (25 pieces)

Size	Order no.	Price
s	900 9406	18.90
M	900 9407	18.90
M long	988 0525	18.90
L	900 9408	18.90
With snap closure, pack of 100		
M	988 2141	49.90



Sterilization bag paper/

clear film without gussets with excellent self-adhesive closure. Applied indicator strips for chemical and steam sterilization, as well as EO sterilization.

Size	Crowd	Order no.	Price
Self-adhesive			
8.3x15.9cm	Pack (200 pieces)	900 2984	11.90
8.9 x 25.4cm	Pack (200 pieces)	900 2983	12.90
13.4 x 28 cm	Pack (200 pieces)	900 2985	19.90
19.0 x 35.6 cm	Pack (200 pieces)	988 2127	26.90
25.4 x 38.1cm	Pack (200 pieces)	900 2986	30.90
Not self-adhesive			
5.0 x 30 cm	Pack (500 pieces)	900 3059	39.90
7.5x30cm	Pack (500 pieces)	900 3060	49.90
10x30cm	Pack (500 pieces)	900 3058	59.90



Steam and gas sterilization packaging

From the role

In a paper-film combination with an indicator already printed on it.
The transparent film corresponds to DIN standard 50049. The paper for medical use corresponds to DIN standard 58953. The welding lines (multiline) are very resistant and arranged at even intervals.

- > For autoclaves
- > Germ-proof
- > Tear-resistant
- > Transparent

Size	Order no.	Price
Roll length: 100 m		
5.0 cm wide	900 2987	9.95
7.5 cm wide	900 2988	14.95
10.0 cm wide	900 2989	19.95
20.0 cm, height 6.5 cm	900 2990	34.95
Roll length: 200 m		
5.0 cm wide	900 2991	20.90
7.5 cm	900 2992	27.90
wide 10.0 cm wide	900 2993	32.90
15.0 cm wide	988 0485	45.90
20.0 cm wide	988 0486	59.90
25.0 cm wide	988 0487	69.90
30.0 cm wide	988 0488	79.90



Container filter paper

With indicator, for single use

Hydrophobic round container filter papers, according to DIN EN 868-2, with indicator, which turns dark brown after sterilization.

Size: Ø 19mm. Color white

Pack (100 pieces)	988 0901	10.90
from 5 packs each		9.90



Sterilization marker

After sterilization, the writing changed from wine red to green. The sterilization marker complies with EN 867-1, Class A and ISO 11140-1, Class 1.

Piece	900 8492	13.50
-------	----------	-------



NEW

piece only

749.00

PRICE-PERFORMANCE

See cleaning agents from
page 30

1 ultrasonic cleaning device Sweep 2000

The Sweep 2000 ultrasonic cleaning device is an innovative device with the latest technology. Particularly noteworthy are the automatically locking lid, the larger drain, the unique and newly designed hanging basket which allows for hand-free draining, and the digital control panel. The new housing makes cleaning quick and easy. The device is ideal for the cleaning process in medical practices, as it standardizes the cleaning cycle thanks to the new hinged and snap-in lid. The new lid also ensures a closed cleaning cycle and very low operating noise.

Technical data:

Power: 43Hz 95 watts; Capacity: 4 liters; Timer: 60 minutes; Heating: 55 watts; Dimensions: 37.5 x 26.7 x 26.7 cm; Tank dimensions: 28.6 x 14 x 8.9 cm; Weight: 6.8kg

Scope of delivery:

1x cleaning unit + lid; 1x drain + hose; 1x power cable; 1x instruction manual

Piece	900 6375	749.00
-------	----------	--------



Accessories set

Contents: 1 hanging basket, 1 perforated lid for 2 cups, 1 perforated lid for 1 Beaker, 2 glass cups (600 ml), 1 drill box and 1 drill tray

set	900 6470	239.00
-----	----------	--------

2 perforated lids

Perforated lid for 1 cup	979 2678	39.90
Perforated lid for 2 cups	979 2679	49.90

3 hanging baskets

Piece	979 2676	129.90
Hanging basket large unit	979 2846	109.90

Replacement lid (not shown)

Piece	979 2677	79.90
-------	----------	-------

4 glass cups

Glass mug 600 ml	979 2849	28.90
------------------	----------	-------

5 rubber band for glass cups

Rubber band for glass cups 600 ml	979 2851	22.90
-----------------------------------	----------	-------

Sieve

6 box for drills	979 2852	49.90
7 basket for drills	979 2853	25.90



Set only

92.90


PRICE-PERFORMANCE

Bowie & Dick test

The safe and economical solution for weekly autoclave checks. The successful steam penetration and the exposure time are confirmed. Consisting of metal test cassette, loading cards and 20 test cards. According to DIN EN ISO 11140-1, Class 2.

Set consisting of test cassette and 20 tests.

Size	Order no.	Price
Metal cassette and 20 test cards		
set	900 8646	92.90
Refill pack test cards		
20 piece	900 8647	44.90
Replacement metal cassette		
Piece	900 8648	76.50



NEW

Wash test test bar

The HS-Wash test test strip is required for the HS test indicator strips. It is made of stainless steel. It can be used multiple times, making it economical and environmentally friendly.

80 x 10 x 28 mm (W x D x H)	Piece	988 2299	51.90
-----------------------------	-------	----------	-------




Piece off

4.49

PRICE-PERFORMANCE

Indicator bands for steam sterilization

Size	Order no.	Price
50mx1.2cm (L x W)	900 1837	4.49
50mx1.9cm (L x W)	900 1838	4.99
50mx 2.5cm (L x W)	900 1839	5.99



Helix test

Test body including 250 test strips for steam sterilization. Helix test validates any sterilization process in a filled autoclave, air removal and no non-condensable gases (NCG). The indication ink is non-toxic and the adhesives do not contain latex. Corresponds to DIN standard EN 867-5, ISO 11140-1, class 2.

Size	Order no.	Price
Set (test body with 250 test strips)	900 8649	89.90
Test body	900 8651	44.90
250 test strips	900 8650	44.90



Pack off
15.90
PRICE-PERFORMANCE

fixing plaster

Skin-colored roll plaster, length: 5 m, inelastic, very tensile cotton fabric, breathable, permeable to water vapor, for fixing any type of bandages, cannulas, etc., 100% cotton.

Size	Crowd	Order no.	Price
5m x 1.25cm	18 pcs	900 8600	15.90
5m x 2.50cm	12 pieces	900 8601	17.90
5m x 5cm	6 pieces	900 8602	19.90



Pack off
15.90
PRICE-PERFORMANCE

Non-woven plaster

Fixation plaster made of white cellulose acetate, breathable, water vapor permeable, for fixing any type of bandages, compresses, cannulas, etc., 100% non-woven fabric, microporous. Color white

Size	Crowd	Order no.	Price
5m x 1.25cm	18 pcs	900 8606	15.90
5m x 2.5cm	12 pieces	900 8607	16.90
5m x 5cm	6 pieces	900 8608	18.90



Pack off
17.90
PRICE-PERFORMANCE

Silk plaster

Fixation plaster made of white cellulose acetate, length: 5 m, breathable, permeable to water vapor, for fixing any type of bandages, cannulas, etc., 100% cotton and cellulose fibers.

Size	Crowd	Order no.	Price
5m x 1.25cm	18 pcs	900 8603	17.90
5m x 2.5cm	12 pieces	900 8604	21.90
5m x 5cm	6 pieces	900 8605	19.90



Pack off
14.90
PRICE-PERFORMANCE

Transparent plaster

Roll plaster made of transparent, supple polyethylene film, breathable, permeable to water vapor, for fixing any type of bandages, cannulas, etc.

Size	Crowd	Order no.	Price
5m x 1.25cm	18 pcs	900 8609	14.90
5m x 2.50cm	12 pieces	900 8610	16.90
5m x 5cm	6 pieces	900 8611	18.90



Pack off
5.49
PRICE-PERFORMANCE

Hypoallergenic

Breathable

Injection patches

To cover the injection site. Hypoallergenic and breathable for sensitive skin. Dimensions: 2 x 6 cm

Pack (200 pieces)	900 8612	6.90
from 3 packs each		5.49



Piece off
5.99
PRICE-PERFORMANCE




Quick wound dressings

Quick wound dressing with a centrally placed, absorbent compress made of 100% viscose, covered with special fleece to prevent it from sticking to the wound, for initial care and the treatment of minor wounds. Breathable, permeable to water vapor.

Size	Crowd	Order no.	Price
1 skin tones			
5m x 4cm	Piece	900 8613	5.99
5m x 6cm	Piece	900 8614	6.99
5m x 8cm	Piece	900 8615	8.99
2 White			
5m x 4cm	Piece	900 8616	5.99
5m x 6cm	Piece	900 8617	6.99
5m x 8cm	Piece	900 8618	8.99

Gauze pad
Made from gauze.

Size	Crowd	Order no.	Price
1 Non-sterile			
hazelnut	500 pieces	900 8739	13.50
walnut	500 pieces	900 8740	16.50
plum	500 pieces	900 8741	22.50
egg	500 pieces	900 8742	34.90
large	250 pieces	900 8743	22.90
extra large	250 pieces	900 8744	32.90
2 Sterile			
pea	500 pieces	900 3192	13.50
hazelnut	500 pieces	900 3193	14.50
walnut	500 pieces	900 3194	16.50
plum	300 pieces	900 3195	12.90
Sterile – packed in packs of 5			
hazelnut	20 x 5 pieces	900 9841	5.99
walnut	20 x 5 pieces	900 9842	6.55
plum	20 x 5 pieces	900 9843	7.55
egg	20 x 5 pieces	900 9844	7.99
large	20 x 5 pieces	900 9845	8.45



Piece off
6.99
PRICE-PERFORMANCE
Gently removable

Fixing fleece

Fully adhesive, slightly elastic fixing fleece. Skin-friendly, permeable to air and water vapor. Does not constrict and adheres non-slip to the skin. High flexibility. 100% polyester fleece, color: white

Size	Crowd	Order no.	Price
5cm x 10m	Piece	900 8619	6.99
10cm x 10m	Piece	900 8620	9.99
15cm x 10m	Piece	900 8621	16.99
20cm x 10m	Piece	900 8622	18.99



100 pieces from
1.99
PRICE-PERFORMANCE

1 gauze pad

Made of gauze with folded cut edges. Laid 8 and 12 times. Absorbent. The non-sterile compresses can be autoclaved.

Size	Crowd	Order no.	Price
Non-sterile, placed 8 times			
5x5cm	100 pieces	900 3180	1.99
7.5x7.5cm	100 pieces	900 3181	3.99
10x10cm	100 pieces	900 3179	4.99
Non-sterile, placed 12 times			
5x5cm	100 pieces	900 3177	2.99
7.5x7.5cm	100 pieces	900 3178	5.99
10x10cm	100 pieces	900 3175	6.49
Sterile, placed 12 times			
5x5cm	50 x 2 pieces	900 3223	4.99
5x5cm	100 pieces	900 3189	5.99
7.5x7.5cm	40 x 2 pieces	900 3224	6.49
7.5x7.5cm	100 pieces	900 3190	8.99
10x10cm	25 x 2 pieces	900 3225	5.99

2 fleece compresses

Versatile compress made from synthetic polyester fabric.


Size	Crowd	Order no.	Price
Non-sterile, placed 4 times			
5x5cm	200 pieces	900 3263	4.49
7.5x7.5cm	200 pieces	900 3308	4.49
10x10cm	200 pieces	900 3309	6.49



cotton balls

Chlorine-free beads with high absorption capacity, ideal for drying cavities and applying medications.

Size	Quantity	Order no.	Price
00, Ø 4.0 mm	4g	900 1879	16.90
0, Ø 4.8 mm	4g	900 1880	12.90
1, Ø 6.3 mm	10g	900 1881	10.90
2, Ø 7.6 mm	10g	900 1882	11.90



Maxima® cotton roll dispenser with button

Pressure dispenser with button for 50 cotton rolls. By briefly pressing the lid, each cotton roll is dispensed individually. All other roles remain in the dispenser. The heavy, non-slip base ensures stability. Top part made of high-quality acrylic glass.

Note: Only suitable for sizes 1 and 2 cotton rolls. Dispenser is delivered empty.

Size	Order no.	Price
blue	900 6050	6.90
white	900 6049	6.90
red	900 6051	6.90



cotton rolls

White, chlorine-free bleached cotton rolls made of cotton and cellulose with high absorption capacity.

Size 1,	Quantity	Order no.	Price
Ø 8 mm	300g	900 1895	9.50
from 5 packs each			8.50
2, Ø 10 mm	300g	900 1896	10.50
from 5 packs of 3			9.50
each, Ø 12	300g	900 1897	11.50
mm from 5 packs each			10.50



Elastic gauze/fixation bandage

Elastic gauze fixation bandage with woven edges for single use. 45% cotton, 55% polyamide. Separately packaged.

Size	Crowd	Order no.	Price
4m x 4cm	20 piece	988 0144	6.99
4m x 6cm	20 piece	988 0145	7.99
4m x 8cm	20 piece	988 0146	8.99
4m x 10cm	20 piece	988 0147	9.99



Elastic gauze/fixation bandage

Elastic fixing bandage according to DIN 64634 with a surface similar to crepe. 70% viscose, 30% polyamide. High elasticity, very breathable, absorbent, soft and non-slip. Stretch approx. 130%. Color white

Size	Crowd	Order no.	Price
4m x 4cm	20 piece	900 8569	5.99
4m x 6cm	20 piece	900 8570	6.99
4m x 6cm	100 pieces	900 8574	29.90
4m x 8cm	20 piece	900 8571	8.09
4m x 8cm	100 pieces	900 8575	39.90
4m x 10cm	20 piece	900 8572	9.99
4m x 10cm	50 pieces	900 8576	25.90
4m x 12cm	20 piece	900 8573	10.90



Cohesive crepe bandage

Longitudinally elastic fixation bandage, non-slip, no final fixation necessary, sterilizable. Stretch approx. 120%, 30% viscose, 30% polyamide, 40% cotton. Color white

Size	Crowd	Order no.	Price
20m x 4cm	Piece	900 8585	6.50
20m x 6cm	Piece	900 8586	6.90
20m x 8cm	Piece	900 8587	7.90
20m x 12cm	Piece	900 8589	9.50



Creped elastic gauze/fixation bandage

High elasticity, resistant to ointments, sweat, oils and fats, very breathable, 40% viscose, 60% polyamide. Absorbent, soft, non-slip, sterilizable, stretch approx. 135%. Color white


Size	Crowd	Order no.	Price
4m x 6cm	100 pieces	900 8582	29.90
4m x 8cm	100 pieces	900 8583	39.90
4m x 10cm	50 pieces	900 8584	25.90



Ideal bandages

Permanently elastic ideal bandages for relief bandages and follow-up treatment of dislocations and fractures. Material: 88% cotton, 11% polyamide and 1% polyurethane. Porous, breathable. High cotton content. Web-edged. Gas and steam sterilizable at 135°C. Color white

Size	Crowd	Order no.	Price
5m x 4cm	10 pieces	900 1604	14.90
5m x 6cm	10 pieces	900 1605	19.90
5m x 8cm	10 pieces	900 0808	23.90
5m x 10cm	10 pieces	900 0809	29.90
5m x 12cm	10 pieces	900 0810	33.90



Piece off
4.49
PRICE-PERFORMANCE

Tape bandages

Inelastic, particularly tensile plaster bandage made from 100% cotton fabric, for treating contusions, distortions, dislocations and fractures as well as reinforcing compression bandages.

Size	Crowd	Order no.	Price
white			
10 mx 2.5 cm	Piece	988 2022	4.99
from 10 pieces each			4.49
10m x 3.8cm	Piece	988 2023	5.99
from 10 pieces each			5.59
10m x 5.0cm	Piece	988 2024	8.99
from 10 pieces each			8.59
blue			
10m x 3.8cm	Piece	988 2025	5.99
red			
10m x 3.8cm	Piece	988 2026	5.99
yellow			
10 mx 3.8 cm	Piece	988 2027	5.99
from 10 pieces per color			5.59



Compression bandages fine

With a long train

Latex-free, permanently elastic, long and soft stretch, length: 7 m, stretch approx. 165%, sterilizable, 81% cotton, 16% polyamide, 3% polyurethane, color: skin

Size	Crowd	Order no.	Price
7m x 6cm	8 pieces	900 1336	33.90
7m x 8cm	8 pieces	900 1338	44.90
7m x 10cm	8 pieces	900 1340	55.90



10 pieces from
38.90
PRICE-PERFORMANCE

Compression bandages

With a very short pull

Latex-free, length: 5 m, stretch approx. 40%, sterilizable, 67% cotton, 33% polyamide, color: skin

Size	Crowd	Order no.	Price
5mx6cm	10 pieces	900 1348	38.90
5mx8cm	10 pieces	900 1350	47.90
5mx10cm	10 pieces	900 1352	58.90



10 pieces from
11.90
PRICE-PERFORMANCE

Universal bandages

63% cotton, 36% polyamide, 1% polyurethane. Permanently elastic and sterilizable. Stretchability approx. 100%, color: white

Size	Crowd	Order no.	Price
5m x 6cm	10 pieces	900 0842	11.90
5m x 8cm	10 pieces	900 0843	13.90
5m x 10cm	10 pieces	900 0844	14.90
5m x 12cm	10 pieces	900 0845	17.90



Compression bandages strong

With a long train

Latex-free, permanently elastic, long stretch, length: 7 m, stretch approx. 175%, sterilizable, 93% cotton, 7% polyamide, color: skin

Size	Crowd	Order no.	Price
7m x 8cm	8 pieces	900 1344	58.90
7m x 10cm	8 pieces	900 1346	62.90



Piece off
6.19
PRICE-PERFORMANCE

Tubular bandage

Longitudinal and transverse elastic. Roll length: 20 m

Size	Crowd	Order no.	Price
Toe and children's bandage			
Gr. 0/1.5cm	Piece	900 2578	6.19
Finger bandage			
Gr. 1/2.5cm	Piece	900 2579	6.19
Arm and child's leg			
Gr. 2 / 6	Piece	900 2580	10.49
cm foot, leg, armpit bandage; Child's head			
Gr. 3/8cm	Piece	900 2581	11.99
face mask; armpit; Leg and head bandage			
Gr. 4/10 cm	Piece	900 2582	15.99
Strong thigh and child's torso			
Gr. 5/12cm	Piece	900 2583	19.99
Child's torso and thighs			
Gr. 6/16 cm	Piece	900 2584	23.99
Desault trunk bandages			
Gr. 7/21cm	Piece	900 2585	34.99



10 pieces from
49.90
PRICE-PERFORMANCE

Zinc glue bandages

Longitudinal elastic with very good compression effect. Easy to model, quick-drying, semi-rigid and skin-friendly. Stretch approx. 60%, 71% viscose, 29% polyamide, color: white. Pack (10 pieces)

Size	Crowd	Order no.	Price
5m x 10cm		900 8594	49.90
7m x 10cm		900 8595	59.90
10m x 10cm		900 8596	69.90



Piece off
0.64
PRICE-PERFORMANCE

Synthetic padding bandages

Made of skin-friendly fleece, breathable, wrinkle-free, for padding and cushioning. 100% polyester unbleached, color: white

Size	Crowd	Order no.	Price
3m x 6cm	108 pieces	900 8597	68.90
3m x 10cm	72 pieces	900 8598	72.90
3m x 15cm	36 pieces	900 8599	76.90



Piece off
6.89
PRICE-PERFORMANCE

Acrylic bandages

Plaster bandages coated with transparent rubber adhesive for permanent compression bandages and functional bandages. Stretch: approx. 70%.

Size	Crowd	Order no.	Price
2.5m x 6cm	Piece	900 0839	6.89
2.5m x 8cm	Piece	900 0840	7.99
2.5m x 10cm	Piece	900 0841	8.89



SUSTAINABLE PRODUCTS
BY HENRY SCHEIN



Practice Green

Do you already know the new, sustainable product line from Henry Schein? It includes both products made from renewable raw materials also biodegradable products. For example:



PAGE 34

Medical rolls made from recycled paper



PAGE 53

Paper cups biodegradable



PAGE 21

EuroSept Xtra Organic Quick disinfectant wipes



PAGE 66

Acclean® Bamboo toothbrush



PAGE 66

Silk dental floss

You can find many other sustainable products in our webshop

<https://henryscheinmed.at/p/practice-green>

Protecting the planet for a healthier future



NEW

Paper cups, biodegradable

Cup made of PLA paper, biodegradable. The cups offer an environmentally friendly alternative to conventional paper cups. These cups were specially marked with the "green" leaf as a reminder that the cups can be composted. The paper cups made from renewable materials are made of cardboard that is considered suitable for commercial composting according to the EN13432 DIN-Certco standard. Advantages: Made from renewable raw materials! Compostable!

Technical: RIM ROLL Ultra-tight roll designed for superior quality: such as sidewall rigidity to ensure a secure grip.

Capacity 180 ml.		
Pack (1,000 pieces)	988 4317	46.90



NEW

3,000 pieces only
67.90
PRICE-PERFORMANCE

Mug

Coloured

Size	Order no.	Price
180 ml content. Pack (3,000 pieces)		
blue	900 6532	67.90
yellow	900 6533	67.90
orange	900 6536	67.90
	900 8845	67.90
pink lime green	988 3011	67.90
turquoise	988 3014	67.90
black	988 3580	67.90

Cups sorted by color

Assorted in different colors (blue, yellow, green, orange, pink, purple). 600 pieces per color. 180 ml content. Pack (3,000 pieces)

Sorted	900 8366	89.00
--------	----------	-------



NEW

HS paper cups recyclable

Mouthwash cups meet EU wood protection regulations. Technical: RIM ROLL Ultra-tight roll designed for superior quality: such as sidewall rigidity to ensure a secure grip. Internal coating: made with food grade PE (polyethylene). Side wall: made with *PEFC-certified, sustainably sourced cardboard. *PEFC is the largest institution for ensuring and marketing sustainable forest management through an independent certification system. Wood and wood products with the PEFC seal demonstrably come from ecologically, economically and socially sustainable forestry.

Note: We recommend not using the cups with sodium hypochlorite or CHX for a longer period of time.

Size	Order no.	Price
180 ml content. Pack (1,000 pieces)		
blue 988 4318 white from 3 packs each		44.90
	988 4319	44.90
		41.90



Mug

Color white

Size	Order no.	Price
150 ml content		
Pack (100 pieces)	988 2019	3.89
180 ml content		
Pack (3,000 pieces)	900 6535	67.90
200 ml content		
Pack (100 pieces)	988 2020	3.99




Cotton swab large
With extra large cotton head
Cotton head: 3 cm x 10 mm (L x Ø)
Wooden stick: 20cm x 3 mm (L x Ø)

Pack (50 pieces)	988 9886	1.99
from 10 packs each		1.39



cotton swab
Cotton head: 1 cm x 5 mm (L x Ø)
Wooden stick: 15cm x 2 mm (L x Ø)

Pack (100 pieces)	988 9887	1.19
from 10 packs each		0.99



cotton swab
With cotton head on both sides. Length: 80mm

Pack (100 pieces)	988 9884	0.99
from 5 packs each		0.89



500 pieces only
5.89
PRICE-PERFORMANCE

Wooden tongue depressor
Spatula made of extremely flexible and break-resistant birch wood, with a smooth, absolutely splinter-free surface.

Pack (500 pieces)

Size	Order no.	Price
For adults, 15 cm	988 9885	5.89
For children, 14 cm	988 9883	5.89



Stool tube
With a straight floor. Volume: 18 ml.
Dimensions: 22 x 63 mm (Ø x L)

Bags (100 pieces)	900 9740	17.90
from 5 bags each		15.90



100 pieces from
10.90
PRICE-PERFORMANCE

Safely protect against
needlestick injuries

Safety blood lancets
Sterile safety blood lancets with a 1.8 mm deep puncture. The needle is not visible before or after use, reducing cross-infection. Pack (100 pieces)

Size	Quantity	Color	Order no.	Price
18G	1.8mm	green	900 7993	13.90
21G	1.8mm	blue	900 7992	13.90
28G	1.8mm	purple	900 7991	13.90
30G	1.8mm	turquoise	570 0676	13.90
21G	2.2mm	orange	900 7980	13.90
21G	2.8mm	pink	900 7982	13.90
from 5 packs per variety				10.90



Piece off

5.59

PRICE-PERFORMANCE

safety goggles

With anti-fog coating for versatile use in the operating room, laboratory or practice.

Size	Order no.	Price
1 safety goggles yellow lenses	979 3109	5.59
2 safety glasses clear lenses	979 3111	5.59
3 safety glasses clear, light, adjustable	979 3112	5.99
4 safety glasses black, light, adjustable	979 3113	5.99
5 safety glasses laboratory clear	979 3114	5.99



Antifog
UV protection glasses

Reduced UV rays, anti-fog coating, particularly large protective lenses, adjustable temples, side ventilated, suitable for people who wear glasses.

Color	Order no.	Price
orange	979 3252	12.90
from 3 packs each		9.90



Glasses frame with 6
protective shields

Glasses frame and 6 visors in the respective color.
Autoclavable up to 121°C. Pack: glasses frame with 6 visors.

Size	Order no.	Price
white pink-transparent	572 2218	29.50
green-transparent	572 2217	29.50
blue-transparent	572 2219	29.50
from 3 packs per variety	572 2216	29.50
from 3 packs per variety		27.90



NEW

Face shield replacement lenses

Comfortable eye protection that allows clear vision. Spring light, suitable for HS glasses frame/frame.

25x16.5cm	25 pieces	900 9630	17.80
-----------	-----------	----------	-------



100 pieces from
12.90
PRICE-PERFORMANCE

Disposable aprons

Disposable apron made of PE. Perforated for individual removal. Color white. Thickness: 0.014mm
Pack (100 pieces)

Length: 75 x 125	988 1842	13.90
cm from 10 packs each		12.90
Length: 75 x 140	988 1843	14.90
cm from 10 packs each		13.90
Length: 75 x 160	988 1844	15.90
cm from 10 packs each		14.90



One Step HCG Urine Test Strips

The result is available in 15 minutes.
Pack (50 test strips)

Pack	988 9889	9.90
of 3 packs each		8.90



One Step Strep A

The Strep A Rapid Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of Strep A antigens from throat swab samples to facilitate the diagnosis of group A streptococcal infection.

Pack (20 test cassettes)

Pack	988 9891	23.90
of 3 packs each		22.50



One Step HCG Urine Test Cassettes

The result is available in 3 minutes.
Pack (40 test cassettes)

Pack	988 9888	14.90
of 3 packs each		13.50



One Step Multi 5 Drug test cassettes

A rapid test for the simultaneous, qualitative detection of multiple drugs and drug metabolites in human urine.
Testing any combination of the following drugs:
ACE/AMP/BAR/BZO/BUP/COC/THC/MTD/MET/MDMA/MOP/MQL/OPI/PCP/PPX/TCA/TML/KET/OXY/
COT/EDDP/FYL/K2/6-MAM/MDA/ETG/CLO/LSD/MPD/ZOL/MEP/ALC/MDPV/DIA/ZOP/MCAT/7-ACL/
CFYL/CAF/CAT/TRO/ALP/PGB/COD/ZAL/MPRD/CNB/GAB/TZD/CAR/ABP/QTP/FLX/UR-144/KRA/
TLD-/PVP/MES

Including sample validity tests (SVT) for: oxidants/PCC, specific gravity, pH, nitrite, glutaraldehyde, creatinine and bleach

Pack (25 test strips)

Pack	988 9890	61.90
of 3 packs each		55.90

It depends on the content -

Affordable 100 piece packs!



Urine test strips Plus

Blood and glucose test field suppressed with ascorbic acid.
> For visual urine analysis
> With leukocyte field, long shelf life

Test strips														Order no.	Crowd	Price	from 10 packs each
Urispec Plus 3	•	•	•						•					900 9960	100 pieces	6.90	5.90
Urispec Plus 5+Leuko	•	•		•		•		•			•			900 9961	100 pieces	7.90	6.90
Urispec Plus 9+Leuko	•	•	•	•	•	•	•	•	•		•			900 9957	100 pieces	9.90	8.90
Urispec Plus 10 SL	•	•	•	•	•	•	•	•		•	•			900 9958	100 pieces	11.90	10.90
Urispec Plus 11 SYS	•	•	•	•	•	•	•	•	•	•	•			900 9959	100 pieces	12.90	11.90



Pediatric stethoscope
For the small patient: Smaller double chest piece with cold protection ring. Latex free.

1 black	Piece	988 0827	11.90
from 3 pieces each			9.90

Infant stethoscope

2 gray	Piece	988 0828	11.90
from 3 pieces each			9.90



Maxima Cardiology Stethoscope
> Extra deep bell for an unmatched sound
Bass reproduction
> Ultra-sensitive diaphragm for greater gain and sharper high frequency transmission
> Non-shower funnel and membrane retaining edge for patient comfort

Piece	988 3146	39.90
-------	----------	-------




Rappaport stethoscope
The Rappaport model is characterized by its excellent hearing quality. The included accessory kit includes 5 different interchangeable chestpiece accessories (3 horns and 2 diaphragms) for capturing a wide range of frequencies in pediatric and adult patients and 2 pairs of eartips. Latex free. Hose length: 55 cm.

Color	Order no.	Price
1 black	988 0815	15.90
2 navy blue	988 0816	15.90
3 fresh orange	988 0817	15.90
4 royal blue	988 0818	15.90
from 5 pieces per color each		13.90
spare Parts		
Membrane, large	988 0832	1.49
Membrane, small	988 0833	1.49
ear olives	988 0834	2.49
Ear tips, transparent	988 0835	2.49




Flat stethoscope
Stethoscope that can be easily slipped under a blood pressure cuff. Latex free. Hose length: 55 cm.

Color	Crowd	Order no.	Price
1 black	Piece	988 0819	5.49
2 gray	Piece	988 0820	5.49
3 royal blue	Piece	988 0821	5.49
4 red	Piece	988 0822	5.49
from 5 pieces per color each			4.49




Double head stethoscope
Comfortable double-head stethoscope with a lightweight chestpiece. Cold protection ring for better patient comfort. Latex free.

Color	Crowd	Order no.	Price
1 black	Piece	988 0823	11.90
2 gray	Piece	988 0824	11.90
3 royal blue	Piece	988 0825	11.90
4 red	Piece	988 0826	11.90
from 5 pieces per color each			9.90



Maxima stethoscope
The professionals' preferred stethoscope at an unrivaled price. Combination chest piece for all frequency ranges: membrane for listening to high-frequency noises – open funnel for listening to low-frequency noises. Combination chest piece and ear hook made of stainless steel. Latex free. The included accessory set includes a replacement membrane and 2 pairs of ear tips.

Size 1	Order no.	Price	spare Parts
black	988 0812	43.90	membrane 988 0829 1.99
burgundy	988 0813 3	43.90	Ear tips, large 988 0830 3.19
navy blue	988 0814	43.90	Ear tips, small 988 0831 3.19



piece only
23.90
PRICE-PERFORMANCE

Blood pressure monitor
For the practice
> Inexpensive hand model with two-hose cuff

Piece	900 9733	23.90
-------	----------	-------



piece only
299.00
PRICE-PERFORMANCE
SpO2 and pulse

Pulse oximeter
Measurement of oxygen saturation SpO2 and pulse PR in adults.
With alarm limit functions (visual and acoustic). Sensor error alarms. Storage function of up to 72 hours of measurement data, battery operation 3AA. Dimensions: 125x60x 30 mm, weight: 141 g. Incl. finger clip sensor for adults and batteries.

	Order no.	Price
Piece (incl. finger clip sensor)	900 8275	299.00
Finger sensor with clip	900 8277	49.90
Soft finger sensor, made of silicone	900 8278	49.90



piece only
97.50
PRICE-PERFORMANCE

Includes bag

Finger pulse oximeter
SpO2 with LCD lighting
With color display and automatic switch-off function. For all standard applications. Includes battery, bag and instruction manual.

Piece	900 6368	97.50
-------	----------	-------



Set included
Cuff only
49.90
PRICE-PERFORMANCE

2 cuff sizes for exact measurement results
Extra large
Color display

Digital blood pressure monitor
Blood pressure monitor with an extra large display. 120 memory locations, average values of the last 3 measurements. Date and time display, alarm clock and reminder function, arrhythmia detection, blood pressure rating scale according to WHO for accurate assessment of the measured values. Mains or battery operation possible with 4 AA batteries. Incl. adult cuff for normal arm circumference (22-32 cm), storage tray and power supply. Order batteries separately.


	Order no.	Price
Set including cuff	900 8181	49.90
Cuff normal arms		
22-32 cm arm circumference	900 8182	14.90
Cuff strong arms		
up to 42 cm arm circumference	900 8183	15.90



Piece off
23.90
PRICE-PERFORMANCE

Ear thermometer
Measurement result after a few seconds. With fever alarm, memory function and battery level indicator. The delivery includes a battery and 15 disposable measuring covers.

Piece	900 6438	28.90
from 10 pieces each		23.90




Piece off
3.99
PRICE-PERFORMANCE

clinical thermometer
Digital clinical thermometer for quick measurements within just 30 seconds.
> Fever alarm
> Waterproof
> Memory function
Delivery including battery.

Piece	900 6437	5.99
from 5 pieces each		4.99
from 10 pieces each		3.99



24 pieces only
21.90
PRICE-PERFORMANCE



10 pieces from
13.90
PRICE-PERFORMANCE



Piece off
4.99
PRICE-PERFORMANCE



Pack off
9.90
PRICE-PERFORMANCE



Pack off
6.65
PRICE-PERFORMANCE

1 fast cold bag
Activation occurs by pressing the middle together and provides the desired cold for around 30 minutes. Skin-friendly covering.

15x19cm		
Pack (24 pieces)	900 8442	21.90

2 cold/warm compresses
Reusable cold/warm compresses provide effective relief from muscle pain and help reduce swelling. Flexible at low temperatures.

13x14cm		
Pack (10 pieces)	900 8294	13.90
12x29cm		
Pack (10 pieces)	900 8295	21.90

3 softgel cold/warm compresses
Like cold/warm compresses but with soft gel and a soft surface.

13x26.5cm		
Pack (10 pieces)	900 8296	33.90
from 5 packs each		29.90

4 sports cooling spray
For cooling and ice treatment of sports injuries without open wounds, e.g. strains, cramps, bruises, muscle pain, sprains and tennis elbow (acute complaints).

200 ml spray can	900 2546	5.99
from 10 pieces each		4.99

ECG electrodes for long-term ECG
ECG electrodes, round with patch and centered push button. Prepared with contact gel for long-term and stress ECG. Self-adhesive.
Pack (50 pieces)

Size	Order no.	Price
Ø50mm		
for adults	900 6699	12.90
from 5 packs each		11.90
from 10 packs each		9.90

ECG one time
Electrodes for resting ECG
Round with liquid gel foam carrier and centered push-button connection, self-adhesive, for resting ECG. Pack (50 pieces)

Size	Order no.	Price
Ø 55mm for adults		
Pack	900 6696	7.99
Ø 38mm for teenagers		
pack	900 6698	7.99
Ø 32mm for children		
Pack	900 6697	7.99
from 5 packs per variety,		7.79
from 10 packs per variety		6.65

50 pieces from
49.90
PRICE-PERFORMANCE



Oxygen goggles
Latex free. Attachment to the ears.
2 meters of break-resistant hose, medical grade PVC.
Pack (50 pieces)

1 oxygen goggle for adults		
curved approach	900 7036	69.90
2 oxygen goggles for adults		
Standard approach	900 7035	69.90
3 oxygen goggles for adults		
Funnel approach	900 7037	69.90
4 oxygen goggles for children		
Funnel approach	900 7034	69.90
from 5 packs per variety,		59.90
from 10 packs per variety		49.90



Cooler spray
Ozone-friendly cold spray, CFC-free. Exit temperature: approx. -40 °C. Propane/butane mixture. The spray head has an extended spray tube of 52 mm.

Fragrance	Crowd	Order no.	Price
Mint	200ml	900 1383	5.99
Orange	200ml	900 1384	5.99
from 5 pieces per variety			4.99

**Laboratory water**

Not sterile; fully desalinated; laboratory tested quality. Corresponds to VDE 0510 and DIN 57510.

Size	Order no.	Price
4 x 5 liter canisters	988 2042	19.90
10 liter canister	900 6280	6.99
from 3 pieces each		6.49
Accessories for 900 2577 and 900 6280		
tap	979 4136	9.99
Canister key	979 4140	9.99

**Electrode gel**

Extremely high electrical conductivity. Remains moisture-stable for a long time, is not greasy and is easy to wipe off.

> Absolutely bubble-free

Size	Order no.	Price
250ml bottle	900 0847	3.49
from 10 pieces each		2.99
5 liter cubic container	900 0848	19.90
from 3 pieces each		18.90
Tap for 900 0848		
Piece	979 4135	9.90

Protective covers for
ultrasound probes
without lubricant

Protective cover to prevent contamination of the ultrasound probes by body fluids. Makes cleaning and disinfection of the probe easier.

Dimensions: 34 x 195 mm (Ø x L)

Pack (144 pieces)

pack	900 0639	10.50
------	----------	-------



5 kg storage
container

12.90

PRICE-PERFORMANCE

Microbubble free

**Ultrasound gel**

- > Absolutely microbubble-free
- > Harmless to transducers
- > Skin and mucous membrane compatible
- > Without perfume or color additives
- > Free from stabilizers and formaldehyde

Size	Order no.	Price
250ml bottle		
Piece	900 0849	2.15
from 5 pieces		1.95
per 1,000 g dispenser		
Piece	900 0850	3.99
5kg canister		
Piece	900 0881	14.99
from 3 pieces		12.99
per 10 kg canister		
Piece	900 0882	25.90
from 3 pieces each		21.90
Dosing pump for 5 kg canisters		
Piece	979 4137	19.90
Dosing pump for 10 kg canisters		
Piece	979 4138	19.90
Canister key for 900 0881 and 900 0882		
Piece	979 4139	9.90
5kg storage container		
Piece	900 0851	14.90
from 4 pieces each		12.90
Tap for 900 0851		
Piece	979 4135	9.90



100% quality

Factory fresh goods

X-ray films PAN-G

> 100% quality meets the requirements of the European Medical Devices Act (MPG) and the X-ray Ordinance (RÖV).

> Each film contains 100 sheets

Green-sensitized

Size	Order no.	Price
18x24cm	900 8175	78.60
24x30cm	900 8176	92.50
30x40cm	900 8177	145.00

**Acclean CHX**

Mouthwash solution without alcohol

Alcohol-free mouthwash solution is used for professional cleaning and care of the mouth and throat area. The mouthwash solution has a bactericidal effect, creates fresh breath and reduces the formation of new plaque.

Acclean CHX prevents tooth decay, gingivitis and periodontitis and can also be used on inflamed and irritated gums. Can be used as a preventive measure before any dental treatment.

Size	Order no.	Price
250ml bottle	900 9663	5.90
1 liter bottle	900 9664	13.90
5 liter canister	900 9665	35.90

**Radiation protection apron**

Frontal

The apron offers frontal protection. With Velcro fastener. Without pockets. Pb equivalent 0.5 mm.

Colour blue

Size M: 65x120cm	Piece	988 1834	199.00
------------------	-------	----------	--------



Protecting the planet for a healthier future



Silk dental floss

Raft made from natural silk. Environmentally friendly packaged in a cardboard box with a stainless steel cutter. Coil tightly rolled and without additional plastic core.

25 m	Piece	979 6823	4.70
------	-------	----------	------



Acclean®

Bamboo toothbrush

Bamboo toothbrush with white filaments. Sustainable dental care. Handle made of environmentally friendly bamboo wood, sustainable production of the filaments (MedexS®) Green line filaments: sustainable filaments, consists of up to 100% renewable raw materials. There is no longer any need to compromise between sustainability and performance! Made from polyamide 10.10 or 6.10 from sustainable bio-based resources: PA 10.10 => 100% bio-based PA 6.10 => bio-based), vegan. Rounded brush head. Hardness medium.

Piece	979 6822	2.95
-------	----------	------



Toothbrushes

Optimal cleaning thanks to the conically tapered brush head. Handy, non-slip handle.

Pack (12 pieces)	714 0480	9.59
------------------	----------	------



Baby and children's toothbrush

Specially made for baby teeth, suitable for children up to 6 years. Attractive design with anti-slip handle, consisting of two components, with practical suction cup stand. Angled brush head with special contour. 27 filaments. Color range: pink/violet, yellow/red, orange/blue, cedro/green (no color selection possible).

Piece	988 3611	1.40
from 10 pieces each		1.29



Youth toothbrush

Modern toothbrush with rubberized thumb stop for children up to 12 years. 32 filaments. With suction cup for a secure stand on the sink. Compact brush head with a shorter handle (17 cm) specially designed for young people, consisting of two components. Color range: pink/white, black/orange, blue/green, violet/pink (no color selection possible).

Piece	988 3612	1.90
from 10 pieces each		1.69



Disposable toothbrush with

toothpaste

Individually packaged disposable toothbrushes with toothpaste. The brushes are ideal for cleaning teeth before dental treatment. Assorted by color.


100 pieces	988 0707	15.80
from 5 packs each		13.80



Travel toothbrush

Unique mechanism. The housing protects the brush head in your luggage while traveling. With angled neck for less jaw strain. Easy to use, the lid cannot be lost during travel. Space saving. To use: 1. Press, 2. Turn, 3. Lock. Color range: blue, red, green, turquoise (no color selection possible).

Piece	988 3613	1.50
from 10 pieces each		1.39



Patient cape

Transparent patient gown on a roll made of PE film with a neckline and straps, suitable for standard patient gown roll dispensers. 1-ply. Roll (200 pieces)

Size	Order no.	Price
white	900 0787	21.90
dark blue	900 2949	21.90
orange	900 2947	21.90
blue	900 0788	21.90
green	900 0789	21.90
from 5 pieces per variety		19.90



1 urine bottle for men with lid

1 liter, graduated pieces	988 0896	2.99
from 10 pieces each		1.99

2 urine bottle holders

Piece	988 0897	4.79
-------	----------	------



Urine bag holder

Hanging device for urine bags.

Piece	988 0898	1.90
-------	----------	------



napkins

Two-layer napkins made of absorbent cellulose fleece and a layer of moisture-impermeable plastic film (PE) for particularly reliable protection against soiling of clothing. Also suitable as a place mat. Dimensions: 32 x47 cm. Pack (500 pieces)

Size	Order no.	Price
white	900 1030	19.90
Light Blue	900 1031	19.90
dark blue	900 1402	19.90
green	900 1032	19.90
turquoise	900 1034	19.90
yellow	900 1039	19.90
purple	900 1038	19.90
peach	900 1036	19.90
pink	900 1033	19.90
lime green	988 2021	19.90
orange	900 2940	19.90
turquoise	988 3144	19.90
from 5 packs per variety		17.90

piece only

4.99
PRICE-PERFORMANCE



Napkin holder

Chrome plated with ball chain. Total length 40cm.


Piece	900 1746	4.99
-------	----------	------



Medicine dispenser

Medicine dispenser, 4 compartments (morning, lunch, evening, night), with transparent sliding lid.

Size	Order no.	Price
white	900 9635	0.89
blue	900 9636	0.89
from 10 pieces per variety,		0.79
from 20 pieces per variety		0.59



SERVICE	
About Henry Schein	69
team	70-71
Webshop	72
Newsletter	73
Practice Green	74-75
Henry Schein own brand	76-77
Technical service/repairs/ Maintenance	78
Device checks	79
Probe analyses	80
Leasing/Financing	81
Practice planning	82
All from a single source	83



WE ARE
HENRY SCHEIN

The story of Henry Schein



About Henry Schein Inc.

Henry Schein, Inc. (Nasdaq: HSIC) offers high-quality solutions for physicians and healthcare professionals through a broad network of people and technologies.

With solutions for operational, clinical, technological and supply chain areas support all team members relaxed dentists and general practitioners, so that these

their high quality health care offer even more effectively can.



Henry was founded Semblance of Henry and Esther Semblance in 1932: She had a small counter pharmacy in Queens, NY.



Use the strength we give you
Offer and get advice

- > Operations and branches in 32 countries and regions
- > Serving over 1 million customers
- > Approximately 2,100 telesales sales employees
- >90 years in business
- > More than 22,000 Team Sham members
- > 300,000 products in stock

- > Shipping from avg. 156,000 boxes per day
- > 48 distribution and production sites
- > More than 3,200 partner suppliers
- > More than 2,100 field service technicians
- > 180,000 additional products as special items
- > More than 3,400 sales consultants in the field
- > 99% of the items ordered will be delivered on Shipped on the business day the order is received



THE HOT WIRE - A STRONG TEAM AWAITS YOU.

If you have any questions about your order, delivery or other concerns, simply contact our customer service. We look forward to your call and will be happy to help.



Thomas Friedl
Managing Director

thomas.friedl@henryschein.com
Telephone: 01 / 718 19-61



Alexandra Takacs
Marketing &
Administration Manager
alexandra.takacs@henryschein.com
Telephone: 0664 / 457 31 78



Sandra Rausch
Executive Assistant &
Customer Service Specialist
sandra.rausch@henryschein.at
Telephone: 01 / 718 19-15



Sandra Kieseewetter
Customer Service
sandra.kieseewetter@henryschein.at
Telephone: 01 / 718 19 61-25



Melanie Prokop
Customer Service
melanie.prokop@henryschein.at
Telephone: 01 / 718 19 61-10



Yvonne Geyer
Customer Service
yvonne.geyer@henryschein.at
Telephone: 01 / 718 19 61-13

Always there for you personally!
Telephone: 01 / 718 19 61 - 99 Monday - Thursday 8:00 a.m. - 5:00 p.m., Friday 8:00 a.m. - 1:00 p.m. Email: bestellung@henryscheinmed.at Postal
address: Henry Schein Medical Austria GmbH, Schönbrunner Straße 297, 1120 Vienna

AT YOUR LOCAL — GET PERSONAL ADVICE.

Simply call and arrange a personal consultation appointment in your practice. The map shows you which area representative looks after your location. We will be happy to assist you with advice and assistance.



Wolfgang Vrablik
sales Manager
wolfgang.vrablik@henryschein.at
0664 / 822 82 67

Areas:
Vienna (105x, 110x-118x, 123x, 13xx)
Lower Austria (23xx-28xx)
Burgenland (70xx-72xx)



Christian Wippl
External sales
christian.wippl@henryschein.at
0664 / 857 17 60

Areas:
Burgenland (73xx-75xx)
Styria
Carinthia




Ivan Filipovic-Cugura
External sales
ivan.filipovic-cugura@henryschein.at
0664 / 857 17 89

Areas:
Vienna (100x, 101x-104x, 106x-109x, 119x, 122x, 14xx)
Lower Austria (20xx-22xx, 30xx-32xx, 34xx-37xx)



Andreas Schwaiger
External sales
andreas.schwaiger@henryschein.at
0664 / 883 31 278

Areas:
Lower Austria (33xx, 38xx-39xx)
Upper Austria
Salzburg (50xx-57xx)
Tyrol Vorarlberg



Areas:

- | | | |
|-----------------|------------------------|-------------------------|
| 1 Vienna | 4 Upper Austria | 7 Burgenland |
| 2 Lower Austria | 5 Salzburg | 8 Styria |
| 3 Lower Austria | 6 Tyrol and Vorarlberg | 9 Carinthia, East Tyrol |



WWW.HENRYSCHMED.AT AN ADDRESS WITH MANY ADVANTAGES!

YOUR ONLINE ADVANTAGES:

- > Access to over 10,000 articles
- > Everything at a glance with your customer account
(Invoice, order history, shipment tracking)
- > Order history can be viewed at any time
- > Delivery tracking
- > View invoice as PDF
- > Create watch lists
- > No waiting time on the phone
- > Product availability directly visible
- > Register returns & contact the person responsible
Find field staff
- > Online weekly promotions
- > Up to 10% discount on our Sun Schein Days
- > Online special offers



www.henryscheinmed.at

NEWSLETTER - DON'T MISS

ANYTHING FROM NOW ON!

With the Henry Schein newsletter you will always be up to date on all the latest news and promotions on
www.henryscheinmed.at.


Look forward to:

- > Current news
- > Attractive offers and
exclusive online promotions
- > Vouchers and discount promotions

Register now and secure your benefits!

It's that easy:

At the end of the starter journey, simply enter your email
address at www.henryscheinmed.at and fill out your details
and click "Subscribe now".





mission

With Practice Green, Henry Schein aims to encourage people to become more environmentally friendly by offering solutions that support, educate and support our customers' efforts to improve ecological sustainability.

central message



Practice Green is intended to improve health enable beings to future of our planet to have a positive impact by reducing the ecological footprint and promoting sustainability - both in everyday life and at work.



Our goal is to support healthcare professionals to strengthen their sustainability efforts while continuing to provide high-quality patient care.



Henry Schein's customers you can rely on us to offer a selection of products, Services and solutions that help to make the practice or laboratory more environmentally friendly design.



We realized that we have to Our global environmental footprint and our unique position within an ecosystem of relationships with suppliers, business partners and customers can be a driving force for sustainability.



As a global supplier of Lö With solutions for the healthcare sector, Henry Schein supports its customers on their way to improving ecological sustainability.



As part of Practice Green Henry Schein offers his customers Solutions that help, practical To pursue paths to greater sustainability while at the same time complying with applicable health and safety regulations to comply with regulations and a to ensure high quality patient care.



We are reducing ours in the Environmental impact downstream by encouraging our customers to improve their practices to make it more environmentally friendly and to work with us towards sustainability.



This program is consistent with our commitment to our effects on the climate to reduce the amount of waste and to use valuable materials to recycle and reuse and to use energy and water resources intelligently and efficiently.



At Henry Schein, we strive to organize our activities in a way that protects the earth. We support our customers' efforts towards ecological sustainability. Together we can create a healthier world for everyone.

SUSTAINABLE PRODUCTS

BY HENRY SCHEIN



Practice Green

Do you already know the new, sustainable product line from Henry Schein? It includes both products made from renewable raw materials also biodegradable products. For example:



You can find many other sustainable products in our webshop

<https://henryscheinmed.at/p/practice-green>



TOP QUALITY AT PERMANENTLY AFFORDABLE
PRICES — THE HENRY SCHEIN OWN BRAND

Discover our Henry Schein own brand!
We are constantly adding new and, above all,
inexpensive products to our range.

You can recognize our own brand products by the
quality seal, which stands for the best selection,
best quality and consistently low prices.
Fulfilling your customer wishes and even more service
are our top priority.

Many doctors rely on our own brand products in
their daily work.
When are you going to put our own brand to the test?



Five advantages
the Henry Sham
Own brand:

1. Unbeatable value for money
2. Highest quality - "Seal of Excellence"
3. Customer-friendly return policy
4. Diversity - over 10,000 products worldwide
5. Always up to date

1.

Unbeatable
Price-performance ratio

Our dedicated private label team stands worldwide
widely in contact with various producers.

As soon as we see a product that seems relevant
discover, we check according to our strict standards
Quality criteria as to whether it is an option for us.

How do the better prices for own brands
come about?

With over 1 million customers worldwide, Henry
Schein has an enormous purchasing volume.
This reduces purchasing costs and ours

Private label products are offered to you at reasonable
prices.

2.

Highest quality -
"Seal of Excellence"

Our goal is to provide the customer with all our
Products to always deliver the highest quality.
This is what we do as a trading company
Market experience from over 90 years ago.
We learn from our daily exchanges with customers
exactly what needs and requirements
our products must meet the requirements.
Based on these experiences, our
Quality requirements are defined, which are the reason
location for a product range of over 10,000
products worldwide

As a sign of our quality
Everyone will promise
Henry Schein products
Own brand with the "Seal
of Excellence".



4.

Diversity with over
10,000 products worldwide

The Henry Schein own brand has...
10,000 products worldwide via one of the
largest product ranges in the industry.
Be inspired and discover a selection of Henry Schein
own-brand products in our online shop.



3.

Customer friendly
Return Policy

The "Seal of Excellence" stands for our company
trust in private label products. Therefore secure
The seal gives you the possibility of an expanded one
return to. Your responsible person will be happy to inform you
Field service advisor personally.

5.

Always on that
latest status

We develop the Henry product range
Appearance private label for you further. It works
in addition to the inclusion of new products in
our catalog, including more environmentally friendly ones
Alternatives for existing products. In our
We have all the latest products in our online shop
Henry Schein own brand for you

summarizes:



TECHNICAL SERVICE -
FAST, SIMPLE, AUSTRIA!

Call us and we will inform you about options: 01 / 718 19 61 - 99

We will inform you about

- > Repairs
- > Rental equipment
- > Metrological controls (MTK)
- > Safety controls (STK)
- > Installation of medical technology devices
- > Assembly/assembly of practice furniture
- > Maintenance, etc



EQUIPMENT TESTS AS SCHEDULED —
WE WILL HELP YOU!



The use and use of medical devices in Austria is regulated by the Medical Devices Act (MPG) and the Medical Device Operator Ordinance (MPBV). A distinction is made between:

Metrological

Controls (MTK)

Metrological controls are for all devices with measuring function, which are listed in Appendix 2 of the MPBV, mandatory:

- > 1 year review period for medical devices acoustic determination of hearing ability (e.g. audiometer, infrared radiation thermometer)
- > Inspection period 2 years for electric thermometers, Medical devices with replaceable temperature sensors, medical devices for non-invasive Blood pressure measurement, medical devices for determining intraocular pressure (e.g. blood pressure monitor, eye tonometer)
- > Verification period according to the Measurement and Calibration Act medical scales

Safety-related

Controls (STK)

Safety controls are for operators

Periodic safety checks required for medical devices that must be carried out.

According to the Medical Device Operator Ordinance (MPBV), the following devices must undergo STK:

- > All active medical devices, all those products that are operated with energy (electricity, compressed air, etc.), including defibrillators, stimulation current devices, ECG devices, sterilizers, ultrasound devices, electro-surgical devices and much more
- > All non-active medical devices where the Manufacturer specifies an STK

Interval: If the manufacturer specifies an interval, this must be adhered to. If the manufacturer does not explicitly specify a STK or an interval, the intervals according to the MPBV must be adhered to. The test intervals are between 1 and a maximum of 3 years.

REPAIRS AND MAINTENANCE

Qualified device management and reliable service are indispensable foundations on which you can rely you need to rely on when testing and maintaining your medical devices. For all questions about maintenance and repairs, simply contact our technical service - we will take care of it.

Your advantages:

- > Shortest response times
- > Affordable rental equipment
- > Reduced practice downtime
- > Longer life cycles for your investments
- > Maintaining the value of your practice inventory
- > Affordable maintenance flat rates
- > Compliance with all legal requirements



We would be happy to monitor the testing deadlines for your medical-technical devices. With our "reminder service" we ensure that you never miss a test appointment. Contact our technical service!



PROBLEMS WITH YOUR PROBE?
OUR TECHNICAL SERVICE IS HERE FOR YOU!

ServiceDirect

Reliably working probes are a basic requirement for your work. Our technical service supports you in all matters. Use our Henry Schein ServiceDirect.

- > Complete probe analysis
- > Probe preparation/repair (test certificate, Final inspection, documentation, guarantee)
- > Available rental probes
- > Exchange probes at attractive conditions
60% savings compared to buying new (determined EPRN average value), independent of manufacturer, Certified according to ISO 13485 + ISO 9001



OUR SERVICES IN DETAIL:

- | | |
|----------------------------------|------------------------|
| > Cable & connector repair | > Probe preparation |
| > Renewal of the kink protection | > Fix crystal failures |
| > Renewal of the coating | > TEE probe repair |
| > Probe analysis | |

Maximum performance at low costs

A probe analysis detects hidden problems in good time. This means that repairs or reprocessing can be carried out before the analysis quality is lost or the probe becomes irreparable.

INDIVIDUAL LEASING OFFERS —
INVEST IN YOUR SUCCESS!

Our leasing service

Regardless of whether you want to finance a device or your entire practice facility through leasing, We will make you an individual offer according to your wishes. Starting from a net worth of You have the opportunity to win 1,000 euros Use leasing service!



Your advantages:

- > Service billing without investment
- > Increased liquidity because your equity capital is not tied up
- > Free choice of leasing property
- > Transparent calculation
- > Low monthly leasing rates at attractive conditions
- > Leasing fee can be deducted as a business expense
- > Flexibility - Individual monthly rates
- > You can choose between 30, 36, 42, 48, 54, 60 or 72 small, comfortable monthly installments!
- > Planning security with consistent monthly leasing rates
- > Business and tax advantages
- > Faster copying options

FINANCING CONCEPTS —
GREATER FINANCIAL ROOM!

Quick and easy:

- > We offer you tailor-made financing concepts for the new establishment of your practice for individual modernization measures.
- > Depending on the task and taking into account To determine the future profitability of your medical practice, we will create a suitable leasing concept for you.
- > After you have decided on your leasing option If you have divorced, we will take care of the entire process. Our competent contacts will be happy to assist you by telephone.

Just call us: 01 / 718 19 61 - 99





PRACTICE FACILITY

WITH HENRY SCHEIN SERVICE DIRECT

A holistic planning concept offers this Orientation and safety for patients and the entire practice team. Furthermore, it ensures a relaxing stay. Your practice will be unique and distinctive with us as your full service Specialists.

Use our Henry Schein ServiceDirect.

- > Analysis of practice/room conditions
- > Including consideration of investment volume, Design and color requests etc.
- > Individual 3D planning
- > Implementation
- > Organization, delivery and construction of the new ones Practice facility



INTERIOR

Medical practices are workspaces that must offer technical perfection, a pleasant climate for patients and ergonomic work for doctors. We have specialized in setting up ordinations for years. With our regional and international partners, we implement your desired solution for you.

Submit work to us

One contact or twenty?

We are happy to help you with everything from practice planning to basic supplies with consumables. With our cooperation partners we offer a wide range of services. Use us as a nerve-saving interface.

No matter whether you have questions, would like detailed and competent advice or would like to order, we are there for you.

We are
there for her!

ALL FROM A SINGLE SOURCE -

SELECTION, ADVICE, SERVICE!

Over 10,000
products

... with us you will find everything you need
successful ordination!



Many well-known ones
Manufacturer

... are with new products and
proven classics represented.



MELAG

ZOLL



HAEBERLE



BOSCH
+SOHN

boso

seca

Präzision für die Gesundheit

DR. MACH

Miele

CARDIOLINE

Vitalograph

the Speichertechnik Partner

PANNOMED

Affordable
alternative

... benefit from our alternative
to branded products.



All prices in € plus the statutory. VAT Subject to errors and price changes - the current price and our terms and conditions are always available online at henryscheinmed.at. Acceptance of the order is subject to the availability of the goods. Partial deliveries may occur.





2023/24

Own brands
range



Data protection information: You are receiving the information sent to you because you are a customer of Henry Schein Medical Austria GmbH or we purchased your address and contact details from Verlag Dieter Göschl GmbH, Hernals Hauptstraße 213, 1170 Vienna, Austria for the purpose of advertising have. You can object to the transmission of this advertising without giving reasons in the future to Henry Schein Medical Austria GmbH, Schönbrunner Straße 297, 1120 Vienna by letter, by email at info@henryscheinmed.at or by fax at 0800/7181961-23. Further information about data protection at Henry Schein can be found on our website henryscheinmed.at/p/datenschutz. Commercial register number 102456x, commercial register court: Vienna Commercial Court