



Product Name	Type	Color	Mil Thickness	Texturing	Extra Features
Criterion® Aloe Advanced	Exam Non-Sterile	Natural White	3.2 mil	Textured	Aloe & Vitamin E Coated
Criterion® Aloe Green	Exam Non-Sterile	Green	5.1 mil	Textured	Aloe Vera Coated
Criterion® CL	Exam Non-Sterile	Natural	4.3 mil	Textured	
Criterion® CL2	Exam Non-Sterile	Natural White	3.14 mil	Textured Fingers	
Criterion® PC	Exam Non-Sterile	White	3.2 mil	Textured	Polymer Coated
Criterion®	Surgical Sterile	Latex	7.1 mil	Microtextured	Polymer Coated
Criterion®	Exam Non-Sterile	Blue	3.6 mil	Textured Fingers	
Criterion® COATS™	Exam Non-Sterile	White	2.8 mil	Textured	Colloidal Oatmeal Coated
Criterion® N200	Exam Non-Sterile	Blue	2.4 mil	Textured	
Criterion® N250	Exam Non-Sterile	Lilac	3.1 mil	Textured	
Criterion® N300	Exam Non-Sterile	Ice Blue	2 mil	Textured Fingers	
Criterion® Nitrile EC	Exam Non-Sterile	Periwinkle	4.7 mil	Textured	Extended Cuff
Criterion® Pure Freedom	Exam Non-Sterile	Slate Blue	2 mil	Textured Fingers	Accelerator Free
Criterion® Radiant Razz	Exam Non-Sterile	Radiant Razz	2.4 mil	Textured Fingers	Aloe & Vitamin E Coated
Criterion® CR	Exam Non-Sterile	Green	4 mil	Textured	Polymer Coated
Criterion® CR 4G	Exam Non-Sterile	Arctic Lime	2.8 mil	Textured Fingers	
Criterion® CR	Surgical Sterile	White	5.9 mil	Textured	Polymer Coated
Criterion®	Exam Non-Sterile	Clear	3.15 mil	Smooth	

For our **HENRY SCHEIN BRAND PRODUCTS CATALOG**, as well as additional product information, visit WWW.HENRYSCHEINDENTAL.COM/HSBRAND.



**VALUE AND
QUALITY**



**PRODUCT
SELECTION**



**RELIABLE AND
TRUSTWORTHY**



HENRY SCHEIN®

CRITERION® Brand Gloves

The Henry Schein Criterion® line of gloves is designed to meet all your dental practice needs, whether exam or surgical, latex or non-latex. All are backed by the Henry Schein Seal of Excellence.

LATEX POWDER-FREE EXAM

NON-STERILE



CRITERION® ALOE ADVANCED

Latex Powder-Free Aloe Exam
Natural white latex color. Coated with an enhanced compound of aloe and vitamin E. Textured and nonchlorinated.

Gloves..... 100/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:
XS (112-5840)
S (112-5841)
M (112-5842)
L (112-5843)
XL** (112-5844)

**90/Box



CRITERION® ALOE GREEN

Latex Powder-Free Aloe Exam
Inside of glove is coated with a layer of aloe to soothe and moisturize hands. Textured and chlorinated. Green.

Gloves..... 100/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:
XS (108-7848)
S (108-7850)
M (108-8553)
L (108-7855)
XL** (108-7856)

**90/Box



CRITERION® CL

Latex Powder-Free Chlorinated Exam
Comfortable fit and low protein levels. Textured and chlorinated.

Gloves..... 100/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:
XS (102-5354)
S (102-8296)
M (102-5355)
L (102-7763)
XL (102-5357)



CRITERION® PC

Latex Powder-Free Polymer-Coated Exam
Synthetic polymer coating eases donning. Textured finish for enhanced grip. Excellent tactile sensitivity. Nonchlorinated.

Gloves..... 100/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:
XS (102-5418)
S (102-5419)
M (102-5421)
L (102-6730)
XL (102-5422)

LATEX POWDER-FREE SURGICAL

STERILE



CRITERION®

Sterile Latex Powder-Free Surgeons
Powder-free to help eliminate powder-related complications and irritation. Advanced polymer coating on the inside surface allows for easy donning when hands are wet or damp. Microtextured surface and optimized tack level provide excellent grip in wet or dry conditions. Supersoft latex film and anatomical design provides exceptional comfort to reduce hand fatigue even during long procedures.

Gloves..... 50 Pair/Box
4 Boxes (Case) @

Specify:
5½ (570-0635) 7½ (570-0639)
6 (570-0636) 8 (570-0640)
6½ (570-0637) 8½ (570-0641)
7 (570-0638) 9 (570-0642)

NITRILE NON-LATEX EXAM

NON-STERILE



CRITERION®

Nitrile Powder-Free Exam
Soft and stretchy for better comfort and fit. Fully textured and chlorinated. 3.6-mil thick. Blue.

Gloves..... 100/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:
XS (111-8533)
S (111-8535)
M (111-8536)
L (111-8537)
XL** (111-8539)

**90/Box

*of the same size

All gloves are gluten-free.

NITRILE NON-LATEX EXAM (Cont'd)

**200
PER BOX**



CRITERION® COATS™

Nitrile Powder-Free

Through a combination of patented and proprietary processes, an all-natural colloidal oatmeal USP coating has been added to our barrier protective glove. Provides continuous transfer and delivery of unique skin protection during typical glove use. Gluten-free. 2.8-mil thick.

Gloves..... 200/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:

XS	(112-6220)
S	(112-6221)
M	(112-6222)
L	(112-6223)
XL	(112-6224)

**200
PER BOX**



CRITERION® N200

Nitrile Powder-Free Exam

Innovative formulation emulates the fit and feel of latex with added strength and puncture resistance. "Real-feel" grip for optimal tactile sensitivity. Textured fingers provide enhanced, superior grip in both wet and dry applications. 2.4-mil thick. Blue.

Gloves..... 200/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:

XS	(900-7437)
S	(900-7438)
M	(900-7439)
L	(900-7440)
XL**	(900-7441)

**180/Box

**250
PER BOX**



CRITERION® N250

Nitrile Powder-Free Exam

Thin nitrile glove that provides added protection and delivers greater comfort and increased tactile sensitivity. Textured fingertips for enhanced grip in a dental setting. Not made from natural latex to help protect practitioners from skin allergies. Economical packaging: 60% less storage space and 50% less packaging waste. Chemotherapy tested. 3.1-mil thick. Lilac.

Gloves..... 250/Box
10 Boxes (Full Case)* @

Specify:

XS	(570-0099)	L	(570-0102)
S	(570-0100)	XL	(570-0103)
M	(570-0101)		

**300
PER BOX**



CRITERION® N300 ULTRA

Nitrile Powder-Free Exam

Ultrathin nitrile gloves that are so comfortable they feel as if you are not wearing gloves at all. Get the protection you need without sacrificing the comfort you want. Natural feel and formfitting. Textured fingertips. Eco-friendly packaging, with 300 gloves in a box, reduces waste and saves space. 2.0-mil thick. Ice Blue.

Gloves..... 300/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:

XS	(112-6829)
S	(112-6830)
M	(112-6831)
L	(112-6832)
XL**	(112-6833)

**180/Box



CRITERION® NITRILE EC

Nitrile Powder-Free Exam, Extended Cuff

12" exam gloves that offer extra protection for high-risk situations, providing protection between the cuff and sleeve. Textured for a secure grip. Offer a high level of nitrile barrier protection with comfort and tactile sensitivity similar to latex. Tested and approved for use with chemotherapy drugs, formaldehyde, gluteraldehyde, and chlorhexidine gluconate. 4.7-mil thick. Periwinkle.

Gloves..... 50/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:

XS	(112-6785)
S	(112-6786)
M	(112-6787)
L	(112-6788)
XL	(112-6790)

**200
PER BOX**



CRITERION® PURE FREEDOM

Nitrile Powder-Free Exam

Powder-free nitrile gloves that are manufactured without the use of chemical accelerators or sulfur, preventing type-IV allergic reactions. They reduce the potential for dermatitis for glove wearers while providing softness, comfort, and maintaining superior protection. Textured fingertips for improved grip. 2.0-mil thick. Slate Blue.

Gloves..... 200/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:

XS	(570-0630)	L	(570-0633)
S	(570-0631)	XL**	(570-0634)
M	(570-0632)		

**180/Box

**200
PER BOX**



CRITERION® RADIANT RAZZ

Nitrile Powder-Free Exam

Extra-soft stretch nitrile glove with soothing aloe and vitamin E lining that provides an enhanced coating to rejuvenate and nourish hands. Thin feel fits like latex and provides improved tactile sensitivity. Textured fingers provide excellent grip in a wet/dry dental setting. "Smart Dispense" technology allows a single-glove delivery with fewer wrinkles. 2.4-mil thick. Radiant Razz color.

Gloves..... 200/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:

XS	(112-7143)	M	(112-7146)
S	(112-7144)	L	(112-7147)

*of the same size

All gloves are gluten-free.

*of the same size

INFECTION CONTROL

CHLOROPRENE NON-LATEX EXAM

NON-STERILE



CRITERION® CR

Chloroprene Powder-Free Exam

Offers the feel and elasticity of a latex glove without the worry of an allergic reaction to latex or nitrile. More natural feel and fit to hand. Higher tensile strength. 4.0-mil thick. Green.

Gloves..... 100/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:

XS (112-6857)
S (112-6858)
M (112-6859)
L (112-6860)
XL (112-6861)



CRITERION® CR 4G

Chloroprene Powder-Free Exam

Third-generation chloroprene formulation, which provides improved comfort and fit, in a thinner construction. Finger texture delivers improved tactile sensitivity. Improved modulus provides a softer, more flexible glove, which results in a snug but comfortable fit that can improve performance, increase dexterity, and reduce hand fatigue. 2.8-mil thick. Arctic Lime color.

Gloves..... 100/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:

XS (112-7306)
S (112-7307)
M (112-7308)
L (112-7309)
XL (112-7310)

SYNTHETIC NON-LATEX SURGICAL

STERILE



CRITERION® CR

Sterile Polychloroprene Powder-Free Surgeons

For the latex-sensitive user, this glove offers a similar fit and feel to latex. Excellent tactile sensitivity with a textured surface for enhanced gripping. Unique formulation allows gloves to be worn for extended periods without hand fatigue. Ideal for all oral surgery, perio, and implant procedures. White.

Gloves..... 50 Pair/Box
4 Boxes (Case) @

Specify:

5½ (112-7079)
6 (112-7080)
6½ (112-7081)
7 (112-7082)
7½ (112-7083)
8 (112-7085)
8½ (112-7086)
9 (112-7088)

VINYL NON-LATEX EXAM

NON-STERILE



CRITERION®

Vinyl Powder-Free Exam

Excellent solution for those at risk of latex sensitivity. Contains no natural latex. Supple touch for maximum tactile sensitivity. Smooth and nonchlorinated. Clear.

Gloves..... 100/Box
10 Boxes (Mix & Match) @
10 Boxes (Full Case)* @

Specify:

XS (570-0087)
S (570-0088)
M (570-0089)
L (570-0090)
XL** (570-0091)

**90/Box

Glove-Sizing Assistant

X-SMALL/ Sterile: 5–5.5


SMALL/ Sterile: 6–6.5

MEDIUM/ Sterile: 7–7.5

LARGE/ Sterile: 8–8.5

X-LARGE/ Sterile: 9

Place the widest part of your hand over the size grid as shown



*of the same size

All gloves are gluten-free.

INFECTION CONTROL

GLOVE DISPENSERS

MAXI-DISPENZ® WIRE GLOVE BOX DISPENSERS

Durable steel wire frame with white-powdered finish can be mounted on wall for easy dispensing of gloves. Selection of 5 different choices meets any office needs. Hardware is included.

A. Dual Box Dispenser

Holds 1 box horizontally or 2 boxes vertically.
10 $\frac{3}{4}$ "W x 4"H x 3 $\frac{1}{16}$ "D
(112-5471)..... Ea.

B. Double Horizontal Box Dispenser

Holds 2 glove boxes horizontally.
11"W x 11 $\frac{1}{2}$ "H x 4"D
(112-5473)..... Ea.

E. Double Vertical Box Dispenser

Holds 2 glove boxes vertically.
11"W x 11 $\frac{1}{2}$ "H x 4"D
(112-5472)..... Ea.

D. Triple Horizontal Box Dispenser

Holds 3 glove boxes horizontally.
16"W x 8"H x 4"D
(112-5475)..... Ea.

C. Triple Vertical Box Dispenser

Holds 3 glove boxes vertically.
11"W x 14"H x 4"D
(112-5474)..... Ea.



MAXI-DISPENZ® ACRYLIC GLOVE BOX DISPENSERS

Clear acrylic glove box holder. Available in 4 sizes to meet any office or facility's needs. Accommodates most brands and sizes. Option to wall mount with predrilled holes (screws NOT included) or place on countertop.

A. Single Box Holder

10 $\frac{1}{16}$ "W x 4 $\frac{3}{4}$ "H x 3 $\frac{1}{2}$ "D
(102-6467)..... Ea.

B. Double Box Holder


12"W x 10"H x 3 $\frac{3}{4}$ "D
(102-6468)..... Ea.

C. Triple Box Holder

11"W x 14 $\frac{1}{2}$ "H x 3 $\frac{1}{2}$ "D
(102-8344)..... Ea.

D. Quadruple Box Holder

11"W x 19 $\frac{1}{4}$ "H x 3 $\frac{1}{2}$ "D
(102-7201)..... Ea.



INFECTION CONTROL

ASTM LEVEL 1

- Oral exams for low-risk patients
- Procedures involving low levels of spray and splatter
- Orthodontic treatment with little or no spray
- Hand instrumentation with little moisture or spray
- Fluoride treatments
- Basic lab work

ASTM LEVEL 2

- Oral exams for known higher-risk patients
- Procedures involving moderate levels of spray, splatter and airborne particles
- High or low-speed handpieces
- Dry, low-aerosol procedures
- Bleaching, orthodontic, endodontic and periodontic procedures
- Lab work with exposure to blood or contaminated materials

ASTM LEVEL 3

- Oral exams for known high risk patients (AIDS, hepatitis B, tuberculosis)
- Procedures involving high levels of spray, splatter and airborne particles
- Ultrasonic scalers and high-speed aerosols
- Air polishers and air-water syringes
- Oral surgery (tooth extractions, maxo/facial reconstruction)

ASTM REQUIREMENTS	ASTM LEVEL 1	ASTM LEVEL 2	ASTM LEVEL 3
Bacterial filtration efficiency (BFE)	≥ 95%	≥ 98%	≥ 98%
Sub-micron particulate efficiency at 0.1 µm (PFE)	≥ 95%	≥ 98%	≥ 98%
Resistant to penetration by synthetic blood	@ 80 mm Hg	@ 120 mm Hg	@ 160 mm Hg
Differential Pressure (ΔP)	<4.0 mm H ₂ O/cm ²	<5.0 mm H ₂ O/cm ²	<5.0 mm H ₂ O/cm ²
Flame Spread	Class 1	Class 1	Class 1

EARLOOP MASKS

PROCEDURE MASKS LEVEL 1

Light, comfortable mask for extended wear. Three-ply, soft, nonirritating inner and outer facing.

- Fluid resistant
- Full-width, flexible nosepiece
- Fiberglass-free
- Not made with natural rubber latex

Earloop Masks..... 50 per Box
6 Boxes (Case) @

Specify:

Arctic Blue (570-0738) Pink (104-3730)
Blue (104-3809) Yellow (104-2849)
Fuchsia (570-0739) White (104-6611)
Lavender (104-8600)

Compare to ValuMax Ultra 3-in-1 Masks

PROCEDURE MASKS LEVEL 2

Spunbond outer layer with white inner facing for enhanced wearer comfort.

- Fluid resistant
- Extra-long aluminum nosepiece
- Made in the USA
- Not made with natural rubber latex

Earloop Masks..... 50 per Box
6 Boxes @

Specify:

Blue (104-7799) Pink (104-9200)
Lavender (104-8071)

MAXIMA® PROCEDURE MASKS LEVEL 3

Same excellent features as the Henry Schein Procedure Mask but with a higher filtration rate and fluid-penetration resistance.

- Fluid resistant
- Full-width, flexible nosepiece provides fit and security
- Fiberglass-free
- Not made with natural rubber latex

Earloop Masks..... 50 per Box
6 Boxes (Case) @

Specify:

Blue (570-1527) Teal (570-1525)
Lavender (570-1530) White (570-1526)
Pink (570-1528) Yellow (570-1529)

MAXI-GARD™ AIR PROCEDURE MASKS WITH AIR SPACE™ TECHNOLOGY LEVEL 1

Feature an innovative, contoured design that enhances air flow space, and provide complete coverage and fit with extra breathing room that allows for a more comfortable, extended wear.

- Fluid resistant
- Adjustable, full-width nosepiece
- Fiberglass-free
- Not made with natural rubber latex

Earloop Masks..... 50 per Box
6 Boxes (Case) @

Specify:

Blue (570-0722) Lavender (570-0724)
White (570-0723)

MAXI-GARD™ AF2 ANTI-FOG PROCEDURE MASKS LEVEL 2

Outer and inner spunbond layers with white, nonwoven antifog strip. Blue.

- Antifog strip for fog reduction and increased comfort
- Fluid-resistant
- Adjustable nosepiece
- Not made with natural rubber latex

Earloop Masks
(112-7252)..... 50 per Box
6 Boxes @

Compare to ValuMax Extra-Safe Sensitive Mask

MAXI-GARD™ SFT MASKS LEVEL 1

Outer and inner spunbond layers with a white inner layer for sensitive skin.

- Fluid resistant
- Extra-long aluminum nosepiece
- 3" aluminum chin piece for a tight fit
- Made in the USA
- Not made with natural rubber latex

Masks—Level 1..... 50 per Box
6 Boxes @

Specify:

Pink (900-7452) Blue (900-7451)

MAXI-GARD™ SFT MASKS LEVEL 2

Outer and inner spunbond layers with a white inner layer for sensitive skin.

- Fluid resistant
- Provides a higher level of protection
- Extra-long aluminum nosepiece
- 3" aluminum chin piece for a tight fit
- Made in the USA
- Not made with natural rubber latex

Earloop Masks..... 50 per Box
6 Boxes @

Specify:

Pink (900-7454) Blue (900-7453)

PROCEDURE PINK RIBBON MASKS LEVEL 1

Same excellent features as the Henry Schein Procedure Mask but with a pink-ribbon design to support breast cancer awareness. Not made with natural rubber latex.

Earloop Masks
(900-4210)..... 50 per Box
6 Boxes (Case) @

ULTRASOFT PROCEDURE MASKS LEVEL 1

Spunbond outer layer with a white inner facing.

- Fluid resistant
- Hypoallergenic
- Will not lint, tear, or shred
- Made in the USA
- Not made with natural rubber latex

Earloop Masks..... 50 per Box
6 Boxes @

Specify:

Blue (104-8072) Pink (104-8073)



TISSUE-LINED PROCEDURE MASKS LEVEL 1

Designed for individuals with sensitive skin.
Soft-tissue inner layer. White.

- Fluid resistant
- Full-width, flexible nosepiece
- Not made with natural rubber latex
- Dye-free

Earloop Masks

(104-4632).....50 per Box
6 Boxes (Case) @

SURGICAL TIE-ON MASKS



PLEATED SURGICAL MASKS LEVEL 1

Pleated masks with cloth-like ties for a better fit.

- Fiberglass-free
- Not made with natural rubber latex

(104-6809).....50 per Box
6 Boxes (Case) @

MOLDED MASKS



MOLDED FACE MASKS

Available in either regular size blue or 10% smaller pink. Not made with natural rubber latex.

Molded Masks

Office Pack..... 50 per Box
5 Boxes @

Specify:

Pink (100-8350) Blue (100-8356)

Molded Masks

Bulk Pack..... 100 per Box
5 Boxes @

Specify:

Pink (100-6723) Blue (100-2686)



ECONOMY MOLDED FACE MASKS

Preformed with an aluminum nose guard. Blue.

- Fiberglass-free
- Not made with natural rubber latex

Molded Masks

(101-6829)..... 50 per Box
5 Boxes @



ADAPTED FACE MASKS

Unique contoured design prevents fogging of preventive eyewear.

- Fiberglass-free
- Not made with natural rubber latex

Molded Masks..... 50 per Box
20 Boxes @

Specify:

Blue (112-5575) Pink (112-5576)

MASK/FACE SHIELD COMBINATIONS



EARLOOP/EYE SHIELD PROCEDURE MASKS LEVEL 3

Soft, comfortable loops and absorbable inner facing make this mask very comfortable for extended wear. Teal.

- Fluid resistant
- Optically clear shield is treated on both sides with antifog material
- Not made with natural rubber latex

(104-7321).....25 per Box
4 Boxes (Case) @

FACE SHIELDS



DISPOSABLE FACE SHIELDS

Optically clear, see-through, full-length visor provides increased protection against splash, spatter, and aerosols. Features a 1½" thick, 9" long foam headband that keeps the face shield off the face to provide maximum comfort and protection. Ultrasonically welded, 13" elastic headband for added strength and increased wearer comfort. Special fog-free coating on inside and outside of visor. Not made with natural rubber latex. Made in the USA.

13"W x 7½"L

(101-2254)..... 24 per Box
4 Boxes @

PROTECTIVE EYEWEAR



MAXI-GARD® PROTECTIVE EYEWEAR

Clear, shatter-resistant polycarbonate lens with scratch-resistant coating. Meets ANSI Z87+ standards. 10-base curve allows for complete coverage. Sure-grip rubber temple insert provides extra comfort.

806 Series.....Ea.
6 @

Specify:

Black (900-4436) Pink (900-4779)
Blue (900-4488) Red (900-4781)



MAXI-GARD® PROTECTIVE EYEWEAR

Clear, shatter-resistant polycarbonate lens with scratch-resistant coating. Meets ANSI Z87+ standards. 9-base curve allows for great wraparound protection and coverage. Stop hinge prevents temple tip from scratching lens. Protective case included. Metal.

811 Series
(900-4437).....Ea.



MAXI-GARD® ROTECTIVE EYEWEAR

Stylish, wraparound protective eyewear for greater comfort and protection. Clear, shatter-resistant polycarbonate lens. 10-base curve for complete coverage. Includes a scratch-resistant, antifog coating and foam inserts for added comfort. Meets ANSI Z87+ standards.

Eyewear.....Ea.

Specify:

Black Frame (570-0434) Pink Frame (570-0431)
Blue Frame (570-0433) Purple Frame (570-0432)



PROTECTIVE EYEWEAR

Lightweight, wraparound, unisex protective eyewear. 10-base curve for complete coverage. Excellent optical quality. Ideal for patient and staff use. Gray lens and temples.
(570-0435).....Ea.



SAFETY GLASSES

Clear plastic lens with side shields for eye protection in the operator or lab. Unisex.

(100-7382).....Ea.
6 @

INFECTION CONTROL

DISPOSABLE GARMENTS

Compare to ValuMax
Extra-Safe Jackets & Lab Coats



MAXI-GARD® JACKETS

Premium protective jacket available in 9 fashionable, attractive colors. Stylish hip-length jacket is made from high-quality SMS fabric that meets OSHA fluid-resistance requirements for protection. Lightweight, breathable, non-latex material for maximum wearer comfort. Detailed finishing includes a 3-pocket design and snap-front closure.

Sizes S–XL.....10 per Pkg.
5 Pkg. (Case)@

Specify:

Small

Black (570-0706)
Blueberry (112-6250)
Ceil Blue (112-6255)
Light Pink (112-6264)
Purple (112-6267)
Raspberry (112-6272)
Sky Blue (112-6277)
Teal (112-6282)
White (112-6249)

Medium

Black (570-0707)
Blueberry (112-6253)
Ceil Blue (112-6258)
Light Pink (112-6265)
Purple (112-6269)
Raspberry (112-6273)
Sky Blue (112-6279)
Teal (112-6284)
White (112-6288)

Large

Black (570-0708)
Blueberry (112-6251)
Ceil Blue (112-6257)
Light Pink (112-6263)
Purple (112-6268)
Raspberry (112-6271)
Sky Blue (112-6275)
Teal (112-6283)
White (112-6286)

X-Large

Black (570-0709)
Blueberry (112-6252)
Ceil Blue (112-6254)
Light Pink (112-6262)
Purple (112-6266)
Raspberry (112-6270)
Sky Blue (112-6274)
Teal (112-6281)
White (112-6285)

Size 2XL–Ceil Blue
(112-6259).....10 per Pkg.
5 Pkg. @

Size 3XL–Ceil Blue
(112-6261).....10 per Pkg.
5 Pkg. @



MAXI-GARD® LAB COATS

Premium protective lab coat available in 9 fashionable, attractive colors. Stylish knee-length coat is made from high-quality SMS fabric that meets OSHA fluid-resistance requirements for protection. Lightweight, breathable, non-latex material for optimal wearer comfort. Detailed finishing includes a 3-pocket design and snap-front closure.

Sizes S–XL.....10 per Pkg.
5 Pkg. (Case)@

Specify:

Small

Black (570-0710)
Blueberry (112-6289)
Ceil Blue (112-6293)
Light Pink (112-6302)
Purple (112-6303)
Raspberry (112-6308)
Sky Blue (112-6311)
Teal (112-6316)
White (112-6323)

Medium

Black (570-0711)
Blueberry (112-6292)
Ceil Blue (112-6295)
Light Pink (112-6301)
Purple (112-6306)
Raspberry (112-6310)
Sky Blue (112-6313)
Teal (112-6318)
White (112-6324)

Large


Black (570-0712)
Blueberry (112-6291)
Ceil Blue (112-6294)
Light Pink (112-6300)
Purple (112-6305)
Raspberry (112-6309)
Sky Blue (112-6314)
Teal (112-6317)
White (112-6321)

X-Large

Black (570-0713)
Blueberry (112-6290)
Ceil Blue (112-6296)
Light Pink (112-6299)
Purple (112-6304)
Raspberry (112-6307)
Sky Blue (112-6312)
Teal (112-6319)
White (112-6320)

Size 2XL–Ceil Blue
(112-6297).....10 per Pkg.
5 Pkg. @

Size 3XL–Ceil Blue
(112-6298).....10 per Pkg.
5 Pkg. @



Compare to ValuMax Dual-Fabric
Surgical Cover Gowns

MAXI-GARD® COVER GOWNS

Knee-length cover gown with complete front coverage. Dual-fabric design offers protection as it maximizes comfort. Sleeves and front of gown are made from quality SMS fabric that meets OSHA fluid-resistance requirements for protection; back of gown is made of a lightweight, breathable, material for wearer comfort. Knitted cuffs and open-back design includes neck and waist ties for individual adjustments. Not made with natural rubber latex. Blue.

Regular
(112-6245)..... 10 per Pkg.
5 Pkg. (Case) @

X-Large
(112-6246)..... 10 per Pkg.
5 Pkg. (Case) @



MAXI-GARD® DISPOSABLE GOWNS

Unique one-piece-design gown with thumb loop on sleeve makes donning gloves easier. Impervious, 100% polyethylene barrier. Easy on/off. Individually wrapped. Convenient dispenser carton for reduced storage requirements. Unisize: one size fits most. Blue.

(104-5806)..... 15 per Box



Compare to ValuMax
Isolation Gowns

MAXI-GARD® ISOLATION GOWNS

Single-use, full-back disposable gowns designed to protect health-care professionals and patients from cross contamination. Made of strong, lightweight, and breathable polypropylene material. Neck ties and belt allow individual adjustments for the perfect fit. Knit cuffs for increased comfort and protection. Not made with natural rubber latex.

Regular.....50 per Pkg.

Specify:
Blue (570-1103) White (570-1109)
Green (570-1108) Yellow (570-1105)
Pink (570-1107)

X-Large..... 50 per Pkg.

Specify:
Blue (570-1104) Yellow (570-1106)

COVER GOWNS

Full-length gown made of fluid-resistant, nonwoven polypropylene material. Full cut, open neck, open back, and waist ties. Sewn neck binder and knit cuffs.

Gowns..... 10 per Pkg.
5 Pkg. (Case) @

Specify:
Blue, Medium/Large (100-5870)
Pink, Small/Medium (100-6784)
Pink, Large (100-9392)

Blue, X-Large
(100-5806)..... 10 per Pkg.
5 Pkg. (Case) @



ISOLATION GOWNS

Durable, one-ply, polypropylene material. Glued construction, no stitch holes. Elastic cuff. Excellent closure.

Regular..... 10 per Pkg.
10 Pkg. (Case) @

Specify:
Blue (101-2585) Yellow (101-1789)

X-Large..... 10 per Pkg.
10 Pkg. (Case) @

Specify:
Blue (101-7475) Yellow (101-4096)



LAB GOWNS

Comfortable, lightweight, snap-front gowns made from fluid-resistant, nonwoven polypropylene. Long sleeves with elastic cuffs. White.

Gowns..... 10 per Pkg.
5 Pkg. (Case) @

Specify:
Small (900-7594) Large (900-7596)
Medium (900-7595) X-Large (900-7597)



POLYCOATED APRONS

Worn over ordinary or other protective apparel, these 26" x 42" full-length aprons are excellent for hygiene as well as added protection. Made of polycoated polypropylene.

(100-7483)..... 25 per Pkg.
4 Pkg. (Case) @

PROTECTIVE GOWNS

Polycoated material protects the wearer from fluid contamination. Long sleeves with elastic cuffs provide a secure closure at wrists. Ties at neck and waist allow total closure. Full cut for comfort. White.

Medium/Large
(100-5849)..... 10 per Pkg.
5 Pkg. (Case) @

X-Large
(100-6865)..... 10 per Pkg.
5 Pkg. (Case) @



PROTECTIVE LAB COATS

Safeguard your clothing from splashes or spills. Made of high-quality polypropylene, each lab coat is lightweight and fluid-resistant, offering the perfect combination of comfort and protection. All coats are full length to provide maximum coverage. Blue, knit cuffs and collar with a snap-front closure.

Blue..... 10 per Pkg.
5 Pkg. (Case) @

Specify:
Small (900-4673) Large (900-4675)
Medium (900-4674) X-Large (900-4676)



CAPS



BOUFFANT CAPS

Full-size cut is sized properly to provide total hair coverage and allows for elastic to be comfortable. Prevents hair loss from contaminating sensitive area or site. Elastic not made with natural rubber latex to minimize allergic reaction.

White
21"
(102-4151)..... 100 per Box
24"
(102-6614)..... 100 per Box
Blue
21"
(104-6032)..... 100 per Box
24"
(104-6058)..... 100 per Box

SHOE COVERS

SHOE COVERS

Non-Sterile
Nonconductive with a nonskid design. One size fits all. Not made with natural rubber latex. Blue
(104-9452)..... 50 Pairs per Box
(104-6350)..... 150 Pairs per Case



HAND SOAPS



ANTIBACTERIAL LIQUID HAND SOAP

Provides a luxurious feel, leaving hands smooth and soft. Mild, gentle and moisturizing. Pleasant cucumber and aloe fragrance. Active ingredient is chloroxylenol.

Pump Bottles
(900-4439)..... 8 oz.
12 x 8 oz. @
(900-4438)..... 16 oz.
12 x 16 oz. (Case)@
Refill
(900-4440)..... Gal.
4 Gal. (Case)@



FOAMING ANTIBACTERIAL HAND SOAP

Antibacterial, foaming hand soap contains a natural-based surfactant that provides conditioning to the skin. Rich lather is naturally moisturizing with pleasant spring water. Clean fragrance. Active ingredient is chloroxylenol.

Pump Bottle
(900-4441)..... 18.6 oz.
12 x 18.6 oz. (Case)@



MAXICLENS® ANTIMICROBIAL HAND WASH

Contains 0.6% PCMX for antimicrobial action and aloe vera to add moisture with each use. The product aids in reducing cross contamination while helping you meet APIC and OSHA recommendations for hand washing.

800-mL Bag for Manual Dispenser
(112-5987)..... Ea.
12 (Case)@
800-mL Bag for Automatic Dispenser
(112-5990)..... Ea.
12 (Case)@
Manual Dispenser
(112-5983)..... Ea.
Automatic Dispenser
(112-5985)..... Ea.
Stand for Dispenser
(112-5989)..... Ea.



MAXICLENS® ANTIMICROBIAL/ANTISEPTIC SKIN CLEANSER

Formulated with detergents, moisturizers, skin protectors, and humectants that cleanse without drying. Fast-acting formula that is proven effective at 15 seconds against a broad spectrum of bacteria. Effective for up to 6 hours after use. Active ingredient: 4% CHG.

(431-0032)..... Qt.
4 Qt. (Case)@
(431-0099)..... Gal.
4 Gal. (Case)@
Hand Pump for Quart Bottle
(431-1478)..... Ea.

WATERLESS HAND GELS & SANITIZERS



MAXICLENS® INSTANT HAND SANITIZER AND DISPENSER

Kills 99.99% of germs in 15 seconds without the use of soap and water. Contains aloe vera to add moisture with every use. Henry Schein's hand sanitizer helps you meet OSHA, APIC, and CDC hand-washing recommendations. Citrus scent. 66.5% ethyl alcohol.

800-mL Bag for Manual Dispenser
(112-5988)..... Ea.
12 (Case)@
800 mL Bag for Automatic Dispenser
(112-5991)..... Ea.
12 (Case)@
Manual Dispenser
(112-5984)..... Ea.
Automatic Dispenser
(112-5986)..... Ea.
Stand for Dispenser
(112-5989)..... Ea.



HAND SANITIZER

Waterless, alcohol-based gel hand sanitizer. Ideal for use where soap and water are not readily available. Rinse-free: no water required. Kills 99.9% of most common disease-causing germs spread by hands. Soothing emollients and a pleasant, fresh, clean scent. Leaves hands soft and clean. 66.5% ethyl alcohol.

Pump Bottles
(900-4443)..... 8 oz.
12 x 8 oz. @
(900-4442)..... 16 oz.
12 x 16 oz. (Case)@

TOWELETTES & WIPES



QUICK WIPES

9" x 5½" premoistened, ready-to-use hand towelettes. Wipes are durable and easy to dispense. Cleans and refreshes without soap and water. Hands will self-dry in seconds, leaving skin soft and smooth.

Pop-up Container
(104-4176).....150 per Cont.
12 Cont. @

HAND CREAMS & LOTIONS



HAND CREAM

Eliminates rough skin and dryness that comes from overexposure to soaps and sanitizers. Contains cucumber and aloe emollients that help to leave hands feeling luxuriously soft.

Pump Bottle
(900-4446)..... 8 oz.
12 x 8 oz. @

STERILANTS

MAXICIDE® PLUS

3.4% glutaraldehyde solution that can be used as a sterilant or high-level disinfectant for up to 28 days. 90 minutes at 77°F (25°C) for high-level disinfection, 10 hours at 77°F (25°C) for sterilization. Contains a surfactant to enhance the disinfection process and a rust inhibitor to reduce corrosion. No dilution required.

(102-5796)..... Qt.
4 Qt. (Case)@

(102-2865)..... Gal.
4 Gal. (Case)@



MAXICIDE®

2.5% buffered glutaraldehyde solution that can be used as a sterilant and high-level disinfectant for up to 28 days. 90 minutes at 77°F (25°C) for high-level disinfection, 10 hours at 77°F (25°C) for sterilization. Contains a surfactant to enhance the disinfection process and a rust inhibitor to reduce corrosion. No dilution required.

(102-6003)..... Qt.
4 Qt. (Case)@

(102-7394)..... Gal.
4 Gal. (Case)@



MAXICIDE® ACTIVATED

Preactivated 2.65% glutaraldehyde solution. Can be used as a sterilant and high-level disinfectant for up to 30 days. Tuberculocidal in 45 minutes. Ready to use with a fresh, pleasant scent.

Bottle
(104-4681)..... Gal.
4 Gal. (Case)@



MAXICIDE® OPA 28

Non-glutaraldehyde 0.575% ortho-phthalaldehyde (OPA) solution. Ready-to-use high-level disinfectant that achieves disinfection in just 10 minutes at room temperature (20°C) and offers a reuse period of up to 28 days. Significantly reduces the risk of occupational hazards including eye, skin, and respiratory tract irritations and sensitization. Test strips should be used to verify the minimum recommended concentration (MRC) prior to each use.

(112-6976)..... Gal.
4 (Case)@

Test Strips
(112-6977)..... 50 per Bottle



INSTRUMENT TRAYS



GERMICIDE TRAYS

For cold sterilization and disinfection of instruments. Tray automatically lifts and drains when cover is opened. Removable folding tray handle. High impact plastic, available with stainless steel or transparent cover. Resistant to all phenolic solutions. Includes bur tray. Approximately a 1-Qt. capacity.

Standard size tray: 3"H x 10 $\frac{3}{4}$ "W x 5"D

Stainless Cover..... Ea.

Specify:

Beige (100-8193) White (100-8485)
Blue (100-7778) Mauve (101-8195)

Opaque Smoke

Plastic Cover..... Ea.

Specify:

Beige (100-6126) Mauve (101-8196)
Blue (100-7564)

AUTOClave BAGS



AUTOClave BAGS

Built-in Indicator/Paper, 2 $\frac{1}{2}$ " Wide, 1 $\frac{1}{2}$ " Opened Sides, Class 1 Indicator For use in steam sterilizers.

6"
(100-0166)..... 1000 per Box
5 Boxes@

8"
(100-0529)..... 1000 per Box
5 Boxes@

10"
(100-1686)..... 1000 per Box
5 Boxes@



SELF-SEAL AUTOCLAVE BAGS

Heavy-duty, surgical-grade kraft paper. Dual indicators for steam or chemical vapor. Its self-sealing feature ensures a high-integrity seal. Class 1 indicator.

Dimensions: 2 $\frac{1}{2}$ " x 1 $\frac{1}{2}$ " x 10 $\frac{1}{2}$ "
(101-5915)..... 1000 per Box

STERILIZATION POUCHES



SELF-SEAL STERILIZATION POUCHES WITH DUAL INDICATORS

Superior-quality pouch with dual indicators. The internal and external indicators are preprinted on same side and change color when processed with steam and EtO. Additional process indicator strips are not required with every use. Meets CDC infection control guidelines for effective sterilization process. Made with medical-grade paper and an easy-to-open thumb notch. Class 1 indicator. Pouches and packaging are lead-free and not made with natural rubber latex.*

2 $\frac{1}{4}$ " x 4"
(112-4859)..... 200 per Box
12 Boxes (Case)@

2 $\frac{3}{4}$ " x 9"
(112-4856)..... 200 per Box
6 Boxes (Case)@

3 $\frac{1}{2}$ " x 5 $\frac{1}{4}$ "
(112-4857)..... 200 per Box
6 Boxes (Case)@

3 $\frac{1}{2}$ " x 9"
(112-4853)..... 200 per Box
6 Boxes (Case)@

3 $\frac{1}{2}$ " x 9"
(112-4855)..... 500 per Box
8 Boxes (Case)@

3 $\frac{3}{4}$ " x 12"
(112-4866)..... 200 per Box
6 Boxes (Case)@

4 $\frac{1}{4}$ " x 11"
(112-4863)..... 200 per Box
6 Boxes (Case)@

5 $\frac{1}{4}$ " x 10"
(112-4854)..... 200 per Box
6 Boxes (Case)@

5 $\frac{1}{4}$ " x 12"
(112-4861)..... 200 per Box
6 Boxes (Case)@

7 $\frac{1}{2}$ " x 13"
(112-4858)..... 200 per Box
6 Boxes (Case)@

8" x 13"
(112-4864)..... 100 per Box
6 Boxes (Case)@

9" x 15"
(112-4860)..... 100 per Box
4 Boxes (Case)@

10 $\frac{1}{2}$ " x 17"
(112-4862)..... 100 per Box
4 Boxes (Case)@

13" x 20"
(112-4865)..... 100 per Box
4 Boxes (Case)@

*The use of weekly biological indicators is recommended.



MAXIMA® MULTIPARAMETER SELF-SEAL STERILIZATION POUCHES WITH SURE-CHECK® TECHNOLOGY

Confirms that all 3 parameters for sterilization have been met: time, temperature, and presence of steam. Internal and external multiparameter indicators are printed directly on the pouch. Eliminates the extra step and cost of adding separate multiparameter strips to each pouch. Medical-grade paper with triple seals; chevron safety seal at base for added protection; thumb notch for easy opening. Not made with natural rubber latex. Made in the USA. Class 1 indicator.

2 1/4" x 4" (900-7838)..... 200 per Box
6 Boxes @

2 3/4" x 9" (900-4602)..... 200 per Box
6 Boxes @

3 1/2" x 5 1/4" (900-4604)..... 200 per Box
6 Boxes @

3 1/2" x 9" (900-4601)..... 200 per Box
6 Boxes @

5 1/4" x 10" (900-4603)..... 200 per Box
6 Boxes @

7 1/2" x 13" (900-4782)..... 200 per Box
5 Boxes (Case) @

12" x 18" (900-4783)..... 100 per Box
5 Boxes (Case) @



DUAL PURPOSE NYLON FILM POUCHES

All nylon film with indicator. Self-sealing. Can be used in steam or dry heat. Class 1 indicator.

2" x 9 1/2" (101-3363)..... 100 per Box
12 Boxes (Case) @

3" x 10" (101-3373)..... 100 per Box
12 Boxes (Case) @

4" x 10" (101-3383)..... 100 per Box
12 Boxes (Case) @

6" x 11" (101-3393)..... 100 per Box
12 Boxes (Case) @

HEAT SEALERS



IMPULSE HEAT SEALER

Seals nylon tubing, then automatically turns off. Seals perfectly within 1-2 seconds. UL certified and durable.

110-Volt Unit (100-3048)..... Ea.

Replacement Kit (101-5802)..... Ea.

Contains: 2 heating elements, 1 upper cloth & 1 lower cloth.

TUBING



DUAL-PURPOSE NYLON STERILIZATION TUBING

Transparent: contents are easy to identify. Sturdy: resists tearing and puncturing. Dual purpose: can be used in steam or dry heat. 100 ft. per roll. Class 1 indicator.

Continuous Clear 2" (101-2224)..... Ea.
12 (Case) @

3" (101-2225)..... Ea.
12 (Case) @

4" (101-2226)..... Ea.
12 (Case) @

6" (101-2227)..... Ea.
12 (Case) @

Continuous with Indicator 2" (101-1111)..... Ea.

3" (101-1110)..... Ea.

4" (101-2222)..... Ea.

6" (101-2223)..... Ea.

Perforated with Indicator 2" x 9 1/2" (101-3323)..... Ea.

3" x 10" (101-3333)..... Ea.

4" x 10" (101-3343)..... Ea.

6" x 11" (101-3353)..... Ea.

INDICATOR TAPE



LEAD-FREE STEAM INDICATOR TAPE

Resilient, pressure-sensitive, and lead-free adhesive tape designed for use on disposable wraps. The backing is a natural- or blue-colored, saturated crepe paper that provides secure wraps throughout sterilization expansion. Indicator lines will change from beige to dark brown/black when proper levels of moisture and temperature have been achieved. ISO 11140-1:2014 Type 1 compliant.

1/2" Roll (570-1511)..... Ea.

3/4" Roll (570-1510)..... Ea.

1" Roll (570-1509)..... Ea.

INDICATOR STRIPS



MAXITEST® CLASS 5 INTEGRATOR

Class 5 sterilization indicator designed to integrate the 3 variables of sterilization: time, temperature, and saturated steam. Immediate readout results. When proper sterilization parameters have been met, the integrator will enter the blue (safe) area. For use with all steam processes: gravity, prevacuum, and flash.

(570-0310)..... 100 per Bag
(570-0311)..... 1000 per Bag



INDICATOR STRIPS

For steam sterilizers. Indicator ink changes from white to black, indicating sterilization process has been completed. Class 1 indicator.

(102-6332)..... 250 per Box



DRY HEAT STERILIZER INDICATOR LABELS

For dry heat sterilizers. Indicator ink changes from green to black, visually indicating that the sterilization process has been completed. Class 1 indicator.

(102-7982)..... 1000 per Box

STERILIZATION CASSETTES

INSTRUMENT CASSETTES

Stainless steel locking tray complete with silicone rails for storage of dental instruments. Fully autoclavable, features a unique design that optimizes the sterilization process. Upper and lower sections of the tray are fully interchangeable. Colored rails for instrument identification and organization.



5-Instrument Cassette
7.1" x 2.8" x 1.3"..... Ea.

Specify:
Blue (112-4963) Red (112-5297)
Green (112-5301) Yellow (112-5442)
Purple (112-5336)

5-Instrument
Replacement Rails.....3 per Bag

Specify:
Blue (112-5455) Red (112-5456)
Green (112-5457) Yellow (112-5459)
Purple (112-5458)



10-Instrument Cassette
7.1" x 5.6" x 1.3"..... Ea.

Specify:
Blue (112-4964) Red (112-5299)
Green (112-5302) Yellow (112-5443)
Purple (112-5439)

10-Instrument
Replacement Rails.....3 per Bag

Specify:
Blue (112-5460) Red (112-5461)
Green (112-5462) Yellow (112-5463)
Purple (112-5464)



16-Instrument Cassette
11.3" x 7.1" x 1.3"..... Ea.

Specify:
Blue (112-4965) Red (112-5300)
Green (112-5303) Yellow (112-5444)
Purple (112-5441)

16-Instrument
Replacement Rails.....6 per Bag

Specify:
Blue (112-5465) Red (112-5466)
Green (112-5468) Yellow (112-5470)
Purple (112-5469)

CSR WRAPS



CSR WRAP

For general-purpose sterilization, providing a barrier against airborne and waterborne bacteria. Conforms easily around various sizes and shapes. Provides a barrier yet allows excellent penetration for various methods of sterilization including steam, ethylene oxide, radiation, and E-beam. Autoclavable 1-ply sheets.

12" x 12"..... 100 per Case
(100-3347)..... 10 Cases@

15" x 15"..... 100 per Case
(102-5193).....

18" x 18"..... 100 per Case
(100-4683).....

24" x 24"..... 100 per Case
(100-6950)..... 6 Cases@

30" x 30"..... 100 per Case
(100-8187).....

36" x 36"..... 100 per Case
(100-2016).....



DRI-GARD® CSR WRAP

All-purpose sterilization overwraps. Provide high protection against air- and fluid-borne particles, yet allow effective penetration for various methods of sterilization, including steam, ethylene oxide, radiation, and E-beam. Resistant to tears and punctures. Conform to various pack and tray sizes and are easily secured with standard adhesives or tapes for full enclosure and assured sterility. Autoclavable, 1-ply sheets.

12" x 12"..... 500 per Case
(900-7447).....

15" x 15"..... 500 per Case
(900-7448).....

20" x 20"..... 500 per Case
(900-7449).....

24" x 24"..... 500 per Case
(900-7450).....

BIOLOGICAL MONITORS



MAXITEST®

Biological Monitoring System
Mail-in service includes 3 bacterial spore test strips (2 test strips, 1 control strip) per test package. Postage prepaid; culturing service included. Monitors all sterilizers including steam, chemical vapor, dry heat, and EtO. Third-party documentation serves as tangible legal document where proof of sterility testing must be provided. Immediate and confidential notification of test failures. Test results are available electronically.

12-Test Package
(101-5830)..... Ea.
Contains: 12 packets, each containing 2 test strips & 1 control strip.

48-Test Package
(101-1253)..... Ea.
Contains: 48 packets, each containing 2 test strips & 1 control strip.

Not for export.



MAXITEST® VALUE TEST PACKS

Biological monitoring system used to test steam, dry heat, chemical vapor, and EtO processes. Mail-in service with results reported quarterly. Immediate and confidential notification of test failures. Test results available electronically. Postage not included.

(101-4715)..... 52 per Box
Contains: 52 test packets, each containing 1 test strip & 1 control strip.

Not for export.



MAXITEST® IN-OFFICE BIOLOGICAL MONITORING SYSTEM

In-office system for monitoring steam sterilizers. Easy to use and interpret. Results in 24-48 hours. Simply process vial in a normal sterilizer cycle. Then activate vial, incubate, and record results. If any spores survive the process, the vial medium will turn yellow. MaxiTest Incubator is 24 volts and holds up to 14 indicators. Indicators are also compatible with the 3M™ Attest™ incubator.

Starter Kit
(102-6076)..... Ea.
Contains: 1 of each: box of 25 indicators, incubator, crusher & record keeper.

Indicator Refills
(102-5804).....25 per Box
(112-6848)..... 100 per Box

Biological Monitors continued
on next page.

MAXITEST® BOWIE-DICK TEST PACK

Designed to evaluate the efficacy of air removal and steam penetration in prevacuum steam sterilizers. Immediate readout results. When proper sterilization conditions are met, test pack changes color. Each test pack contains frequency of use, instructions, and storage/shelf life information preprinted on the box for ease of use. Use daily.

30-Test Package
(570-0309).....Ea.



MAXITEST® STEAM BI TEST PACK

In-office biological indicator (BI) system for gravity and prevacuum steam sterilizers. Allows for immediate release of processed loads through the use of MaxiTest® Class 5 Steam Integrator, which is permanently attached to a load record card. Test pack complies with standards for routine sterilizer efficacy monitoring and biologically testing implant loads. Meets AAMI recommendations for implant loads. Media growth promotion is validated per USP.

25-Test Packet
(570-0312).....Ea.



MAXITEST® STERILIZER TEST PACK

Designed to monitor sterilization conditions in both gravity displacement and prevacuum sterilizers. Immediate readout results require no incubation. Load record card provides permanent, physical proof of sterilization conditions.

25-Test Packet
(570-0313).....Ea.

ULTRASONIC SOLUTIONS



SUPER CONCENTRATED ULTRASONIC CLEANING SOLUTIONS

Fast-acting solutions in small, easy-to-store, 16 oz. bottles. Metered-dose dispenser top ensures exact dose. Available in 3 formulas: General Purpose, for very fast cleaning of blood, mucus, and tissue; Tartar & Stain Remover's citric acid formula is noncorrosive and cleans quickly; and Enzymatic Solution is effective against blood, tissue, and mucus found on instruments.

Metered Dose/Bottles..... 16 oz.
6 x 16 oz. (Case)@

Specify:

General Purpose Cleaner
(Makes 16 Gal.) (101-4412)

Tartar and Stain Remover
(Makes 96 oz.) (101-0827)

Enzymatic Cleaning Solution
(Makes 16 Gal.) (102-6472)

ULTRASONIC CLEANING SOLUTIONS

General Purpose

A nonammoniated concentrate. Dilute 1:10 for use. Effectively cleans such items as dried blood, tissue, oxide compounds, and more. Removes ions from instrument's surface to enhance the effectiveness of the ultrasonic cleaning unit.

(100-3836)..... Gal.
4 Gal. (Case)@



Tartar and Stain Remover

Used in a beaker in the ultrasonic unit to remove tartar and stains from dentures. Ready to use, so it is easy and convenient. Potent stain remover to eliminate stains on dentures, crowns, and bridges.

(100-4263)..... Gal.
4 Gal. (Case)@



Temporary Cement Remover ♦

Ready-to-use formula, so no mixing is required. Strong, fast-acting formula removes cement quickly and completely.

(100-5168).....Gal.
4 Gal. (Case)@



Plaster and Stone Remover

(101-3277)..... Gal.
4 Gal. (Case)@



♦ Hazardous material



MAXITABS™ ULTRASONIC CLEANING TABLETS

General Purpose Cleaner

Provide safe and effective cleaning properties in multiple cleaning cycles throughout the day. Does not leave a residue on instruments or inside tank. One tablet makes 1 gallon of solution; 1 box makes 32 gallons.

(900-4088)..... 32 per Box

Tartar & Stain Remover

Offers safe and effective cleaning efficacy for removal of calculus, tartar, tobacco, and food stains. Specifically designed for dentures, bridges, and dental appliances. 1 tablet makes 5 oz. of solution.

(900-4089)..... 32 per Box



MAXIZYME®

Low-foaming enzymatic detergent for use as a presoak manual detergent, ultrasonic detergent, and evacuator cleaner. Dissolves blood, protein, tissue and soil. Noncorrosive, nonflammable, and rust-inhibiting. Concentrated and economical.

(101-9031)..... Gal.
4 Gal. (Case)@

MAXIZYME® TABS

Enzyme cleaning tablets for ultrasonic cleaning, evacuation-system cleaning, and instrument presoak. 2 tablets dissolve easily in 1 gallon of warm tap water. Easy to use and requires less storage space.

(900-0475)..... 64 per Box



ULTRA CLEAN

General-purpose, nonammoniated ultrasonic cleaner. Removes dried blood, tissue, soil, plaster, stone, pumice, and other light contaminants. Nonflammable and nonfoaming. Supplied as powder in container. Concentrated solution: dilute 1 part to 9 parts of water.

1-gal Bottle
(102-2255).....Gal.
4 Gal. (Case)@





ULTRA CLEAN-L

General-purpose, nonammoniated, and super concentrated. Easier to dilute than powders. Extra cleaning power and pleasant scent. Add water to fill container; dilute 1:10 for working-strength solution.

(100-3185).....Gal.
4 Gal. (Case)@



UNIT DOSE ULTRASONIC CLEANING SOLUTIONS

Single-dose, super-concentrated powders for use in ultrasonic cleaning units, auxiliary pans, and beakers. Eliminates lifting and storing heavy gallon-size containers. Cleans instruments directly in ultrasonic cleaning tank.

General Purpose Cleaner

For cleaning instruments directly in the ultrasonic cleaning tank. One packet of powder makes 1 gal. of solution.

(101-0377).....24 per Box

Tartar & Light Stain Remover

For removing tartar and light stains from dentures, partials, and bridges. Dissolve 1 packet of powder into a beaker of 400 mL warm tap water. Place beaker in ultrasonic cleaner.

(101-6767).....24 per Box



ORTHO ULTRASONIC SOLUTION

General-purpose cleaner for the removal of blood, serum, and debris from instruments. Also cleans buffing compounds, rouge, tripoli, oxidizers, burs, glass slabs, metal parts, gold, and porcelain. Special aluminum-safe, no-rinse formula with corrosion inhibitors. Concentrated, 40:1 dilution.

(104-9207).....Gal.
4 Gal. (Case)@

INSTRA-CLEAN

Formulated especially for cleaning all dental instruments. Neutral pH: very mild and kind to your hands. Built-in antirust compound. Concentrated lemon scent and convenient dispenser container.

8-oz Bottle
(100-4250).....8 oz.
4 x 8 oz. (Case)@

STAIN KLEEN

For tartar, light stains, and permanent cement removal. Also used for solder, flux, rust, oxides, and investments. Supplied as a powder in container. Add water for ready-to-use solution.

(100-4179).....Gal.
4 Gal. (Case)@

SURFACE BARRIERS



SURFACE BARRIERS

Disposable barrier products to aid in the prevention of cross contamination.

Light Handle Covers—4"W x 5¾"L
For "T" bar handles.

(570-1316).....500 per Box

Syringe Covers—2½"W x 10"L

For 3-way air/water syringes.

(570-1320).....500 per Box

Syringe Covers with Opening

2½"W x 10"L

For 3-way air/water syringes.

(104-7423).....500 per Box
4 Boxes@

Chair Covers—27½"W x 24"L

(570-1933).....225 per Box
4 Boxes@

Cover Film

Disposable sheets with an adhesive back.

4"W x 6"L.....1200 per Roll

Specify:

Clear (570-1315) Blue (570-1318)

Curing Light Sleeves

(112-7093).....100 per Box
5 Boxes@

Tray Barriers—10½"W x 14"L

(570-1321).....500 per Box

X-ray Covers—15"W x 26"L

(570-1319).....250 per Roll

DISPOSABLE SYRINGE TIPS



SURE TIP®

Disposable air/water syringe tip. Durable, all-plastic construction. Dry air every time, no water leaking. Clear plastic allows visibility of water flow. No adapter or converter kit needed with original equipment. Will fit most models of air/water syringes if not previously converted to disposable tips.

Syringe Tips

(102-8981).....250 per Pkg.
4 Pkg. @

Bulk Package—Clear

(570-0104).....1000 per Pkg.

Conversion Kit for ADEC and Screw-In Types

(642-1941).....Ea.



Compare to Crystal Tips

SURE TIP® COLORS

Disposable air/water syringe tips in attractive colors. Durable, all-plastic construction prevents cross contamination. Dry air every time: no leaking water. Four air channels enhance spray and effectively separate water and air flow. Interchangeable with most metal tips: no conversion is necessary.

Syringe Tips.....250 per Bag
4 Bags @

Specify:

Assorted (900-7022) Blue (900-7021)

Bulk Package.....1000 per Bag

Specify:

Assorted (900-7972) Blue (900-7973)



SURE TIP® PLUS

Disposable air/water syringe tip. Five-air-channel design guarantees a perfect spray and instant air/water separation. Plastic exterior with a metal core for strength. Interchangeable with most metal tips without conversion. Locking indicator and locking groove ensure tip is seated securely.

White

(900-7030).....150 per Pkg.
6 Pkg. @

(900-7031).....1600 per Pkg.



Compare to Starz Tips

SURE TIP® FOCUS

Disposable, durable, 6-air-channel design enhances air/water spray pattern, as compared to metal tips, while preventing cross contamination. Easily bendable to create custom angles and gentle on soft tissue to provide a more comfortable patient experience. No adapter or converter kit is required with original equipment. Includes a nonscratching end and locking indicator that is gentle on the O-ring in the syringe head, preventing water crossover. Available in an attractive assortment of colors.

Syringe Tips

(112-7182).....150 per Bag
4 Bags @

Bulk Package

(112-7274).....1000 per Pkg.

SURFACE CLEANERS & DISINFECTANTS



MAXISPRAY® PLUS

Broad-spectrum, tuberculocidal, disinfectant/decontaminant cleaner. Effective against viruses (hydrophilic and lipophilic), bacteria (including VRE), and fungi. Kills MRSA and hepatitis B and C in 2 minutes and TB in 3 minutes. May be used as an immersion solution for noncritical instruments and as a surface disinfectant. Can be used on plastics, stainless steel, chrome, glass, Plexiglas®, and painted surfaces. Contains no phenols. Nontoxic and nonirritating. Easy, ready-to-use formula. No messy mixing, measuring, or rinsing required. Intermediate level.

Pump Spray
(101-4114)..... 24 oz.
12 x 24 oz. (Case)@

Refill
(101-0285)..... Gal.
4 Gal. (Case)@

TOWELETTES



MAXIWIPE® LA DISINFECTING TOWELETTES

Low-alcohol disinfectant towelettes that clean and disinfect hard, nonporous surfaces and patient-care equipment. These fast-action wipes kill TB in 3 minutes and MRSA, HIV-1, HBV, and HCV in 2 minutes. The durable, non-woven towelettes will not bunch up during use. The deep-well lid and easy-access tab provide quick opening and tight closures between wipes. Towelette is an all-in-one cleaner and disinfectant. The XL wipe is 3X the size of the regular wipe. Wipes contain 17.2% alcohol. Intermediate level.

Large—6" x 6¾"
(112-6843)..... 160 per Canister
12 Canisters (Case)@

X-Large—9" x 12"
(112-6844)..... 65 per Canister
6 Canisters@

MAXIWIPE® GERMICIDAL CLOTHS

Disinfectant wipe that kills TB and RSV in 1 minute; HBV, HCV, HIV-1, rhinovirus, MRSA, VRE, influenza A (H1N1) virus, *Acinetobacter baumannii*, and *Campylobacter jejuni* in 2 minutes. Thick, absorbent, hexagonal-embossed material helps remove soil from surfaces without tearing. Pliable to get into hard corners of equipment. Unique deep-well lid on canister seals securely to prevent moisture loss. Cleans and disinfects hard, nonporous surfaces and patient-care equipment with a fast kill time of 2 minutes. Bactericidal, tuberculocidal, and virucidal. Contains 55% alcohol. Intermediate level.

Large—6" x 6¾"
(106-6794)..... 160 per Canister
12 Canisters (Case)@

X-Large—7½" x 15"
(112-5305)..... 65 per Canister
6 Canisters@



WASTE BAGS



INFECTIOUS WASTE COLLECTION BAGS

Constructed of high-grade polyethylene, for safe collection, handling and disposal of infectious wastes destined for incineration or a landfill. Red, labeled with biohazard symbol.

Small—10 Gallons, 15"W x 9"H x 24"D
(100-9305)..... 100 per Box
Medium—18 Gallons, 15"W x 9"H x 31"D
(100-8878)..... 100 per Box
Large—33 Gallons, 23"W x 8"H x 41"D
(100-3416)..... 100 per Box

SHARPS DISPOSAL



SHARPS CONTAINER

Shatter-resistant, puncture-resistant, and transparent lid allows for easy viewing. Ideal for safe disposal of needles, needles with syringes, vacutainers, lancets, and blades (e.g., razors and scalpels). Counterbalanced lid is designed to prevent infection/cross contamination by limiting access to container contents and facilitating one-handed operation. Constructed of puncture-resistant polypropylene. Safe to autoclave or incinerate and meets or exceeds FDA standards, including OSHA and EPA requirements, for local biohazard waste-collection services.

1 Qt.
(112-7149)..... Ea.
5 Qt.
(112-7152)..... Ea.
5.4 Qt.
(112-7151)..... Ea.
2 Gal.
(112-7150)..... Ea.

SHARPS CONTAINER WITH ROTARY LID

Rotating lid allows for maximum storage and safe disposal of larger sharps. High-impact polypropylene material. Fracture- and puncture-resistant. Lid locks for final disposal.

2.2 Qt.
(570-0213)..... Ea.
3 Gal.
(570-0212)..... Ea.



PHARMACEUTICAL DISPOSAL



MEDSAFE EASY-PAK

Mail-Back Pharmaceutical Returns Service. Unique reverse pharmaceutical distribution service for expired pharmaceuticals that ensures compliant destruction. "Pay as you go" system is FDA, DEA, and EPA compliant nationwide. Place all of the following pharmaceuticals in one box: Schedule II (CII), Schedule III–V (CVs), and all noncontrolled pharmaceuticals (Rx, OTC, hazardous, and compound).^{*} Includes prepaid return shipping label.

Envelope
(570-0387)..... Ea.
2 Gal.
(570-0390)..... Ea.
5 Gal.
(570-0391)..... Ea.
10 Gal.
(570-0389)..... Ea.
15 Gal.
(570-0388)..... Ea.

^{*}Vaccines and in-date products are not accepted and should go through the general return process.

AUTOCLAVE CLEANERS



AUTOCLAVE CLEANER

A mild, slightly alkaline liquid that acts as a powerful, concentrated cleaner to remove grease, grime, and accumulation of scale deposits. Prolongs the life of the autoclave with routine maintenance. Mix 4 oz. of cleaner into 64 oz. of distilled water. Excellent for use with steam. Cost-effective 32 oz. bottle makes up to 544 oz. (4.2 gal.) of cleaning solution.

Bottle
(112-5839)..... 32 oz.
4 x 32 oz. (Case)@