


SURGICRYL[®]

910

POLYGLACTINE



Type	Braided and coated multifilament
Composition	Polyglactine 910, a copolymer made of 90% glycolide and 10% L-lactide
Coating	Poly-glycolide-co-L-lactide and calcium stearate (<1%)
Colour	Violet
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption between 56 and 70 days. Approximately 50% of tensile strength remain after 21 days.
Presentation	Box with 12 sutures; Cassette



*For quality
&
safety*

smi

www.sutures.be



CHARACTERISTICS

- ☒ High tensile strength
- ☒ Synthetic
- ☒ Breaks down by hydrolysis
- ☒ Braided structure

- ☒ Coated multifilament

BENEFITS

Excellent tissue support

Minimal tissue reaction

Predictable absorption rate

Excellent handling properties
High flexibility
High tensile strength
Excellent knot security

Soft passage through the tissues with minimal sawing,
tissue drag and trauma
Good knot tie down / adjustability
Reduced capillary action

INDICATIONS General soft tissue closing
and/or ligation


especially in General surgery
Skin closure
Gastrointestinal surgery
Plastic surgery
Ophthalmic surgery
Gynaecology
Obstetrics
Orthopaedics
Urology

CONTRAINDICATIONS Cardiovascular surgery
Neurosurgery


Being absorbable not to use when an
extended wound support is required.

Caution in patients with delayed
wound healing (e.g. compromised
immune system).

TENSILE STRENGTH



ABSORPTION





Without needle

Metric	USP	Length		Order number
1.5	4/0	150 cm	<i>violet</i>	15150015
2	3/0	150 cm	<i>violet</i>	15200015
3	2/0	150 cm	<i>violet</i>	15300015
3.5	0	150 cm	<i>violet</i>	15350015
4	1	150 cm	<i>violet</i>	15400015
5	2	150 cm	<i>violet</i>	15500015
3	2/0	3 x 45 cm	<i>violet</i>	15300045-3
3.5	0	3 x 45 cm	<i>violet</i>	15350045-3
1.5	4/0	6 x 45 cm	<i>violet</i>	15150045-6
2	3/0	6 x 45 cm	<i>violet</i>	15200045-6
3	2/0	6 x 45 cm	<i>violet</i>	15300045-6
3.5	0	6 x 45 cm	<i>violet</i>	15350045-6
1.5	4/0	12 x 45 cm	<i>violet</i>	15150045
2	3/0	12 x 45 cm	<i>violet</i>	15200045
3	2/0	12 x 45 cm	<i>violet</i>	15300045
3.5	0	12 x 45 cm	<i>violet</i>	15350045
4	1	12 x 45 cm	<i>violet</i>	15400045
5	2	12 x 45 cm	<i>violet</i>	15500045
2	3/0	2 x 70 cm	<i>violet</i>	15200070-2
3	2/0	2 x 70 cm	<i>violet</i>	15300070-2
3.5	0	2 x 70 cm	<i>violet</i>	15350070-2
4	1	2 x 70 cm	<i>violet</i>	15400070-2
5	2	2 x 70 cm	<i>violet</i>	15500070-2
4	1	7 x 75 cm	<i>violet</i>	15400075



Without needle

Ligapack



Metric	USP	Length		Order number
2	3/0	250 cm	<i>violet</i>	15200025
3	2/0	250 cm	<i>violet</i>	15300025
3.5	0	250 cm	<i>violet</i>	15350025
4	1	250 cm	<i>violet</i>	15400025
5	2	250 cm	<i>violet</i>	15500025



Round bodied taper point

1/2 circle



		Metric	USP	Length		Order number
1/2 circle	HR 13 mm		0.7	6/0	75 cm	<i>violet</i> 15070113
			1	5/0	75 cm	<i>violet</i> 15100113
			1.5	4/0	75 cm	<i>violet</i> 15150113
	HR 17 mm		0.7	6/0	75 cm	<i>violet</i> 15070117
			1	5/0	75 cm	<i>violet</i> 15100117
			1.5	4/0	75 cm	<i>violet</i> 15150117
			2	3/0	75 cm	<i>violet</i> 15200117
			3	2/0	75 cm	<i>violet</i> 15300117
	HR 22 mm		1.5	4/0	75 cm	<i>violet</i> 15150122
			2	3/0	75 cm	<i>violet</i> 15200122
			3	2/0	75 cm	<i>violet</i> 15300122
			4	1	75 cm	<i>violet</i> 15400122
	HR 26 mm		1	5/0	75 cm	<i>violet</i> 15100126
			1.5	4/0	75 cm	<i>violet</i> 15150126
			2	3/0	75 cm	<i>violet</i> 15200126
			3	2/0	75 cm	<i>violet</i> 15300126
			3.5	0	75 cm	<i>violet</i> 15350126
			4	1	75 cm	<i>violet</i> 15400126
			5	2	75 cm	<i>violet</i> 15500126
	HR 30 mm		2	3/0	75 cm	<i>violet</i> 15200130
			3	2/0	75 cm	<i>violet</i> 15300130
			3.5	0	75 cm	<i>violet</i> 15350130
			4	1	75 cm	<i>violet</i> 15400130
			5	2	90 cm	<i>violet</i> 15500130
	HR 36 mm		2	3/0	75 cm	<i>violet</i> 15200136
			3	2/0	75 cm	<i>violet</i> 15300136
			3	2/0	75 cm	<i>undyed</i> 15300136U
			3.5	0	90 cm	<i>violet</i> 15350136
			4	1	90 cm	<i>violet</i> 15400136
			5	2	90 cm	<i>violet</i> 15500136



Round bodied taper point

1/2 circle



		Metric	USP	Length		Order number
HR 40 mm		2	3/0	75 cm	<i>violet</i>	15200140
		3	2/0	75 cm	<i>violet</i>	15300140
		3.5	0	75 cm	<i>violet</i>	15350140-75
		3.5	0	90 cm	<i>violet</i>	15350140
		4	1	75 cm	<i>violet</i>	15400140-75
		4	1	90 cm	<i>violet</i>	15400140
		5	2	75 cm	<i>violet</i>	15500140-75
		5	2	90 cm	<i>violet</i>	15500140
HR 45 mm		4	1	90 cm	<i>violet</i>	15400145
		5	2	90 cm	<i>violet</i>	15500145
HR 48 mm		3	2/0	90 cm	<i>violet</i>	15300148
		3.5	0	90 cm	<i>violet</i>	15350148
		4	1	90 cm	<i>violet</i>	15400148
		5	2	90 cm	<i>violet</i>	15500148
HR 60 mm		4	1	75 cm	<i>violet</i>	15400160
		5	2	90 cm	<i>violet</i>	15500160
HR 65 mm		3	2/0	90 cm	<i>violet</i>	15300165
		3.5	0	90 cm	<i>violet</i>	15350165
HR 80 mm		4	1	100 cm	<i>violet</i>	15400180
		5	2	100 cm	<i>violet</i>	15500180




Round bodied taper point double armed

1/2 circle



		Metric	USP	Length		Order number
2 x HR 26 mm		2	3/0	75 cm	<i>violet</i>	15200126-2
		3	2/0	75 cm	<i>violet</i>	15300126-2

 Round bodied taper point loop	Metric	USP	Length	Order number		
1/2 circle 	HR 45 mm 	3.5	0	150 cm	<i>violet</i>	15350945
	HR 65 mm 	5	2	240 cm	<i>violet</i>	15500965

 Round bodied taper point			Metric	USP	Length		Order number
 5/8 circle	FR 26 mm		2	3/0	75 cm	violet	15200326
	FR 27 mm		3	2/0	75 cm	violet	15300327
			3.5	0	75 cm	violet	15350327

 Round bodied taper point	Metric	USP	Length	Order number		
3/8 circle 	DR 17 mm 	1.5	4/0	45 cm	<i>violet</i>	15150217
		2	3/0	75 cm	<i>violet</i>	15200217
	DR 22 mm 	2	3/0	75 cm	<i>violet</i>	15200222
	DR 30 mm 	4	1	75 cm	<i>violet</i>	15400230
	DR 32 mm 	3	2/0	90 cm	<i>violet</i>	15300232
		3.5	0	75 cm	<i>violet</i>	15350232

 Reverse cutting			Metric	USP	Length		Order number
1/2 circle 	HS 12 mm		0.7	6/0	45 cm	<i>violet</i>	15071412
	HS 17 mm		1	5/0	75 cm	<i>violet</i>	15101417
			1.5	4/0	75 cm	<i>violet</i>	15151417
	HS 22 mm		2	3/0	75 cm	<i>violet</i>	15201422
	HS 23 mm		5	2	90 cm	<i>violet</i>	15501423
	HS 26 mm		2	3/0	75 cm	<i>violet</i>	15201426
			3	2/0	75 cm	<i>violet</i>	15301426
			3.5	0	75 cm	<i>violet</i>	15351426
			4	1	75 cm	<i>violet</i>	15401426
	HS 30 mm		2	3/0	75 cm	<i>violet</i>	15201430
	HS 36 mm		3	2/0	75 cm	<i>violet</i>	15301436
			3.5	0	90 cm	<i>violet</i>	15351436
			4	1	90 cm	<i>violet</i>	15401436
			5	2	90 cm	<i>violet</i>	15501436
	HS 40 mm		2	3/0	90 cm	<i>violet</i>	15201440
			3	2/0	90 cm	<i>violet</i>	15301440
			3.5	0	90 cm	<i>violet</i>	15351440
			4	1	90 cm	<i>violet</i>	15401440
			5	2	90 cm	<i>violet</i>	15501440
	HS 48 mm		3.5	0	90 cm	<i>violet</i>	15351448
			4	1	90 cm	<i>violet</i>	15401448
			5	2	90 cm	<i>violet</i>	15501448

 Reverse cutting			Metric	USP	Length		Order number
3/8 circle 	DS 12 mm		0.7	6/0	45 cm	<i>violet</i>	15071512
			1	5/0	45 cm	<i>violet</i>	15101512
			1.5	4/0	45 cm	<i>violet</i>	15151512
	DS 16 mm		0.7	6/0	45 cm	<i>violet</i>	15071516
			1	5/0	75 cm	<i>violet</i>	15101516
			1.5	4/0	75 cm	<i>violet</i>	15151516
			2	3/0	45 cm	<i>violet</i>	15201516-45
	DS 16 mm Easy Pass		2	3/0	50 cm	<i>undyed</i>	15201516UP



Reverse cutting

3/8 circle



		Metric	USP	Length		Order number
DS 16 mm black		0.7	6/0	45 cm	<i>violet</i>	15071516VB
DS 19 mm		1	5/0	75 cm	<i>violet</i>	15101519
		1.5	4/0	75 cm	<i>violet</i>	15151519
		2	3/0	75 cm	<i>violet</i>	15201519
DS 19 mm Easy Pass		1.5	4/0	50 cm	<i>undyed</i>	15151519UP
		2	3/0	50 cm	<i>undyed</i>	15201519UP
DS 19 mm black		1	5/0	75 cm	<i>violet</i>	15101519VB
		1.5	4/0	75 cm	<i>violet</i>	15151519VB
DS 24 mm		1.5	4/0	75 cm	<i>violet</i>	15151524
		2	3/0	75 cm	<i>violet</i>	15201524
		3	2/0	75 cm	<i>violet</i>	15301524
		3.5	0	75 cm	<i>violet</i>	15351524
DS 24 mm Easy Pass		2	3/0	50 cm	<i>undyed</i>	15201524UP
DS 26 mm		2	3/0	75 cm	<i>violet</i>	15201526
		3	2/0	75 cm	<i>violet</i>	15301526
DS 26 mm Easy Pass		2	3/0	75 cm	<i>undyed</i>	15201526UP
		3	2/0	75 cm	<i>undyed</i>	15301526UP-75
DS 30 mm		2	3/0	75 cm	<i>violet</i>	15201530
		3	2/0	75 cm	<i>violet</i>	15301530
		3	2/0	90 cm	<i>violet</i>	15301530-90
		3.5	0	75 cm	<i>violet</i>	15351530
		3.5	0	90 cm	<i>violet</i>	15351530-90
		4	1	90 cm	<i>violet</i>	15401530-90
DS 30 mm Easy Pass		3	2/0	75 cm	<i>undyed</i>	15301530UP-75
DS 36 mm		3	2/0	75 cm	<i>violet</i>	15301536
		3.5	0	75 cm	<i>violet</i>	15351536
		4	1	75 cm	<i>violet</i>	15401536
DS 36 mm Easy Pass		2	3/0	75 cm	<i>violet</i>	15201536P
DS 40 mm		3.5	0	75 cm	<i>violet</i>	15351540
		4	1	75 cm	<i>violet</i>	15401540
DS 45 mm		3	2/0	75 cm	<i>violet</i>	15301545
DS 80 mm		7	5	150 cm	<i>violet</i>	15701580



Round bodied taper cutting point

Metric

USP

Length

Order number

1/2 circle



HRT 13 mm		1.5	4/0	75 cm	<i>violet</i>	15150513
HRT 17 mm		1.5	4/0	75 cm	<i>violet</i>	15150517
		2	3/0	75 cm	<i>violet</i>	15200517
HRT 26 mm		1.5	4/0	75 cm	<i>violet</i>	15150526
		2	3/0	75 cm	<i>violet</i>	15200526
		3	2/0	75 cm	<i>violet</i>	15300526
HRT 30 mm		3	2/0	75 cm	<i>violet</i>	15300530
		3.5	0	75 cm	<i>violet</i>	15350530
		4	1	75 cm	<i>violet</i>	15400530
HRT 36 mm		3	2/0	90 cm	<i>violet</i>	15300536
		3.5	0	90 cm	<i>violet</i>	15350536
		4	1	90 cm	<i>violet</i>	15400536
HRT 40 mm		3.5	0	90 cm	<i>violet</i>	15350540
		4	1	90 cm	<i>violet</i>	15400540
		5	2	90 cm	<i>violet</i>	15500540
HRT 48 mm		8	6	150 cm	<i>violet</i>	15800548



Round bodied taper cutting point

Metric

USP

Length

Order number

3/8 circle



DRT 17 mm		1.5	4/0	75 cm	<i>violet</i>	15150417
		2	3/0	75 cm	<i>violet</i>	15200417



Round bodied taper cutting point

Metric

USP

Length

Order number

3/8 circle




1/2 circle



DRT 13 mm		1	5/0	75 cm	<i>violet</i>	15100530-0413
HRT 30 mm						

 Cutting		Metric	USP	Length		Order number
Straight 	GS 60 mm 	2	3/0	75 cm	<i>violet</i>	15202160
		3	2/0	75 cm	<i>violet</i>	15302160

 Blunt		Metric	USP	Length		Order number
1/2 circle 	HRB 36 mm 	3.5	0	90 cm	<i>violet</i>	15351836
		4	1	90 cm	<i>violet</i>	15401836
	HRB 40 mm 	4	1	90 cm	<i>violet</i>	15401840
	HRB 45 mm 	5	2	90 cm	<i>violet</i>	15501845

 Blunt		Metric	USP	Length		Order number
3/8 circle 	DRB 63 mm 	4	1	100 cm	<i>violet</i>	15401963
	DRB 85 mm 	5	2	75 cm	<i>violet</i>	15501985

IN CASSETTES

Use within 4 months after first opening

Metric	USP	Length	Order number
1.5	4/0	15 m	17151
2	3/0	15 m	17201
2	3/0	25 m	17202
3	2/0	15 m	17301
3	2/0	25 m	17302
3	2/0	50 m	17305
3.5	0	15 m	17351
3.5	0	25 m	17352
3.5	0	50 m	17355
4	1	15 m	17401
4	1	25 m	17402
4	1	50 m	17405
5	2	15 m	17501
5	2	25 m	17502
5	2	50 m	17505
5	2	100 m	175010
6	3&4	15 m	17601
6	3&4	25 m	17602
6	3&4	50 m	17605
7	5	15 m	17701
7	5	25 m	17702
7	5	50 m	17705
8	6	15 m	17801
8	6	25 m	17802
8	6	50 m	17805

